1-Instruction for use	DEDICAM Products (prosthetics, healing caps and impression posts)	p7 to p60
2-Istruzioni per l'uso	Prodotti DEDICAM (protesi, cappette di guarigione e transfer da impronta)	p61 to p118
3-Brugsanvisning	DEDICAM-produkter (proteser, helingshætter og aftryksstifte)	p119 to p173
4-Gebrauchsanweisung	DEDICAM Produkte (Prothetik, Gingivaformer und Abformpfosten)	p174 to p250

This instruction for use contains information on / Le presenti istruzioni per l'uso contengono informazioni su / Denne brugsanvisning indeholder oplysninger om/ Diese Gebrauchsanweisung enthält Informationen zu:

1	General notes on product design	7
2	Product descriptions, including intended purpose/intended use, indications/contraindications and application	7
2.1	One-piece abutments and healing caps made of titanium for Camlog and BioHorizons Implant systems.	7
2.2	One-piece abutments and healing caps for other implant systems	11
2.4	One-piece healing caps and impression posts made of PEEK for Camlog and BioHorizons Implant systems	15
2.5	Mesostructures and crowns for bonding bases (on single implants).....	18
2.6	Bridges, bars, and implant base structures at implant level (screw-retained)	22
2.7	Bridges, bars and implant base structures at abutment level (screw-retained)..	26
2.8	Bridges, bars, and implant base structures for bonding bases (screw-retained)	30
2.9	Suprastructures for bars	35
2.10	Crowns and bridges (bonded)	36
2.11	Double crowns	39
2.12	Veneers, inlays, onlays and partial crowns	41
3	Compatible products for DEDICAM Implant prosthetics.....	42
3.1	Scanbodies for Camlog implant systems	42
3.2	Scanbodies for other implant systems and Multi-unit abutments	43
3.3	Abutment screws.....	43
3.4	Prosthetic screws for CAMLOG and CONELOG Bar abutments.....	44
3.5	Prosthetic screw for Camlog and BioHorizons Multi-unit abutments.....	45
3.6	Prosthetic screws for Nobel Biocare Multi-unit abutments	45
3.7	Lab screws for Camlog implant systems.....	45
3.8	Screwdriver for abutment, prosthetic and lab screws for Camlog implant systems.....	46
3.9	Torque wrench for screws for Camlog implant systems.....	46

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Instruction for use – DEDICAM Products (prosthetics, healing caps and impression posts)

1 General notes on product design

- All DEDICAM Products described in these instructions for use are custom-made devices designed virtually using suitable dental CAD software. They are manufactured individually for each patient in the Biotech Dental Digital Manufacturing center and are therefore intended for single use only.
- Custom-made devices are made under the responsibility of a health care professional and are manufactured in accordance with their written prescription. The prescription provides patient specific design characteristics.
- If not otherwise noted, for implant-supported structures the CAD system specific DEDICAM CAD library must be used. Those structures can only be used for the implants, prefabricated abutments or bonding bases assigned via the CAD library.
- The design of healing caps, impression posts, abutments, bars and crowns/bridges on implants must be fit adequate to the emergence profile to avoid too much pressure of the gingiva.

The following applies to all DEDICAM Products as a matter of principle:

- The prosthetic and physiological principles must be observed in the design. In particular, the lower mechanical load-bearing capacity of diameter-reduced implants must be taken into account.
- Scanning a wax-up as a design template is possible, but the final product does not offer the same surface quality as a purely virtual design.

2 Product descriptions, including intended purpose/intended use, indications/contraindications and application

Note: Gingiva former and healing abutment are synonyms to the term healing cap used in this document.

2.1 One-piece abutments and healing caps made of titanium for Camlog and BioHorizons Implant systems.

2.1.1 Product description

DEDICAM One-piece abutments for CAMLOG, CONELOG, iSy implants and BioHorizons Internal Hex implant-abutment connections serve as customized connection elements for cemented or adhesive-bonded crowns with a dental implant and bridges with two or more dental implants.

DEDICAM One-piece healing caps for CAMLOG, CONELOG, iSy implants and BioHorizons Internal Hex implant-abutment connections serve as customized elements to shape the implant surrounding soft

tissue.

The abutments and healing caps are manufactured from a titanium alloy blank with a prefabricated implant-abutment interface and - where applicable - color-coded by partial anodization according to the color schema of the dedicated implant system.

The design requires a scan using the CAMLOG, CONELOG, iSy or BioHorizons Scanbodies (see 3.1 “Scanbodies for Camlog implant systems” or 3.2 “Scanbodies for BioHorizons Internal Hex implants) or other scanbodies stored in the DEDICAM CAD library.

For the design of double crowns and for the simultaneous design of abutments and crowns or bridges (“file-splitting”) on iSy Implants, the iSy Scanbody must be used on the implant shoulder during the scan.

For CAMLOG Implants with article number K10xx.xxxx, the Platform Switching option is available via the DEDICAM CAD libraries. Abutments and healing caps designed in this way are reduced in diameter to the level of the implant shoulder compared to products without Platform Switching.

Each abutment and healing cap includes a corresponding abutment screw packaged separately. These can also be reordered individually (see 3.4 “Abutment screws”).

Depending on the CAD system used, a tapped hole M1.4, compatible with the titanium screw M1.4 (REF 33000700) of bredent medical GmbH & Co. KG, can be placed in the abutment for a horizontally screw-retained suprastructure. The tapped hole in the DEDICAM CAD library must be selected for this purpose.

The abutments, healing caps and abutment screws are non-sterile and intended for single use. Abutments and healing caps are available finely milled or polished. For abutments, the polished surface is limited to the subgingival section.

Material:

Titanium alloy (Ti-6Al-4V ELI):

CTE (25 to 500°C): $\sim 9 \times 10^{-6}/K$

2.1.2 Intended purpose/intended use.

2.1.2.1 Abutments

Abutments are intended for the functional and/or esthetic oral rehabilitation of partially or fully edentulous patients in combination with endosseous implants in the maxilla and/or mandible.

2.1.2.2 Healing caps

Healing caps are intended for the conditioning of the soft tissue in combination with endosseous implants in the maxilla and/or mandible.

2.1.3 Indications

DEDICAM One-piece abutments and healing caps are components compatible with dental implants for the restoration of partially or fully edentulous jaws with prosthetic restorations such as implant-supported single crowns, bridges and full dentures. Specifically, for:

- single-tooth gaps,
- partially edentulous jaws with several missing teeth or
- edentulous jaws.

2.1.4 Contraindications

The contraindications of the respective implant systems must be considered. These can be found in the respective instructions for use of the implants.

2.1.4.1 Abutments

- Clinical situation which does not permit compliance with design specifications (see 2.1.5.1 “Design”)
- Reuse
- CAMLOG Abutments: Platform Switching on CAMLOG Implants with article number

J10xx.xxxx(no longer available and succeeded by implants with article number K10xx.xxxx)

- Restorations with a length relationship to the implant length greater than 1:1.25
- Primary splinted abutments
- Double crowns on implants with diameter 3.3mm
- Single-tooth restorations with free-end pontic
- Extension attachments
- Parafunctions, e.g. bruxism
- Direct laser welding
- Cast-on technique.
- Modification of the implant-abutment connection incl. the implant shoulder support surface
- Direct full ceramic veneer
- Single crown restorations for implant diameter 3.0 and 3.3mm outside the region of the upper lateral incisors and lower central and lateral incisors

2.1.4.2 Healing caps

- Clinical situation which does not permit compliance with design specifications (see 2.1.5.1 “Design”)
- Reuse
- Modification of the implant-abutment connection incl. the implant shoulder support surface
- CAMLOG Healing caps: Platform Switching on CAMLOG Implants with article number J10xx.xxxx(no longer available and succeeded by implants with article number K10xx.xxxx)

2.1.5 Application

2.1.5.1 Design

The following requirements must be met when designing the abutment or healing cap:

- Angle of the prosthetic axis to the screw axis for abutments up to max. 30°
- Wall thickness in the load-bearing regions is min. 0.6mm
- Radius of curves is min. 0.5mm
- Construction height above the screw head for the abutment is min. 0.2mm and min. 1.0mm for the healing cap.
- Apical minimum diameter according to Table 1


Figure 1: Design limits abutment


Figure 2: Design limits healing cap

Caution:

- To prevent notching effects and the associated risk of fracture of the abutment or healing cap, sharp transitions on the individual geometry must be avoided. If necessary, the abutment or healing cap must be discarded or shapes that may cause notching must be eliminated by polishing under microscopic examination.
- The lower mechanical load-bearing capacity of diameter-reduced implants must be taken into account.

Minimum diameter of the implant support area to be considered for the option of Platform Switching:

CAMLOG Implant Ø [mm]	Minimum shoulder Ø of the healing cap for Platform Switching [mm] (One must not fall below this value. If the diameter does not meet the minimum requirement, the gingiva may be squashed when inserting an abutment for Platform Switching)	Minimum shoulder Ø of the abutments for Platform Switching [mm] (One must not fall below this value. If the diameter does not meet the minimum requirement, the internal configuration of the implant is no longer. covered)
3.3	Platform Switching not possible as the minimum shoulder Ø is 3.3mm	
3.8	3.25	3.2
4.3	3.75	3.7
5.0	4.35	4.3
6.0	5.05	5.0

Table 1: Minimum diameter of the implant support area for the option of Platform Switching

Caution: The option of Platform Switching may only be used in conjunction with CAMLOG Implants with article number K10xx.xxxx.

2.1.6 Post-processing

Caution:

- The contact surfaces between the abutment or healing cap to the implant must not be abrasive blasted or processed mechanically. In this context, please note that discolorations of the titanium alloy can already occur at temperatures above 300°C. Heating to above 800°C leads to structural changes with a reduction in strength.
- The surface quality and wall thickness may not be changed or reduced in the load-bearing regions.

If post-processing of the abutment or the healing cap is required, please note that:

- the implant-abutment connection incl. the implant shoulder support surface is to be protected using an implant analog.
- with Camlog implant systems fixation is performed with the corresponding lab screw (see 3.8 “Labscrews for the Camlog implant systems”).
- that the design specifications are complied with (see 2.1.5.1 “Design”)

2.1.7 Handover to clinician

The post-processed abutment or healing cap must be completely cleaned of processing residues and other contamination using a suitable method.

The specific abutment screw included with the abutment and healing cap must be clearly assigned and passed to the recipient in its original unopened packaging.

2.1.8 Insertion

Caution: To avoid injuries to the gingiva, the following points must be observed during further treatment:

- When using a CAMLOG Healing cap for Platform Switching, only platform switching components (e.g. PS impression posts, PS abutments) must be used.
- When using a customized CONELOG or iSy Healing cap, the use of a scanbody, an impression post or a directly screw-retained bridge, bar or implant base structure is only possible if the implant shoulder is not covered by tissue. For customized BioHorizons Healing caps, this information must be observed for directly screw-retained bridges, bars or implant base structures.

Definitive insertion of the abutment or healing cap in the implant is performed after cleaning, disinfection and recommended sterilization (necessary in the case of healing caps and abutments if not set after completed healing phase, or if set subcrestal) as follows:

1. Remove implant cover screw or healing cap.
2. Clean the inner configuration of the implant and the implant shoulder in the case of CAMLOG ^{impls}.
3. Ensure that the contact surfaces from implant to abutment or healing cap are absolutely clean and, in particular, that bone remains, and soft tissue are completely removed.
4. Insert the abutment or healing cap in the implant.
5. Screw in the new, unused abutment screw supplied with the abutment or healing cap. Use the screwdriver of the corresponding implant system for this purpose.
6. Tighten the abutment screws for fixation of the abutments with a definitive torque of 20Ncm (Camlog ^{impls}) or 30Ncm (BioHorizons Implants).

Warning:

- Lab screws must not be used on the patient.
- The final fixation of the abutment must be performed with new unused abutment screws.
- Do not mix up abutment screws from different implant systems.

Caution:

- If the torque exceeds the recommended value, this may result in damage/destruction of the screw and/or damage to the abutment and/or implant.
 - If the torque lies below the recommended value, this can lead to loosening of the abutment.
 - Only hand-tightened abutment screws for fixation of healing caps.
7. Make sure of a proper fit in the implant. A control x-ray is advisable to make sure the abutment respectively the healing cap is correctly seated on the implant. No soft tissue may be trapped and there must be no visible gap between the abutment or healing cap and the implant.
 8. Retighten with the same torque after at least 5 minutes.

Warning:

- For healing caps, the screw channel must be sealed with a suitable and removable material to secure the only hand-tightened abutment screw.
 - Beforehand, cover the screw head with a completely removable material such as cotton or gutta-percha.
 - This procedure is to be observed for hygienic reasons in the case of abutments.
9. In the case of prosthetic immediate restoration of the implant: the wound margins are closed tightly with atraumatic suture material. Do not tie the sutures too tightly. The gingiva must adapt tightly against the abutment or healing cap in transgingival healing.

2.2 One-piece abutments and healing caps for other implant systems

2.2.1 Product description

DEDICAM One-piece abutments for implant systems of other selected manufacturers serve as customized connection elements for cemented or adhesive-bonded crowns with a dental implant and bridges with two or more dental implants.

DEDICAM One-piece healing caps for implant systems of other selected manufacturers serve as customized elements to shape the implant surrounding soft tissue.

They can only be used for the implant assigned via the CAD library. The selection of compatible implant systems is listed in Table 4 under 3.4 "Abutment screws".

The abutments and healing caps are manufactured from a titanium alloy blank with prefabricated implant-abutment interface and are not color-coded.

The design requires a scan using the DEDICAM Scanbodies (see 3.3 "Scanbodies for other implant
Version 04

systems”) or other scanbodies stored in the DEDICAM CAD library.

Each abutment and healing cap includes a corresponding abutment screw packaged separately. These are compatible with the corresponding screwdriver of the implant manufacturer and can also be ordered separately (see 3.4 “Abutment screws”).

Depending on the CAD system used, a tapped hole M1.4, compatible with the titanium screw M1.4 (REF 33000700) of bredent medical GmbH & Co. KG, can be placed in the abutment for a horizontally screw-retained suprastructure. The tapped hole in the DEDICAM CAD library must be selected for this purpose.

The abutments, healing caps and abutment screws are non-sterile and intended for single use. Abutments and healing caps are available finely milled or polished. For abutments, the polished surface is limited to the subgingival section.

Material:

Titanium alloy (Ti-6Al-4V ELI)

CTE (25 to 500°C): $\sim 10.3 \times 10^{-6}/K$

Depending on the implant system and implant diameter, some abutment screws have a carbon (DLC) or titanium nitride (TiN) coating

2.2.2 Intended purpose/intended use.

2.2.2.1 Abutments

Abutments are intended for the functional and/or esthetic oral rehabilitation of partially or fully edentulous patients in combination with endosseous implants in the maxilla and/or mandible.

2.2.2.2 Healing caps

Healing caps are intended for the conditioning of the soft tissue in combination with endosseous implants in the maxilla and/or mandible.

2.2.3 Indications

DEDICAM One-piece abutments and healing caps are components compatible with dental implants for the restoration of partially or fully edentulous jaws with prosthetic restorations such as implant-supported single crowns, bridges, and full dentures. Specifically, for:

- single-tooth gaps,
- partially edentulous jaws with several missing teeth or
- edentulous jaws.

2.2.4 Contraindications

Caution: The contraindications of the implant to be restored in the respective manufacturer’s instructions must be observed in addition to those listed here.

2.2.4.1 Abutments

- Clinical situation which does not permit compliance with design specifications (see 2.3.5.1 “Design”)
- Reuse
- Use on an incompatible implant.
- Restorations with a length relationship to the implant length greater than 1:1.25
- Primary splinted abutments
- Double crowns on diameter-reduced implants
- Single-tooth restorations with free-end pontic
- Parafunctions, e.g., bruxism
- Direct laser welding
- Cast-on technique.
- Modification of the implant-abutment connection incl. the implant shoulder support surface
- Ceramic veneers

2.2.4.2 Healing caps

- Clinical situation which does not permit compliance with design specifications (see 2.3.5.1 “Design”)
- Reuse
- Use on an incompatible implant.
- Modification of the implant-abutment connection incl. the implant shoulder support surface

Caution:

To avoid injuries to the gingiva, the following point must be observed during further treatment with conical implant-abutment connections:

When using a customized healing cap, the use of a scanbody, an impression post or a directly screw-retained bridge, bar or implant base structure is only possible if the implant shoulder is not covered by tissue.

2.2.5 Application

2.2.5.1 Design

The following requirements must be met when designing the abutment or healing cap:

- Angle of the prosthetic axis to the screw axis for abutments up to max. 30°
- Wall thickness in the load-bearing regions is min. 0.6mm
- Radius of curves is min. 0.5mm
- Construction height above the screw head for the abutment is min. 0.2mm and min. 1.0mm for the healing cap.


Figure 3: Design limits abutment


Figure 4: Design limits healing cap

Caution:

To prevent notching effects and the associated risk of fracture of the abutment or healing cap, sharp transitions on the individual geometry must be avoided. If necessary, the abutment or healing cap must be discarded or shapes that may cause notches must be eliminated by polishing under microscopic examination.

The lower mechanical load-bearing capacity of diameter-reduced implants must be taken into account.

2.2.5.2 Post-processing

Caution:

- The contact surfaces between the abutment or healing cap to the implant must not be abrasive blasted or processed mechanically. In this context, please note that discolorations of the titanium alloy can already occur at temperatures above 300°C. Heating to above 800°C leads to structural changes with a reduction in strength.
- The surface quality and wall thickness may not be changed or reduced in the load-bearing regions.

If post-processing of the abutment or the healing cap is required, please note that:

- the implant-abutment connection incl. the implant shoulder support surface is to be protected using an implant analog
- fixation is performed with the corresponding lab screw.
- that the design specifications are complied with (see 2.3.5.1 “Design”)

2.2.5.3 Handover to clinician

The post-processed abutment or healing cap must be completely cleaned of processing residues and other contamination using a suitable method.

The specific abutment screw included with the abutment and healing cap must be clearly assigned and passed to the recipient in its original unopened packaging.

2.2.5.4 Insertion

Definitive insertion of the abutment or healing cap in the implant is performed after cleaning, disinfection and recommended sterilization (necessary in the case of healing caps and abutments if not set after completed healing phase, or if set subcrestal) as follows:

1. Remove implant cover screw or healing cap.
2. Clean the inner configuration of the implant and - in case of butt-joint implant-abutment connections - the implant shoulder.
3. Ensure that the contact surfaces from implant to abutment or healing cap are absolutely clean and, in particular, that bone remains and soft tissue are completely removed.
4. Insert the abutment or healing cap in the implant.
5. Screw in the new, unused DEDICAM Abutment screw supplied with the abutment or healing cap. Use the screwdriver of the corresponding implant manufacturer for this purpose.
6. Tighten the abutment screw with the specified tightening torque from Table 4 under 3.4 “Abutment screws”, whereby the value is to be compared with the respective specification in the implant manufacturer's user information and reduced if necessary.

Caution:

- If the torque exceeds the recommended value, this may result in damage/destruction of the screw and/or damage to the abutment and/or implant.
 - If the torque lies below the recommended value, this can lead to loosening of the abutment.
7. Only hand-tighten abutment screws for fixation of healing caps.
 8. Make sure of a proper fit in the implant. A control x-ray is advisable to make sure the abutment resp. the healing cap is correctly seated on the implant. No soft tissue may be trapped and there must be no visible gap between the abutment or healing cap and the implant.
 9. Retighten with the same torque after at least 5 minutes.

Warning:

- For healing caps, the screw channel must be sealed with a suitable and removable

material to secure the only hand-tightened abutment screw. Beforehand, cover the screw head with a completely removable material such as cotton or gutta-percha. This procedure is recommended for hygienic reasons in the case of abutments.

10. In the case of prosthetic immediate restoration of the implant: the wound margins are closed tightly with atraumatic suture material. Do not tie the sutures too tightly. The gingiva must adapt tightly against the abutment or healing cap in transgingival healing.

2.4 One-piece healing caps and impression posts made of PEEK for Camlog and BioHorizons Implant systems

2.4.1 Product description

DEDICAM One-piece healing caps made of PEEK for CAMLOG, CONELOG, iSy, BioHorizons Tapered internal and BioHorizons Internal implants serve as customized elements to shape the implant surrounding soft tissue.

DEDICAM Impression posts made of PEEK for CAMLOG, CONELOG, iSy, BioHorizons Tapered internal and BioHorizons Internal implants serve as tools to transfer the oral implant position and the profile of the previously formed peri-implant soft tissue to the master cast. During impression taking, the impression post supports the peri-implant soft tissue to avoid its collapse.

The design requires a scan using the CAMLOG, CONELOG, iSy or BioHorizons Scanbodies (see 3.1 “Scanbodies for the Camlog implant systems” or 3.2. “Scanbodies for BioHorizons Internal implants”) or other scanbodies stored in the DEDICAM CAD library.

The healing caps and impression posts sit on the implant shoulder in each case.

Impression posts are available for both the open tray and (not for BioHorizons Implants) the closed tray.

For CAMLOG Implants with article number K10xx.xxxx, the Platform Switching option is available via the DEDICAM CAD libraries. Healing caps and impression posts designed in this way are reduced in diameter to the level of the implant shoulder compared to products without Platform Switching.

Each healing cap and each impression post includes a corresponding screw packaged separately. Impression posts for closed trays are supplied with a separately packaged impression cap.

The healing caps, impression posts, corresponding screws and impression caps are non-sterile and intended for single use.

Material:

PEEK (polyetheretherketone)

2.4.2 Intended purpose/intended use

2.4.2.1 Healing caps

Healing caps are intended for the conditioning of the soft tissue in combination with endosseous implants in the maxilla and/or mandible.

2.4.2.2 Impression posts

Impression posts are intended for the three-dimensional registration of implant positions in the oral cavity.

2.4.3 Indications

DEDICAM One-piece healing caps and impression posts are components compatible with dental

implants for the restoration of partially or fully edentulous jaws with prosthetic restorations such as implant-supported single crowns, bridges and full dentures. Specifically for:

- single-tooth gaps,
- partially edentulous jaws with several missing teeth or
- edentulous jaws.

2.4.4 Contraindications

2.4.4.1 Healing caps

- Clinical situation which does not permit compliance with design specifications (see 2.4.5.1 “Design”)
- Reuse
- Modification of the implant-abutment connection incl. the implant shoulder support surface
- CAMLOG healing caps: Platform Switching on CAMLOG Implants with article number J10xx.xxxx(no longer available and succeeded by implants with article number K10xx.xxxx)
- A wear period of over 180 days

Caution: To avoid injuries to the gingiva, the following points must be observed during further treatment:

- When using a CAMLOG Healing cap for Platform Switching, only platform switching components(e.g. PS impression posts, PS abutments) must be used.
- When using a customized CONELOG or iSy Healing cap, the use of a scanbody, a standard impression post or a directly screw-retained bridge, bar or implant base structure is only possible if the implant shoulder is not covered by tissue. For customized BioHorizons Healing caps, this information must be observed for directly screw-retained bridges, bars and implant base structures.

2.4.4.2 Impression posts

- Reuse
- Modification of the implant-abutment connection incl. the implant shoulder support surface

Caution: To avoid injuries to the gingiva, the following point must be observed during further treatment: When using a customized CONELOG, iSy or BioHorizons Impression post no stock healing cap should have been used previously as they do not cover the implant shoulder.

2.4.5 Application

2.4.5.1 Design

The following requirements must be met when designing the healing cap:

- The wall thickness is at least 0.6mm
- Avoid sharp edges in the design
- The construction height above the screw head for the healing cap is at least 1.0mm
- Apical minimum diameter according to Table 2
- The individual portion of the healing cap and impression post is maximum 7mm in height
- In the case of the impression post, the individualized portion with a height of less than 7mm is automatically supplemented with a cylinder. The fixed portion is also automatically added to the individualized and, if necessary, cylindrically extended portion to a total height of 12mm.
- A groove on the occlusal surface in alignment with the one characteristic element of the implant-abutment interface, e. g. a cam is recommended for easier intraoral placement.


Figure 5: Construction height limit healing caps
Version 04


Figure 6: Minimum thicknesses of healing caps


Figure 7: Fixed (a) and individualized (b) proportion of impression posts for open trays


Figure 8: Fixed (a) and individualized (b) proportion of impression posts for closed tray

Minimum diameter of the implant support area to be considered for the option of Platform Switching:

CAMLOG Implant Ø [mm]	Minimum shoulder Ø of the healing cap for Platform Switching [mm] (One must not fall below this value. If the diameter does not meet the minimum requirement, the gingiva may be squashed when inserting an abutment for Platform Switching)
3.3	Platform Switching not possible as the minimum shoulder Ø is 3.3mm
3.8	3.2
4.3	3.7
5.0	4.3
6.0	5.0

Table 2: Minimum diameter of the implant support area for the option of Platform Switching

Caution: The option of Platform Switching may only be used in conjunction with CAMLOG Implants with article number K10xx.xxxx.

2.4.5.2 Post-processing

Caution:

- The contact surfaces between the healing cap or impression post to the implant must not be abrasive blasted or processed mechanically.
- The surface quality and wall thickness may not be changed or reduced in the load-bearing regions.

If post-processing of the healing cap or the impression post is required, please note that:

- the implant-abutment connection incl. the implant shoulder support surface is to be protected using an implant analog
- with Camlog implant systems fixation is performed with the corresponding lab screw or impression post screw (see 3.8 "Lab screws for Camlog implant systems")
- that the design specifications are complied with (see 2.4.5.1 "Design")
- that cross-toothed milling cutters suitable for machining resins are used
- in the case of a closed impression, the contact surfaces of the retention portion for the impression caps must not be modified so as not to impair the precision of the impression or the repositioning of the impression post

In the case of additive shape correction of the contact surfaces to the gingiva, please note that:

- only suitable materials such as composites/PMMA are used
- the implant-abutment connection incl. the implant shoulder support surface is to be protected using an implant analog
- with Camlog implant systems fixation is performed with the corresponding lab screw (see 3.8 "Lab screws for the Camlog implant systems") or impression post screw

2.4.5.3 Handover to clinician

The post-processed healing cap or impression post must be completely cleaned of processing residues and other contamination using a suitable method.

The specific screw included with the healing cap or impression post must be clearly assigned and passed to the recipient in its original unopened packaging.

2.4.5.4 Insertion

The healing cap or impression post is inserted into the implant after cleaning, disinfection and sterilization as follows:

1. Remove implant cover screw or implant abutment.
2. Clean the inner configuration of the implant and the implant shoulder in the case of CAMLOG Implants.
3. Ensure that the contact surfaces from implant to the healing cap or impression post are absolutely clean and, in particular, that bone remains, and soft tissue are completely removed.
4. Insert the healing cap or impression post in the implant.
5. Healing cap: screw in the new, unused abutment screw supplied with the healing cap. Use the screwdriver of the corresponding implant system for this purpose.
6. Impression post: screw in the new, unused impression post screw supplied with the impression cap. Use the screwdriver of the corresponding implant system for this purpose or just manually.
7. Only hand-tighten the screws for fixation of healing caps and impression posts.
8. Make sure of a proper fit in the implant.

Warning:

- Impression posts for open trays: the BioHorizons Impression post screw is not secured against dropping out.
- Impression posts for closed trays: the CAMLOG, CONELOG and iSy Impression post screw is not secured against dropping out.
- For healing caps, the screw channel must be sealed with a suitable and removable material to secure the only hand-tightened abutment screw. Beforehand, cover the screw head with a completely removable material such as for example, sterile Teflon tape or gutta-percha for protection and then close the screw channel with a commercially available composite.

2.5 Mesostructures and crowns for bonding bases (on single implants)

2.5.1 Product description

DEDICAM Mesostructures and crowns for bonding bases serve as replacement for lost tooth substance on a single implant. Mesostructures can also be used as bridge supporting pillars. They are bonded with a:

- CAMLOG, CONELOG or iSy Titanium base CAD/CAM, crown, or
- CAMLOG, CONELOG or iSy Titanium base CAD/CAM free, crown, or
- iSy Implant base, or
- BioHorizons Internal Hybrid Abutment Base, hexed, or
- bonding bases from other selected

suppliers to form so-called “hybrid abutment” or “hybrid crown.”

The structures are fitted with a continuous screw hole (also available without screw hole for the iSy Implant base by ordering via “3Shape Inbox”) and an inner geometry corresponding to a bonding base assigned in the CAD system.

CAMLOG, CONELOG or iSy Titanium base CAD/CAM free, crown:

Their design enables a vertical shift of the screw channel in the mesostructure or crown of up to 25° in

relation to the longitudinal axis. Exception: CONELOG and iSy Titanium base with gingival height 2.0mm is limited up to 15°. In addition, a horizontal deviation of the screw channel of up to 35° relative to the vertical cutout edges for the rotation lock is permissible.

Caution: Violation of the horizontal deviation will render the antirotational feature ineffective.

BioHorizons Internal Hybrid Abutment Base, hexed:

Their design enables a vertical shift of the screw channel in the mesostructure or crown of up to 15° in relation to the longitudinal axis. A horizontal deviation of the screw channel relative to the sloped cutout edges is not possible. For vertical shifts between 1 – 15° the use of the special CAD library “DEDICAM BioHorizons Internal hybrid abutment base angulated screw channel, hexed” is mandatory.

The mesostructures and crowns are designed virtually with the aid of suitable CAD software and the DEDICAM CAD library or freely available, unencrypted CAD libraries. They can only be used for the bonding base assigned via the CAD library.

For the design of double crowns and for the simultaneous design of abutments and crowns or bridges (“file-splitting”) on iSy Implants, the iSy Scanbody must be used on the implant shoulder during the scan.

Mesostructures and crowns for bonding bases are supplied non-sterile. Structures on a CAMLOG, CONELOG or iSy Titanium base (regular) or a BioHorizons Internal Hybrid Abutment Base, hexed, are supplied together with the titanium base and the abutment screw. Structures, designed on BioHorizons Internal Hybrid Abutment Base, hexed and with an angulated screw channel, are supplied with the BioHorizons Precision Angled Screw in addition.

The structures are available for individual characterization in different tooth shades. The zirconia “IPS e.max ZirCAD” is also available in various degrees of translucency and optionally with a natural color gradient from dentin to incisal.

Temporary restorations up to a maximum wearing time of 12 months can be made from a tooth-coloring plastic (polymethylmethacrylate).

Mesostructures and crowns for bonding bases are supplied already densely sintered for zirconia and in the “blue” metasilicate phase for lithium disilicate. Lithium disilicate structures are to be subjected to crystallization firing.

Crowns for bonding bases are particularly suited if the screw channel is located occlusally or palatal/lingual.

Material:

IPS e.max ZirCAD (zirconia)	MT and MT Multi	CTE (25 to 500°C): $10.4 \pm 0.5 \times 10^{-6}/K$
	LT	CTE (25 to 500°C): $10.5 \pm 0.5 \times 10^{-6}/K$
IPS e.max CAD (lithium disilicate)*		CTE (100 to 400°C): $\sim 10.2 \times 10^{-6}/K$
		CTE (100 to 500°C): $\sim 10.5 \times 10^{-6}/K$
Telio CAD (polymethylmethacrylate)		CTE: not applicable

* Only for selected bonding bases – see current DEDICAM Catalog.

2.5.2 Intended purpose/intended use

Mesostructures and crowns are intended for the functional and/or esthetic oral rehabilitation of partially or fully edentulous patients as dental prosthesis elements in combination with endosseous implants in the maxilla and/or mandible.

2.5.3 Indications

DEDICAM Mesostructures and crowns are components of a prosthetic solution for the restoration of partially or fully edentulous jaws with prosthetic restorations such as implant-supported single crowns,

bridges and full dentures. Specifically for:

- single-tooth gaps,
- partially edentulous jaws with several missing teeth or
- edentulous jaws.

2.5.4 Contraindications

- Clinical situation which does not permit compliance with design specifications (see 2.5.5.1 “Design”)
- Reuse
- Primary splints
- Restorations with a length relationship to the implant length greater than 1:1.25
- Restorations on bonding bases designed for multi-unit restorations and other incompatible third-party bonding bases
- Double crowns on diameter-reduced implants
- Single-tooth restorations with free-end pontic
- Extension attachments on mesostructures and crowns
- Parafunctions, e.g. bruxism
- A wear period of over 12 months for restorations made of Telio CAD

Caution: The contraindications of the implant to be restored and the bonding bases in the respective manufacturer's instructions must be observed.

2.5.5 Application

2.5.5.1 Design

The following requirements must be met when designing the mesostructures and crowns for bonding bases:

- Angle of the prosthetic axis to screw axis up to 20°
- CAMLOG, CONELOG and iSy Titanium base CAD/CAM free, crown:
 - Vertical shift of the screw channel up to 25° in relation to the longitudinal axis. Exception: CONELOG and iSy Titanium base with gingival height 2.0mm is limited up to 15°.
 - Horizontal shift of the screw channel of up to 35° relative to the vertical cutout edges for the rotation lock
- BioHorizons Internal Hybrid Abutment Base, hexed:
 - Vertical shift of the screw channel up to 15° in relation to the longitudinal axis.
 - A horizontal shift of the screw channel relative to the sloped cutout edges is not possible.
- Wall thicknesses:

	Mesostructure/ crown IPS e.max ZirCAD	Mesostructure IPS e.max CAD	Crown: IPS e.max CAD	Crown: Telio CAD
Minimum wall thickness – circular around screw channel [mm]	0.7 (LT) 1.0 (MT + MT Multi)	0.5	1.5	0.8
Minimum wall thickness – circular on bonding base shoulder [mm]	0.5	0.5	0.5	0.5
Minimum wall thickness – incisal/occlusal [mm]	0.7 (LT) 1.0 (MT + MT Multi)	--	1.5	1.5

Maximum wall thickness – circular around screw channel [mm]	--	6.0	6.0	--
---	----	-----	-----	----

Table 3: Wall thicknesses of mesostructures and crowns for bonding bases

The lower mechanical load-bearing capacity of diameter-reduced implants must be taken into account.

2.5.5.2 Post-processing

Mesostructures and crowns for bonding bases made of lithium disilicate must be crystallized before or after shade customization and try-in, whereby attention should be paid to sufficient support for the structure.

Crowns made of PMMA for bonding bases can be modified to optimize the emergence profile by applying suitable composite or – taking into account the minimum permissible wall thicknesses – by grinding.

Caution: The surface quality and wall thickness may not be changed or reduced in the load-bearing regions.

The processing instructions of the material manufacturers must be observed, these are available under:

- IPS e.max ZirCAD (zirconia): www.ivoclar.com/en_gb/eifu
- IPS e.max CAD (lithium disilicate): www.ivoclar.com/en_gb/eifu
- Telio CAD (polymethylmethacrylate): www.ivoclar.com/en_gb/eifu

After surface conditioning according to the respective specifications of the material and mounting composite manufacturers, the structures are permanently bonded extra orally to the bonding bases prepared for bonding, taking into account the correct alignment and, if necessary, using bonding aids (included component of CAMLOG, CONELOG and iSy Titanium bases CAD/CAM, crown; only suitable for straight screw channels). If a CAMLOG, CONELOG or iSy Titanium base CAD/CAM free, crown, is used in the procedure as given in the “Instruction for use – CAMLOG Titanium bases CAD/CAM” (document ID: J8000.0445), “Instruction for use – CONELOG Titanium bases CAD/CAM” (document ID: J8000.0446) or “Instruction for use – iSy Titanium bases CAD/CAM” (document ID: J8000.0472) respectively then these must be observed. Particular attention should be paid to this: when used with an angulated screw access the new unused abutment screw must be placed before bonding.

For bonding, the structure is rotated on the titanium base until the anti-rotational features engage, then gently pressed onto the titanium base. Excess bonding material must be removed immediately.

Caution: Cement contact with gingiva has been known to cause inflammation.

2.5.5.3 Suitable surface conditioners and fixation composites

To connect the mesostructures or crown to a bonding base, we recommend the following products from Ivoclar Vivadent AG:

- Bonding agent for PMMA build-up: “SR Connect”
- Primer for bonding base and ceramic build-up: “Monobond Plus”
- Fixation composites: “Multilink Hybrid Abutment”

Observe the manufacturer's processing instructions. Biotech Dental Digital SAS accepts no responsibility for the durability of the bond.

2.5.5.4 Handover to clinician

The mesostructure or crown for the bonding base must be completely cleaned of processing residues and other impurities using a suitable procedure after bonding with the bonding base as a hybrid abutment or hybrid crown.

The specific abutment screw supplied with a CAMLOG, CONELOG or iSy Titanium base CAD/CAM, crown, or a BioHorizons Internal Hybrid Abutment Base, hexed, must be clearly assigned and passed to the recipient in its original unopened packaging (not applicable for pre-assembled (“captured”) abutment screws in case of angulated screw channels).

2.5.5.5 Insertion

Definitive insertion of the hybrid abutment or the hybrid crown in the implant is performed after cleaning, disinfection, and recommended sterilization (necessary if not set after completed healing phase, or if set subcrestal) as follows:

1. Remove implant cover screw or healing cap.
2. Clean the inner configuration of the implant and - in case of butt-joint implant-abutment connections - the implant shoulder.
3. Ensure that the contact surfaces from implant to the hybrid abutment or hybrid crown are absolutely clean and, in particular, that bone remains, and soft tissue are completely removed.
4. Insert the hybrid abutment or hybrid crown into the implant.
5. Screw in the new, unused abutment screw corresponding to the implant system and the implant diameter using the corresponding screwdriver. Note: an angulated screw channel requires use of a specific screwdriver. It is recommended to use the regular screwdriver for straight screw channels.
6. Screw the abutment screw tight with the torque specified by the manufacturers of the bonding base and implant.

Caution:

- If the torque exceeds the recommended value, this may result in damage/destruction of the screw and/or damage to the hybrid abutment or hybrid crown and/or implant.
 - If the torque lies below the recommended value, this can lead to loosening of the hybrid abutment or hybrid crown.
7. Make sure of a proper fit in the implant. A control X-ray is recommended to check the correct fit of the hybrid abutment or the hybrid crown in the implant. No soft tissue may be trapped and there must be no visible gap between the abutment and the implant.
 8. Hybrid abutment: for hygiene reasons, the screw channel is to be sealed with suitable removable materials. Beforehand, cover the screw head with a completely removable material such as cotton or gutta-percha.
 9. Hybrid crown: cover the screw head with a suitable and completely removable material. Then seal the screw channel tightly with suitable materials, e.g., a composite.

2.6 Bridges, bars, and implant base structures at implant level (screw-retained)

2.6.1 Product description

DEDICAM Bridges at implant level are customized artificial restorative elements and serve as replacement for lost tooth substance and are screw-retained on two or more implants.

DEDICAM Bars at implant level serve as customized connection elements for removable prosthetic tooth restorations and are screw-retained on two or more implants.

DEDICAM Implant base structures at implant level serve as customized metallic support structures for prosthetic tooth restorations made of prosthetic plastic and prefabricated denture teeth made of plastic or ceramic and are screw-retained on two or more implants.

The structures are screwed directly to the implants without additional abutments. They are available for Camlog implant systems, BioHorizons Implants with the Internal Hex implant-abutment connections as well as for implants from other selected suppliers (see Table 4 under 3.4 "Abutment screws").

In the case of conical implant-abutment connections, the bridges, bars, and implant base structures at implant level are also supported on the implant shoulder. Butt-joint (e. g. CAMLOG) and conical (e. g. CONELOG, iSy or BioHorizons Internal) implant platforms can be combined.

The anatomical bridges and bridge frameworks at implant level can consist of 2 to 16 elements. Bars and implant base structures at implant level are available for 2 up to 10 implants. With the iSy Implant system, bridges and implant base structures must be supported on at least 6 implants, bars on at least 4 implants.

The design requires a scan using the CAMLOG, CONELOG, iSy, BioHorizons or DEDICAM Scanbodies (see 3.1 “Scanbodies for Camlog implant systems,” 3.2 “Scanbodies for BioHorizons Internal implants” or 3.3 “Scanbodies for other implant systems”). The use of the iSy Multifunctional cap during scanning is not recommended for the construction of multi-unit structures at implant level due to the high accuracy required for this purpose.

For CAMLOG Implants with article number K10xx.xxxx, the Platform Switching option is available via the DEDICAM CAD libraries. Structures designed in this way are reduced in diameter to the level of the implant shoulder compared to products without Platform Switching.

The structures are non-sterile and intended for single use. Bridges and implant base structures have a milled surface, bars are also available polished.

The bridges, bars and implant base structures at implant level are supplied together with the corresponding abutment screws, packaged separately, and in a quantity corresponding to the implants to be restored with the structure. The abutment screws can also be reordered individually (see 3.4 “Abutment screws”). Bridges, bars, and implant base structures at implant level are fixated to the implants using the supplied abutment screws and observing the specified torque (see 3.4 “Abutment screws”).

The provided CAD files of bridges, bars and implant base structures at implant level are analyzed for angulation of the implant axes in manufacturing (CAM software) among other things. Hereby, the implant connections are individually optimized to the overall design and the insertion direction so that the guidance and positioning is given the necessary freedom on the one hand, and stability is guaranteed on the other. Due to this method of detecting interference contacts and their subsequent removal through milling, divergences in the implant axes can be compensated by up to 50° depending on the situation and implant system.

The anatomically reduced bridge framework can be veneered with suitable layering materials or with CAD/CAM fabricated suprastructures.

With the aid of the DEDICAM CAD library, bridges and bars can be supplemented with attachments for the fixation of removable restorations and bars can also be designed with different profiles.

A suitable bar profile is available in the DEDICAM CAD library for implant base structures.

Suitable bar supraconstructions are available for DEDICAM Bars specifically designed for this purpose. The design of the bars and bar supraconstructions is performed by the DEDICAM Scan & Design Service.

Material:

Titanium alloy (Ti-6Al-4V)
CoCr alloy

CTE (20 to 600°C): $\sim 10.3 \times 10^{-6}/K^{-1}$
CTE (20 to 500°C): $\sim 14.0 \times 10^{-6}/K^{-1}$

2.6.2 Intended purpose/intended use

Bridges, bars and implant base structures are intended for the functional and/or esthetic oral rehabilitation of partially or fully edentulous patients as dental prosthesis elements in combination with endosseous implants in the maxilla and/or mandible.

2.6.3 Indications

DEDICAM Bridges, bars and implant base structures at implant level are components compatible with dental implants for the restoration of partially or fully edentulous jaws with prosthetic restorations such as implant-supported bridges and full dentures. Specifically for:

- partially edentulous jaws with several missing teeth or

- edentulous jaws.

2.6.4 Contraindications

- Clinical situation which does not permit compliance with design specifications (see 2.2.5.1 “Design”)
- Reuse
- partially edentulous jaws with several missing teeth or
- edentulous jaws.

2.6.5 Contraindications

- Clinical situation which does not permit compliance with design specifications (see 2.6.5.1 “Design”)
- Reuse
- Use on an incompatible implant
- Platform Switching on CAMLOG Implants with article number J10xx.xxxx (no longer available and succeeded by implants with article number K10xx.xxxx)
- Joint support on tooth (stump) and implant
- Parafunctions, e.g. bruxism

Caution:

- The contraindications given in the respective manufacturer's instructions for the implant to be restored must be observed. For the iSy Implant System, please note that bars are to be planned on at least 4 implants, bridges and implant base structures on at least 6 implants, without mesial or distal extensions and with an angulation correction of not more than 20° to the implant axis. Adequate polygonal support needs to be provided (cross arch).
- The use of different alloy groups in the patient's mouth can lead to galvanic effects.

2.6.6 Application

2.6.6.1 Design

General notes on directly screw-retained vs. abutment-supported implant prosthetics:

We recommend checking the individual cases when designing construction. Both restoration variants have advantages and disadvantages which should be assessed depending on the situation.

In the case of the direct screw retention (design-dependent, modified implant-abutment connection), the restoration is only supported via the implant shoulder, and the fixing screw is used for a considerable part of the overall force transmission. Long load arms and uniaxial distribution of the implants should be regarded critically. For this purpose, angulations up to 50° are possible depending on the situation and the implant system, and this saves costs for construction elements.

For directly screw-retained bridges, bars and implant base structures Biotech Dental Digital reserves the right to refuse manufacturing in case of:

- An angulation between two implants greater than 50° for Camlog implants, BioHorizons Tapered internal and Internal implants and depending on the platform, 30° for implant systems of other suppliers.
- A unilateral construction with implants generally inclined (in particular vestibular or buccal).

The mechanical stability of abutment-supported implant prostheses is regarded as being more beneficial.

The following requirements must be met when designing bridges, bars and implant base structures at implant level:

- Minimum wall thickness: 0.5mm
- Minimum connector cross sections:
 - Bridge element/splinted crowns, anterior tooth: 6mm²
 - Bridge element/splinted crowns, posterior tooth: 9mm²
 - Bars with an individualized profile and implant base structures: 8mm²

- Bridges made of titanium alloy: maximum of 3 connected pontics; in the case of sole support on 2 implants, the maximum distance between the implant axes is 30mm
- Bridges made of CoCr alloy: maximum of 3 connected pontics in the lateral region and 4 in the anterior region; with sole support on 2 implants, the maximum distance between the implant axes is 30mm
- Bridges: maximum length of extensions; width of a molar, however, not more than 10mm
- Bars and implant base structures: maximum distance between two implant axes: 30mm
- Bars and implant base structures: maximum length of extensions: half the length of the adjacent inter-implant segment, however, not more than 10mm
- Indications and contraindications for the implants according to the manufacturer's instructions are to be taken into consideration
- With the iSy Implant system, bridges and implant base structures must be supported on at least 6 implants, bars on at least 4 implants

The lower mechanical load-bearing capacity of diameter-reduced implants must be taken into account. To facilitate the implementation of veneered bridges from the model in the sintering furnace, lateral retaining pins with a minimum diameter of 1mm and a maximum length of 5mm can be designed on the CAD model.

Various bar profiles and attachments are available as design options in the DEDICAM CAD library.

2.6.6.2 Post-processing

If no information on the surface conditioning of the bridge frameworks/implant base structures for the selected veneering/construction material is available from the respective manufacturer, bridge frameworks and implant base structures made of CoCr alloy are to be blasted with pure aluminum oxide with a grain diameter of minimum 110µm – 250µm and with a blasting pressure of 2 to 4 bar, bridge frameworks and implant base structures made of titanium alloy with pure aluminum oxide with a grain diameter of 110µm – 250µm and a blasting pressure of 2 to 3 bar.

Warning:

- During machining or thermal processing of CoCr material (e.g., milling, grinding) dust and fumes can be produced which could entail a health risk if inhaled. Therefore, suitable measures must be applied to reduce dust and vapor exposure, such as the use of extraction and ventilation systems or, for short-term use, suitable personal protective equipment for respiratory air filtration.

The bridge frameworks and implant base structures must then be washed under running water, cleaned with a steam jet and degreased with ethyl alcohol (for titanium alloy never use for cleaning hydrofluoric acid). Do not contaminate the structure afterwards, e. g. touch it with unprotected fingers. Bridge frameworks are veneered with the layer material suitable for the material used (titanium or CoCr alloy) in each case or with CAD/CAM fabricated suprastructures. The manufacturer's processing instructions are to be observed here. In this context, please consider that discolorations of the titanium alloy can already occur at temperatures above 300°C.

Caution:

- The contact surfaces between the bridge, bar or implant base structure to the implant must not be abrasive blasted or processed.
- Heating of titanium to above 800°C leads to structural changes with a reduction in strength.

The processing instructions of the material manufacturers must be observed, these are available under: CoCr alloy (Starbond CoS Disc Basic):

Titanium alloy (Starbond Ti5 Disc):

Please contact via email at service@scheftner.dental

Please contact via email at service@scheftner.dental

To avoid aesthetic disadvantages, CoCr restorations veneered with a layering material must have the

oxide layer removed from the surfaces that come into contact with the soft tissue.

2.6.6.3 Handover to clinician

The post-processed structure must be completely cleaned of processing residues and other contamination using a suitable method.

The specific abutment screws supplied together with the bridges, bars and implant base structures at implant level must be passed on to the recipient in the unopened original packaging and clearly assignable.

2.6.6.4 Insertion

Bridges, bars and implant base structures at implant level are definitively inserted into the implant after cleaning, disinfection and recommended sterilization (necessary if not set after completed healing phase, or if set subcrestal) as follows:

1. Remove implant cover screw or healing cap.
2. Clean the inner configuration of the implant and the implant shoulder.
3. Ensure that the contact surfaces from implant to bridge, bar or implant base structure are absolutely clean and, in particular, that bone remains and soft tissue are completely removed.
4. Insert bridge, bar or implant base structure into implants.
5. Screw in new, unused abutment screws supplied with the structure. Use the screwdriver of the corresponding implant system for this purpose.
6. Tighten the abutment screws cross-wise with the specified tightening torque from Table 4 under 3.4 "Abutment screws" whereby the value is to be compared with the respective specification in the implant manufacturer's user information and reduced if necessary.

Caution:

- If the torque exceeds the recommended value, this may result in damage/destruction of the screw and/or damage to the bridge, the bar or implant base structure and/or the implant.
 - If the torque lies below the recommended value, this can lead to loosening of the bridge, the bar or the implant base structure.
7. Make sure of a proper fit in the implant. A control X-ray is recommended to check the correct fit of the bridge, the bar or the implant base structure in the implant. No soft tissue may be trapped and there must be no visible gap between the bridge, bar or implant base structure and the implant.
 8. Retighten with the same torque after at least 5 minutes.
 9. Bridges and implant base structures: cover the screw head with a completely removable material such as cotton or gutta-percha for protection. Then seal the screw channel tightly with suitable materials, e.g. a composite.
 10. Bars: for hygiene reasons, the screw channels are to be sealed with suitable and removable materials. Beforehand, cover the screw head with a completely removable material such as cotton or gutta-percha.

2.7 Bridges, bars and implant base structures at abutment level (screw-retained)

2.7.1 Product description

DEDICAM Bridges at abutment level are customized artificial restorative elements and serve as replacement for lost tooth substance and are screw-retained on two or more abutments.

DEDICAM Bars at abutment level serve as customized connection elements for removable prosthetic tooth restorations and are screw-retained on two or more abutments.

DEDICAM Implant base structures at abutment level serve as customized metallic support structures for prosthetic tooth restorations made of prosthetic plastic and prefabricated denture teeth made of plastic or ceramic and are screw-retained on two or more abutments.

The structures are screw-retained occlusally on prefabricated abutments offered specifically for this application. No intermediate elements requiring bonding are used. They are available for CAMLOG and CONELOG Bar abutments as well as for iSy Multi-unit abutments from Camlog or BioHorizons and Nobel Biocare.

In the case of conical implant-abutment connections, the bridges, bars and implant base structures at implant level are also supported on the implant shoulder. Butt-joint (e. g. CAMLOG) and conical (e. g. CONELOG, iSy or BioHorizons Internal) implant platforms can be combined.

The anatomical bridges and bridge frameworks at abutment level can consist of 2 to 16 elements. Bars and implant base structures at abutment level are available for 2 to 10 abutments.

The design of bridges, bars and implant base structures at abutment level requires performing the scan at abutment level. Scanning caps for bar abutments for CAMLOG and CONELOG Bar abutments (see **3.1.3** “Scanning caps for bar abutments”) are available for this purpose, and abutment-specific DEDICAM Scanbodies (see **3.1.5** “Scanbodies for other implant systems and Multi-unit abutments”) for the Camlog, BioHorizons and Nobel Biocare Multi-unit abutments.

The anatomically reduced bridge framework can be veneered with suitable layering materials or with CAD/CAM fabricated suprastructures.

With the aid of the DEDICAM CAD library, bridges and bars can be supplemented with attachments for the fixation of removable restorations and bars can also be designed with different profiles. A suitable bar profile is available in the DEDICAM CAD library for implant base structures.

The structures are non-sterile and intended for single use.

Bridges and implant base structures have a milled surface, bars are also available polished.

Optional titanium bonding bases for bar abutments are available for the CAMLOG and CONELOG Implant System as well as for Camlog and BioHorizons Multi-unit abutments (see 2.9 "Bridges, bars and implant base structures for bonding bases (screw-retained)"). In this option, the bridge, the bar or the implant base structure are bonded intraorally to the titanium bonding bases screw-retained to the bar abutments or Multi-unit abutments. This allows compensating possible fit imprecisions (Passive Fit) which could lead to tension.

For bridges, bars and implant base structures at abutment level, the corresponding prosthetic screws are supplied separately packed, but can also be ordered separately (see **3.4** “Prosthetic screws for CAMLOG and CONELOG Bar abutments” , **3.6** “Prosthetic screws for Camlog and BioHorizons Multi-unit abutments” or **3.7** “Prosthetic screws for Nobel Biocare Multi-unit abutments”). The fixation of the structures on the abutments is performed under observation of the specified torque for the prosthetic screws.

Suitable bar supraconstructions are available for DEDICAM Bars specifically designed for this purpose. The design of the bars and bar supraconstructions is performed by the DEDICAM Scan & Design Service.

Material:

Titanium alloy (Ti-6Al-4V)

CTE (20 to 600°C): $\sim 10.3 \times 10^{-6}/K^{-1}$

CoCr alloy

CTE (20 to 500°C): $\sim 14.0 \times 10^{-6}/K^{-1}$

2.7.2 Intended purpose/intended use

Bridges, bars and implant base structures are intended for the functional and/or esthetic oral rehabilitation of partially or fully edentulous patients as dental prosthesis elements in combination with endosseous implants the maxilla and/or mandible.

2.7.3 Indications

DEDICAM Bridges, bars and implant base structures at abutment level are components compatible with dental implants for the restoration of partially or fully edentulous jaws with prosthetic restorations such as implant-supported bridges and full dentures. Specifically for:

- partially edentulous jaws with several missing teeth or

- edentulous jaws.

2.7.4 Contraindications

- Clinical situation which does not permit compliance with design specifications (see 2.8.5.1 “Design”)
- Reuse
- Use on an incompatible abutment
- Joint support on tooth (stump) and abutment
- Parafunctions, e.g. bruxism

Caution:

- The contraindications given in the respective manufacturer's instructions for the implant to be restored and the abutment used must be observed.
- The use of different alloy groups in the patient's mouth can lead to galvanic effects.

2.7.5 Application

2.7.5.1 Design

The following requirements must be met when designing bridges, bars and implant base structures at abutment level:

- Minimum wall thickness: 0.5mm
- Minimum connector cross sections:
 - Bridge element/splinted crowns, anterior tooth: 6mm²
 - Bridge element/splinted crowns, posterior tooth: 9mm²
 - Bars with an individualized profile and implant base structures: 8mm²
- Bridges made of titanium alloy: maximum of 3 connected pontics
- Bridges made of CoCr alloy: maximum of 3 connected pontics in the lateral region and 4 in the anterior region
- Bridges: maximum length of extensions; width of a molar, however, not more than 10mm
- Bars and implant base structures: maximum distance between two implant axes: 30mm
- Bars and implant base structures: maximum length of extensions: half the length of the adjacent inter-implant segment, however, not more than 10mm
- Indications and contraindications for the abutments according to the manufacturer's instructions are to be taken into consideration

The lower mechanical load-bearing capacity of diameter-reduced implants must be taken into account.

To facilitate the implementation of veneered bridges from the model in the sintering furnace, lateral retaining pins with a minimum diameter of 1mm and a maximum length of 5mm can be designed on the CAD model.

Various bar profiles and attachments are available as design options in the DEDICAM CAD library.

2.7.5.2 Post-processing

If no information on the surface conditioning of the bridge frameworks or implant base structures for the selected veneering/construction material is available from the respective manufacturer, bridge frameworks and implant base structures made of CoCr alloy are to be blasted with pure aluminum oxide with a grain diameter of minimum 110µm – 250µm and with a blasting pressure of 2 to 4bar, bridge frameworks and implant base structures made of titanium alloy with pure aluminum oxide with a grain diameter of 110µm – 250µm and a blasting pressure of 2 to 3 bar.

Warning: During machining or thermal processing of CoCr material (e.g. milling, grinding) dust and fumes can be produced which could entail a health risk if inhaled. Therefore, suitable measures must be applied to reduce dust and vapor exposure, such as the use of extraction and ventilation systems or, for short-term use, suitable personal protective equipment for respiratory air filtration.

The bridge frameworks and implant base structures must then be washed under running water, cleaned

with a steam jet and degreased with ethyl alcohol for titanium alloy never use for cleaning hydrofluoric acid). Do not contaminate the structure afterwards, e. g. touch it with unprotected fingers. Bridge frameworks are veneered with the layer material suitable for the material used (titanium or CoCr alloy) in each case or with CAD/CAM fabricated suprastructures. The manufacturer's processing instructions are to be observed here. In this context, please consider that discolorations of the titanium alloy can already occur at temperatures above 300°C.

Caution:

- The surface quality and wall thickness may not be changed or reduced in the load-bearing regions.
- Heating of titanium to above 800°C leads to structural changes with a reduction in strength. The processing instructions of the material manufacturers must be observed, these are available under:

CoCr alloy (Starbond CoS Disc Basic):	Please contact via email at service@scheftner.dental
Titanium alloy (Starbond Ti5 Disc):	Please contact via email at service@scheftner.dental

To avoid aesthetic disadvantages, CoCr restorations veneered with a layering material must have the oxide layer removed from the surfaces that come into contact with the soft tissue.

2.7.5.3 Handover to clinician

The post-processed structure must be completely cleaned of processing residues and other contamination using a suitable method.

The specific prosthetic screws supplied together with the bridges, bars and implant base structures at abutment level must be passed on to the recipient in the unopened original packaging and clearly assignable.

2.7.5.4 Insertion

Bridges, bars and implant base structures at abutment level are definitively screw-retained on the abutment after cleaning, disinfection and recommended sterilization (necessary if not set after completed healing phase, or if set subcrestal) as follows:

1. Remove protective caps or temporary prostheses.
2. Ensure that the contact surfaces from implant abutment to bridge, bar or implant base structure are absolutely clean and, in particular, that bone remains and soft tissue are completely removed.
3. Place bridge, bar or implant base structure on the definitively fixated abutment in the implant.
4. Screw in the new, unused prosthetic screw corresponding to the implant system and the implant diameter using the corresponding screwdriver of the implant manufacturer.
5. Tighten the prosthetic screws cross-wise with the specified torque of 15Ncm.

Caution:

- If the torque exceeds the recommended value, this may result in damage/destruction of the screw and/or damage to the abutment and/or the bridge, the bar or implant base structure.
 - If the torque lies below the recommended value, this can lead to loosening of the bridge, the bar or the implant base structure.
6. Check the proper seating on the abutment. No soft tissue may be trapped and there must be no visible gap between the bridge, bar or implant base structure and the abutment.
 7. Retighten with the same torque after at least 5 minutes.
 8. Bridges and implant base structures: cover the screw head with a completely removable material such as cotton or gutta-percha for protection. Then seal the screw channel tightly with suitable materials, e.g. a composite.
 9. Bars: for hygiene reasons, the screw channels are to be sealed with suitable and removable

materials. Beforehand, cover the screw head with a completely removable material such as cotton or gutta-percha.

2.8 Bridges, bars, and implant base structures for bonding bases (screw-retained)

2.8.1 Product description

DEDICAM Bridges for bonding bases are customized artificial restorative elements and serve as replacement for lost tooth substance and are screw-retained on two or more implants/abutments.

DEDICAM Bars for bonding bases serve as customized connection elements for removable prosthetic tooth restorations and are screw-retained on two or more implants/abutments.

DEDICAM Implant base structures for bonding bases serve as customized metallic support structures for prosthetic tooth restorations made of prosthetic plastic and prefabricated denture teeth made of plastic or ceramic and are screw-retained on two or more implants/abutments.

The structures for bonding bases are screw-retained occlusally on implants on prefabricated bonding bases offered specifically for this application, on CAMLOG/CONELOG Bar abutments or Multi-unit abutments. Bridges, bars and implant base structures made of titanium and CoCr alloy are only available for bonding bases from Camlog and BioHorizons.

The anatomical bridges and bridge frameworks for bonding bases can consist of 2 to 16 elements. Bars and implant base structures for bonding bases are available for 2 up to 10 implants. With the iSy Implant system, bridges and implant base structures must be supported on at least 6 implants, bars on at least 4 implants.

Bridges for bonding bases made of zirconia and plastic are available for individual characterization in different tooth shades. The zirconia "IPS e.max ZirCAD" is also available in various degrees of translucency and optionally with a natural color gradient from dentin to incisal. If hard and/or soft tissue is lost, bridges made of zirconia and plastic can be designed with a gingival structure component as an option.

Temporary restorations up to a maximum wearing time of 12 months can be made from a tooth-coloring plastic (polymethylmethacrylate).

Zirconia structures are already supplied densely sintered.

The anatomically reduced bridge framework can be veneered with suitable layering materials or with CAD/CAM fabricated suprastructures. Implant base structures serve as a metallic support structure for reinforcing restorations made of denture resin and for supporting prefabricated denture teeth.

With the aid of the DEDICAM CAD library, bridges and bars made of titanium or CoCr alloys can be supplemented with attachments for the fixation of removable restorations and bars can also be designed with different profiles. A suitable bar profile is available in the DEDICAM CAD library for implant base structures.

The structures are non-sterile and intended for single use.

Bridges and implant base structures have a milled surface, bars are also available polished.

Bridges, bars and implant base structures for CAMLOG, CONELOG, iSy or BioHorizons Bondingbases are delivered together with the titanium bases.

Suitable bar supraconstructions are available for DEDICAM Bars specifically designed for this purpose. The design of the bars and bar supraconstructions is performed by the DEDICAM Scan & Design Service.

2.8.1.1 CAMLOG, CONELOG and iSy Titanium bases CAD/CAM, bridge

Titanium bases CAD/CAM, bridge, with a considerably shortened implant-abutment connection and without antirotational mechanism are available for the CAMLOG, CONELOG and the iSy Implant

System. These are supplied complete with an abutment screw and are packaged separately with the structure but can also be reordered separately (see 3.12 "CAMLOG, CONELOG and iSy Titanium bases CAD/CAM, bridge").

2.8.1.2 *Titanium bonding bases and titanium caps for CAMLOG and CONELOG Bar abutments*

Titanium bonding bases for bar abutments and titanium caps without an antirotational mechanism are available for CAMLOG and CONELOG Bar abutments (see 3.14 "Titanium bonding "Passive-Fit" bases and titanium caps without retention for bridges for CAMLOG and CONELOG Bar abutments"). In this option, the bridge, bar or implant base structure is bonded intraorally to the titanium bonding bases or titanium caps screw-retained to the bar abutments. This allows compensating possible fit imprecisions (Passive Fit) which could lead to tension. The screw seat is on the titanium bonding base or the titanium cap.

A single bar abutment can be used as position reference for bridges, bars and implant bases structures made of titanium or CoCr alloy. To this purpose, the structure is screw-retained on this bar abutment without a bonding base. For all other bar abutments, the titanium bonding bases, or titanium caps must then be used to ensure freedom of tension.

In addition to the prosthetic screws of the structure, the corresponding "Passive Fit" titanium bonding bases are supplied separately packed but can also be reordered separately. Titanium caps without retention for bridge must be ordered separately (subject to change).

2.8.1.3 *BioHorizons Internal hybrid titanium bases, without hexagon*

Internal hybrid titanium bases are available for the BioHorizons Tapered internal and Internal implants, without hexagon, with greatly shortened implant-abutment connection and without antirotational mechanism. These are supplied complete with abutment screw and are packaged separately with the structure but can also be reordered separately (see 3.5 "BioHorizons Internal hybrid titanium bases, non-hexed").

2.8.1.4 *BioHorizons Bonding bases and titanium copings for Multi-unit abutments*

Bonding bases for the Passive Fit technique and titanium copings are available for Camlog and BioHorizons Multi-unit abutments (see [3.17 "BioHorizons Bonding bases and titanium copings for Multi-unit abutments"](#)). In this option, the bridge, bar or implant base structure is bonded intraorally to the bonding bases, Passive Fit screw-retained to the Multi-unit abutments. Alternatively, plastic bridges can be bonded to the titanium copings. This allows compensating possible fit imprecisions (Passive Fit) which could lead to tension. The screw seat is in the bonding base/titanium coping.

Material:

Titanium alloy (Ti-6Al-4V)		CTE (20 to 600°C): $\sim 10.3 \times 10^{-6}/K^{-1}$
CoCr alloy		CTE (20 to 500°C): $\sim 14.0 \times 10^{-6}/K^{-1}$
IPS e.max ZirCAD (zirconia)	MT and MT Multi	CTE (25 to 500°C): $10.4 \pm 0.5 \times 10^{-6}/K$
	LT	CTE (25 to 500°C): $10.5 \pm 0.5 \times 10^{-6}/K$
Telio CAD (polymethylmethacrylate)		CTE: not applicable

2.8.2 Intended purpose/intended use.

Bridges, bars and implant base structures are intended for the functional and/or esthetic oral rehabilitation of partially or fully edentulous patients as dental prosthesis elements in combination with endosseous implants in the maxilla and/or mandible.

2.8.3 Indications

DEDICAM Bridges, bars and implant base structures for bonding bases are components of a prosthetic solution for the restoration of partially or fully edentulous jaws with prosthetic restorations such as implant-supported bridges and full dentures. Specifically, for:

- partially edentulous jaws with several missing teeth or
- edentulous jaws.

2.8.4 Contraindications

- Clinical situation which does not permit compliance with design specifications (see 2.9.5.1 “Design”)
- Reuse
- Use with an incompatible bonding base.
- Joint support on tooth (stump) and implant or abutment
- Parafunctions, e.g. bruxism
- Restorations made of titanium alloy:
 - Bridge restorations with more than 3 connected pontics
- Restorations made of CoCr alloy:
 - Bridge restorations with more than 3 connected pontics in the posterior region and 4 in the anterior tooth region
- Restorations made of IPS e.max ZirCAD:
 - Bridge constructions with more than 2 connected pontics
 - For material types "MT" and "MT Multi": bridge constructions consisting of more than 3 elements.
- Restorations made of Telio CAD:
 - Wearing period over 12 months
 - Bridge constructions with more than 2 connected pontics

Caution:

- The contraindications of the implant to be restored and the bonding bases in the respective manufacturer's instructions must be observed. For the iSy Implant System, please note that bars are to be planned on at least 4 implants, bridges and implant base structures on at least 6 implants, without mesial or distal extensions and with an angulation correction of not more than 20° to the implant axis. Adequate polygonal support needs to be provided (cross arch).
- The use of different alloy groups in the patient's mouth can lead to galvanic effects.

2.8.5 Application

2.8.5.1 Design

The DEDICAM CAD material libraries and the “DEDICAM Parameter overview for CAD software” contain detailed indication and material-specific design data. Registered customers can obtain the library and the document from their local Camlog distributor.

The lower mechanical load-bearing capacity of diameter-reduced implants must be taken into account. To facilitate the implementation of veneered bridges from the model in the sintering furnace, lateral retaining pins with a minimum diameter of 1mm and a maximum length of 5mm can be designed on the CAD model. The following requirements must be met when designing bridges, bars and implant base structures for bonding bases:

- Minimum wall thickness: 0.5mm
- Minimum connector cross sections:
 - Bridge element/splinted crowns, anterior tooth:

▪ Titanium alloy and CoCr:	6mm ²
▪ IPS e.max ZirCAD LT:	9mm ²
▪ IPS e.max ZirCAD MT/MT Multi:	12mm ²
 - Bridge element/splinted crowns, posterior tooth:

▪ Titanium alloy and CoCr:	9mm ²
▪ IPS e.max ZirCAD LT:	12mm ²
▪ IPS e.max ZirCAD MT/MT Multi:	16mm ²
 - Bars with an individualized profile and implant base structures: 8mm²
- Bridges made of titanium alloy: maximum of 3 connected pontics; in the case of sole support on 2 implants, the maximum distance between the implant axes is 30mm
- Bridges made of CoCr alloy: maximum of 3 connected pontics in the lateral region and 4 in the anterior region; with sole support on 2 implants, the maximum distance between the implant

axes is 30mm

- Bridges made of zirconia or plastic: maximum of 2 connected pontics (IPS e.max ZirCAD MT/MTMulti: 1 pontic)
- Bridges: maximum length of extensions; width of a molar, however, not more than 10mm
- Bars and implant base structures: maximum distance between two implant axes: 30mm
- Bars and implant base structures: maximum length of extensions: half the length of the adjacent inter-implant segment, however, not more than 10mm
- Indications and contraindications for the bonding bases according to the manufacturer's instructions are to be taken into consideration
- With the iSy Implant system, bridges and implant base structures must be supported on at least 6 implants, bars on at least 4 implants

The lower mechanical load-bearing capacity of diameter-reduced implants must be taken into account. To facilitate the implementation of veneered bridges from the model in the sintering furnace, lateral retaining pins with a minimum diameter of 1mm and a maximum length of 5mm can be designed on the CAD model.

Various bar profiles and attachments are available as design options in the DEDICAM CAD library.

2.8.5.2 Post-processing

If no information on the surface conditioning of the bridge frameworks or implant base structures for the selected veneering/construction material is available from the respective manufacturer, bridge frameworks and implant base structures made of CoCr alloy are to be blasted with pure aluminum oxide with a grain diameter of minimum 110µm and with a blasting pressure of 3 to 4bar, bridge frameworks and implant base structures made of titanium alloy with pure aluminum oxide with a grain diameter of 180µm and a blasting pressure of 2 to 3bar.

Warning: During machining or thermal processing of CoCr material (e.g. milling, grinding) dust and fumes can be produced which could entail a health risk if inhaled. Therefore, suitable measures must be applied to reduce dust and vapor exposure, such as the use of extraction and ventilation systems or, for short-term use, suitable personal protective equipment for respiratory air filtration.

The bridge frameworks or implant base structures must then be washed under running water, cleaned with a steam jet and degreased with ethyl alcohol. Do not contaminate the structure afterwards, e. g. touch it with unprotected fingers.

Bridge frameworks are veneered with the layer material suitable for the material used (titanium alloy, CoCr alloy or IPS e.max ZirCAD) in each case or with CAD/CAM fabricated suprastructures. The manufacturer's processing instructions are to be observed here. In this context, please consider that discolorations of the titanium alloy can already occur at temperatures above 300°C.

Bridges made of Telio CAD can be modified to optimize the emergence profile by applying suitable composite or – taking into account the minimum permissible wall thicknesses – by grinding.

Caution:

- The surface quality and wall thickness may not be changed or reduced in the load-bearing regions.
- Heating of titanium to above 800 °C leads to structural changes with a reduction in strength.

The processing instructions of the material manufacturers must be observed, these are available under:

CoCr alloy (Starbond CoS Disc Basic):	Please contact via email at service@scheftner.dental
Titanium alloy (Starbond Ti5 Disc):	Please contact via email at service@scheftner.dental
IPS e.max ZirCAD (zirconia):	www.ivoclar.com/en_gb/eifu
IPS e.max CAD (lithium disilicate):	www.ivoclar.com/en_gb/eifu

To avoid aesthetic disadvantages, CoCr restorations veneered with a layering material must have the oxide layer removed from the surfaces that come into contact with the soft tissue.

2.8.5.3 Bonding on the model (not for bonding bases for CAMLOG or CONELOG Bar abutments or Camlog and BioHorizons Multi-unit abutments)

Before bonding with the CAD/CAM structure, the adhesive surface of the bonding bases must be conditioned according to the specifications of the material and fixation composite manufacturers. The selected and prepared bonding bases, which have a common insertion axis, are screw-retained to the implant analogs in the model with the bonding aid. Then the bridge, bar or implant base structure are bonded extra orally to the screw-retained bonding bases. If the titanium bases have different insertion axes, each joint insertion axis must be bonded to the structure separately. Due to the flexibility of the bonding aids, these can be inserted into diverging screw channels by bending.

Suitable surface conditioners and fixation composites

To connect the bridge, bar or implant base structure to the bonding bases, we recommend the following products from Ivoclar Vivadent:

- Bonding agent for PMMA abutment: "SR Connect"
- Primer for bonding base and zirconia abutment: "Monobond Plus"
- Fixation composites: "Multilink Hybrid Abutment"

Observe the manufacturer's processing instructions. Biotech Dental Digital SAS accepts no responsibility for the durability of the bond.

2.8.5.4 Handover to clinician

The post-processed structure must be completely cleaned of processing residues and other contamination using a suitable method.

The following are to be forwarded to the recipient together with the bridge, bar or the implant base structure in the unopened, original primary packaging and clearly assigned:

- Structures with bonded CAMLOG, CONELOG or iSy Titanium bases CAD/CAM, bridge or BioHorizons Internal Hybrid Abutment Base, non-hexed: the specific abutment screws supplied by Camlog.
- Structures for titanium bonding bases or titanium caps for CAMLOG or CONELOG Bar abutments or BioHorizons Bonding bases, Passive Fit/titanium caps for Multi-unit abutments: the specific bonding bases and prosthetic screws supplied by Camlog

2.8.5.5 Insertion

Before the final fixation of the structure, the manufacturer's instructions in the user information for the bonding bases and the fixation composite must be observed for intraoral bonding of the titanium bonding bases or titanium caps in bridges, bars or implant base structures.

Suitable surface conditioners and fixation composites

To connect the bridge, bar or implant base structure to the bonding bases, we recommend the following products from Ivoclar Vivadent:

- Bonding agent for PMMA abutment "SR Connect"
- Primer for zirconia structures: "Monobond Plus"
- Fixation composite for structures made of zirconia, CoCr and titanium alloy: "Multilink Automix"
- Fixation composite for structures made of PMMA: "Multilink Automix" or (not sterilizable) „Telio CS Cem Implant“

Observe the manufacturer's processing instructions. Biotech Dental Digital SAS accepts no responsibility for the durability of the bond.

Bridges, bars and implant base structures for bonding bases are definitively screw-retained on the implants, bar abutments or Multi-unit abutments after cleaning, disinfection and recommended sterilization (necessary if not set after completed healing phase, or if set subcrestal) as follows:

1. At implant level: remove implant cover screw or implant abutment and clean the inner configuration and shoulder of the implant.
2. At abutment level: remove protective caps or temporary prostheses.

3. Ensure that the contact surfaces from implant or bar abutment or Multi-unit abutment to the bridge, bar or implant base structure are absolutely clean and, in particular, that bone remains and soft tissue are completely removed.
4. Place bridge, bar or implant base structure on the implants or definitively fixated bar abutments or Multi-unit abutments on the implants.
5. Screw in the new, unused abutment or prosthetic screw using the corresponding screwdriver.
6. Tighten the screws cross-wise with the specified torque.

Caution:

- If the torque exceeds the recommended value, this may result in damage/destruction of the screw and/or damage to the implant or abutment and/or the bridge, the bar or implant base structure.
 - If the torque lies below the recommended value, this can lead to loosening of the bridge, the bar or the implant base structure.
7. Check for proper seating on the implants or bar abutments or Multi-unit abutments. No soft tissue may be trapped and there must be no visible gap between the bridge, bar or implant base structure and the implant or bar abutment or Multi-unit abutment.
 8. Retighten with the same torque after at least 5 minutes.
 9. Bridges and implant base structures: cover the screw head with a completely removable material such as cotton or gutta-percha for protection. Then seal the screw channel tightly with suitable materials, e.g. a composite.
 10. Bars: for hygiene reasons, the screw channels are to be sealed with suitable and removable materials. Beforehand, cover the screw head with a completely removable material such as cotton or gutta-percha.

2.9 Suprastructures for bars

2.9.1 Product description

DEDICAM Suprastructures for bars are customized support elements for removable prostheses that are supported directly or on holding elements on a specially designed DEDICAM Bar.

The suprastructures for bars and corresponding bars are designed virtually and fabricated customized in the DEDICAM Manufacturing Center.

Optionally, the bar can be tried in on the patient before the suprastructure for bars is constructed.

The milled structures have a congruent inner geometry to the bar and are finished with prefabricated teeth and denture resin in the dental laboratory to create a dental prosthesis.

These attachments as retaining elements for fixing the suprastructure for bars on the bar are available:

- Preci-Vertex/Preci-Horix
- Mount for MK1
- Locator abutment with M2.0 thread

Male parts for Preci matrices are available in different lengths for the respective anatomical situation. If the matrices (female parts) are to be shortened, it must be ensured that they only rest on the parallel surface of the male part so that a complete overlap is possible.

Retention pins can be attached for the fixation of prefabricated denture teeth. Retention grids are used for the construction of free-end saddles in the upper jaw, and in the lower jaw, hole retentions are used in the covering with denture resin.

The structures are non-sterile, intended for single use and available in both titanium and CoCr alloy materials with finely milled or grained surfaces.

Material:

Titanium alloy (Ti-6Al-4V)
CoCr alloy

CTE (20 to 600°C): $\sim 10.3 \times 10^{-6}/K^{-1}$
CTE (20 to 500°C): $\sim 14.0 \times 10^{-6}/K^{-1}$

2.9.2 Intended purpose/intended use

Suprastructures for bars are intended for the functional and/or esthetic oral rehabilitation of partially or fully edentulous patients as dental prosthesis elements in combination with endosseous implants in the maxilla and/or mandible.

2.9.3 Indications

DEDICAM Suprastructures for bars are components of a prosthetic solution for the restoration of partially or fully edentulous jaws with prosthetic restorations such as implant-supported bridges and full dentures. Specifically for:

- partially edentulous jaws with several missing teeth or
- edentulous jaws.

2.9.4 Contraindications

- Reuse
- Parafunctions, e.g. bruxism

Caution: The use of different alloy groups in the patient's mouth can lead to galvanic effects.

2.9.5 Application

2.9.5.1 Design

The design is performed exclusively by specifically trained employees in the DEDICAM Scan and Design Service.

2.9.5.2 Further processing

The bar supraconstructions are finished with prefabricated teeth and denture resin in the dental laboratory to create a dental prosthesis according to conventional methods. For this purpose, the processing instructions for the respective materials are to be observed. These are available under:

CoCr alloy (Starbond CoS Disc Basic):

Please contact via email at service@scheftner.dental

Titanium alloy (Starbond Ti5 Disc):

Please contact via email at service@scheftner.dental

Warning: During machining or thermal processing of CoCr material (e.g. milling, grinding) dust and fumes can be produced which could entail a health risk if inhaled. Therefore, suitable measures must be applied to reduce dust and vapor exposure, such as the use of extraction and ventilation systems or, for short-term use, suitable personal protective equipment for respiratory air filtration.

2.9.5.3 Handover to clinician

The bar supraconstruction incorporated in the prosthesis must be completely cleaned of processing residues and other contamination using a suitable method.

2.10 Crowns and bridges (bonded)

2.10.1 Product description

DEDICAM Crowns, bridges, crown and bridge frameworks are bonded to tooth stumps or customized abutments. They are designed virtually with the aid of a suitable CAD software and – if compatible with the software used – the DEDICAM CAD material library.

Depending on the selected material, the bridges or bridge frameworks can consist of 2 to 16 elements and can be constructed with a maximum of 3 connected pontics (4 in the anterior region when using CoCr

alloy, 2 when using zirconia, 1 when using lithium disilicate or IPS e.max ZirCAD MT/MT Multi). Appropriate connection cross sections are to be taken into consideration.

DEDICAM Crowns, bridges, crown and bridge frameworks are available for selected materials for individual characterization in various tooth shades and in various degrees of translucency, depending on the indication and material. Zirconia (IPS e.max ZirCAD MT Multi) are also available with a natural shade gradient from dentin to the incisal edge.

Crowns, bridges, crown and bridge frameworks made of zirconia are supplied already densely sintered and in the “blue” metasilicate phase for lithium disilicate. Lithium disilicate structures are to be subjected to crystallization firing.

The anatomically reduced crown or bridge framework can be veneered with suitable layering materials or with CAD/CAM fabricated suprastructures. If hard and/or soft tissue is lost, bridges and bridge frameworks can be designed with a gingival structure component as an option.

Material (Telio CAD are not suitable for frameworks):

Titanium alloy (Ti-6Al-4V)		CTE (20 to 600°C): $\sim 10.3 \times 10^{-6}/K^{-1}$
CoCr alloy		CTE (20 to 500°C): $\sim 14.0 \times 10^{-6}/K^{-1}$
IPS e.max ZirCAD (zirconia)	MT and MT Multi	CTE (25 to 500°C): $10.4 \pm 0.5 \times 10^{-6}/K$
	LT	CTE (25 to 500°C): $10.5 \pm 0.5 \times 10^{-6}/K$
IPS e.max CAD (lithium disilicate)		CTE (100 to 400°C): $\sim 10.2 \times 10^{-6}/K$ CTE (100 to 500°C): $\sim 10.5 \times 10^{-6}/K$
Telio CAD (polymethylmethacrylate)		CTE: not applicable

2.10.2 Intended purpose/intended use

Crowns, bridges, crown and bridge frameworks, bonded to tooth stumps or customized abutments are intended for the functional and/or esthetic oral rehabilitation of partially or fully edentulous patients or patients with compromised teeth as dental prosthesis elements in the maxilla and/or mandible.

2.10.3 Indications

DEDICAM Crowns, bridges, crown and bridge frameworks, bonded to tooth stumps or customized abutments are components of a prosthetic solution for the restoration of partially or fully edentulous jaws or tooth defects with prosthetic restorations.

Restorations made of IPS e.max ZirCAD have different indications depending on their translucency, as shown in the following diagram:

Translucency level	Indications				
	Full-contour crowns	Full-contour 3-unit bridges	Full-contour, 4- to multi-unit bridges with max. 2 pontics	Crown frameworks	3- to multi-unit bridge frameworks with max. 2 pontics
MT Multi Medium Translucency with shade gradation	✓	✓			
MT Medium Translucency	✓	✓			
LT Low Translucency	✓	✓	✓	✓	✓

Figure 9: Indication overview IPS e.max ZirCAD (Reference: Ivoclar Vivadent AG, FL-Schaan)
Version 04

2.10.4 Contraindications

- Clinical situation which does not permit compliance with design specifications (see 2.11.5.1 “Design”)
 - Reuse
 - Abutment-supported restorations: Parafunctions, e.g. bruxism
 - Restorations made of titanium alloy:
 - Bridge restorations with more than 3 connected pontics
 - Restorations made of CoCr alloy:
 - Bridge restorations with more than 3 connected pontics in the posterior region and 4 in the anterior tooth region
 - Restorations made of IPS e.max CAD ZirCAD, LT:
 - Bruxism, in case of veneered restorations
 - Bridge constructions with more than 2 connected pontics or 2 or more connected extension elements
 - Temporary insertion
 - Restorations made of IPS e.max ZirCAD, MT and MT Multi:
 - Bruxism
 - Bridge restorations with more than 1 intermediate pontic or more than 3 elements
 - Temporary insertion
 - Restorations made of IPS e.max CAD:
 - Very deep subgingival preparations
 - Pronounced reduced residual dentition
 - Parafunctions, e.g. bruxism
 - Bridges with more than 3 elements
 - Bridges with a more than 11mm (anterior tooth) or 9mm (premolar) wide pontic
 - Cross-section of connectors smaller than 16mm²
 - Bridges with screw hole
 - Bridge restorations, distal to the 2nd premolar
 - Temporary insertion
 - Restorations made of Telio CAD:
 - Wearing period over 12 months
 - Bridge constructions with more than 2 connected pontics
 - Parafunctions, such as bruxism

Caution:

- For abutment-supported restorations: the contraindications given in the respective manufacturer's instructions for the implant to be restored and the abutment used must be observed.
- The use of different alloy groups in the patient's mouth can lead to galvanic effects.

2.10.5 Application

2.10.5.1 Design

The DEDICAM CAD material libraries and the “DEDICAM Parameter overview for CAD software” contain detailed indication and material-specific design data. Registered customers can obtain the library and the document from their local Camlog distributor.

In particular, the specified minimum wall thickness must be considered when designing. Any change in the cement thickness specified in the material libraries or in the parameter chart may result in a misfit of the restoration.

To obtain perfect results for tooth restorations, the general and material-specific preparation guidelines are to be observed:

- For restorations on tooth stumps: remove sufficient hard tooth substance
- Blocking out undercuts
- Preparation of grooves
- Avoid sharp transitions

- Observe material-specific transitions

To facilitate the implementation of veneered bridges and crowns from the model in the sintering furnace, lateral retaining pins with a minimum diameter of 1mm and a maximum length of 5mm can be designed on the CAD model.

2.10.5.2 Post-processing

If no information on the surface conditioning of the bridge and crown frameworks for the selected veneering material is available from the respective manufacturer, frameworks made of CoCr alloy are to be blasted with pure aluminum oxide with a grain diameter of minimum 110 µm and a blasting pressure of 3 to 4 bar, frameworks made of titanium alloy with pure aluminum oxide with a grain diameter of 180 µm and a blasting pressure of 2 to 3 bar.

Warning: During machining or thermal processing of CoCr material (e.g. milling, grinding) dust and fumes can be produced which could entail a health risk if inhaled. Therefore, suitable measures must be applied to reduce dust and vapor exposure, such as the use of extraction and ventilation systems or, for short-term use, suitable personal protective equipment for respiratory air filtration.

The frameworks must then be washed under running water, cleaned with a steam jet and degreased with ethyl alcohol. Do not contaminate the structure afterwards, e. g. touch it with unprotected fingers.

The frameworks are veneered with the layer material suitable for the material used (titanium, CoCr alloy or IPS e.max ZirCAD, MT and MT Multi) in each case or with CAD/CAM fabricated suprastructures. The manufacturer's processing instructions are to be observed here. In this context, please consider that discolorations of the titanium alloy can already occur at temperatures above 300°C.

Caution:

- The surface quality and wall thickness may not be changed or reduced in the load-bearing regions.
- Heating of titanium to above 800°C leads to structural changes with a reduction in strength.

The processing instructions of the manufacturers for the following materials must be observed, these are available under:

CoCr alloy (Starbond CoS Disc Basic):	Please contact via email at service@scheftner.dental
Titanium alloy (Starbond Ti5 Disc):	Please contact via email at service@scheftner.dental
IPS e.max ZirCAD (zirconia):	www.ivoclar.com/en_us/eifu
IPS e.max CAD (lithium disilicate):	www.ivoclar.com/en_us/eifu
Telio CAD (polymethylmethacrylate):	www.ivoclar.com/en_us/eifu

To avoid aesthetic disadvantages, CoCr restorations veneered with a layering material must have the oxide layer removed from the surfaces that come into contact with the soft tissue.

2.10.5.3 Handover to clinician

The post-processed structures must be completely cleaned of processing residues and other contamination using a suitable method.

2.11 Double crowns

2.11.1 Product description

DEDICAM Primary components (telescopes) for double crown restorations serve as customized supporting elements for a removable buttress-congruent, self-adhesive, removable prosthetic tooth restorations in the jaw.

They are bonded to individual tooth stumps or abutments. Self-locking of the double crowns in their final position, due to static friction, is achieved by fitting the primary and secondary parts in the laboratory, including necessary polishing of the contact surfaces. The degree of self-locking can be determined by the fit, size, angle and roughness of the contact surfaces.

The sliding friction to be achieved with parallel-walled contact surfaces when the double crowns are superimposed is defined by the fit dimensions, size and roughness of the surfaces.

Material:

Titanium alloy (Ti-6Al-4V)	CTE (20 to 600°C): $\sim 10.3 \times 10^{-6}/K^{-1}$
CoCr alloy	CTE (20 to 500°C): $\sim 14.0 \times 10^{-6}/K^{-1}$

2.11.2 Intended purpose/intended use

Crowns, bridges, crown and bridge frameworks are intended for the functional and/or esthetic oral rehabilitation of partially or fully edentulous patients or patients with compromised teeth as dental prosthesis elements in the maxilla and/or mandible.

2.11.3 Indications

DEDICAM Primary components (telescopes) for double crown restorations are components of a prosthetic solution for the restoration of partially or fully edentulous jaws with prosthetic restorations.

2.11.4 Contraindications

- Clinical situation which does not permit compliance with design specifications (see 2.12.5.1 “Design”)
- Reuse
- Primary splints
- Abutment-supported restorations: Parafunctions, e.g. bruxism

Caution:

- For abutment-supported restorations: the contraindications given in the respective manufacturer's instructions for the implant to be restored and the abutment used must be observed.
- The use of different alloy groups in the patient's mouth can lead to galvanic effects.

2.11.5 Application

2.11.5.1 Design

After defining the common insertion axis, the design of the primary components in the CAD software is performed after selection of the corresponding indication, taking into account the current DEDICAM design instructions for CAD software (registered customers receive these from their responsible DEDICAM contact person) and the following factors:

- The best possible compromise between anatomical tooth shape and – for a non-tilting and abrasion-resistant fit – sufficient contact surface area
- Minimum wall thickness: 0.5mm
- Minimum cervical shoulder radius: 0.6mm
- Minimum upper shoulder radius: 0.8mm
- Avoiding undercuts at the contact surface areas
- Parallelism of all contact surface areas in accordance with the insertion axis
- As equal as possible length of the contact surface areas, approx. 5mm is recommended
- Chamfering or rounding of the incisal/occlusal geometry to the contact surface (for easier insertion of the secondary component)

Additional retention elements cannot be added.

2.11.5.2 Parallel milling of the primary component

Warning: During machining or thermal processing of CoCr material (e.g. milling, grinding) dust and fumes can be produced which could entail a health risk if inhaled. Therefore, suitable measures must be applied to reduce dust and vapor exposure, such as the use of extraction and ventilation systems or, for short-term use, suitable personal protective equipment for respiratory air filtration.

Note: A trial-in of the primary component and an impression of the situation with a suitable impression material and subsequent creation of a master cast prior to post-processing has proved successful.

The DEDICAM Primary component is attached to the milling model, the contact surfaces are aligned in the insertion axis and milled in parallel in the milling machine. The milling cutter used should have a radius of at least 0.6mm.

Then the contact surfaces are polished.

The secondary components (not part of the DEDICAM portfolio) are bonded to the removable denture in their final position and by placing them onto the already fixed primary components.

2.11.5.3 Handover to clinician

The post-processed double crowns, together with the restoration, must be completely cleaned of processing residues and other contamination using a suitable method.

2.12 Veneers, inlays, onlays and partial crowns

2.12.1 Product description

DEDICAM Veneers, inlays, onlays and partial crowns are bonded to teeth and serve to replace tooth substance. They are available for the individual characterization in various tooth shades and in two degrees of translucency. Lithium disilicate structures are supplied in the “blue” metasilicate phase and are to be subjected to crystallization firing.

Material:

IPS e.max CAD (lithium disilicate)*

CTE (100 to 400°C): $\sim 10.2 \times 10^{-6}/K$

CTE (100 to 500°C): $\sim 10.5 \times 10^{-6}/K$

2.12.2 Intended purpose/intended use

Veneers, inlays, onlays and partial crowns are intended for the functional and/or esthetic oral rehabilitation as restorative elements of compromised teeth in the maxilla and/or mandible.

2.12.3 Indication

DEDICAM Veneers, inlays, onlays and partial crowns are components for the restoration of compromised teeth with missing tooth substance in the maxilla and/or mandible.

2.12.4 Contraindications

- Reuse
- Restorations made of IPS e.max CAD:
 - Very deep subgingival preparations
 - Pronounced reduced residual dentition
 - Parafunctions, e.g. bruxism
 - Temporary insertion or trial wearing

2.12.5 Application

2.12.5.1 Design

To obtain perfect results for tooth restorations, the general and material-specific preparation guidelines are to be observed:

- Removal of sufficient hard tooth substance
- Blocking out undercuts
- Preparation of grooves
- Avoid sharp transitions
- Observe material-specific transitions

The DEDICAM CAD material libraries and the “DEDICAM Parameter overview for CAD software” contain detailed indication and material-specific design data. Registered customers can obtain the library and the

document from their local Camlog distributor.

2.12.5.2 Post-processing

The processing instructions of the material manufacturer must be observed, these are available under:

IPS e.max CAD (lithium disilicate): www.ivoclar.com/en_gb/eifu

2.12.5.3 Handover to clinician

The post-processed structures must be completely cleaned of processing residues and other contamination using a suitable method.

3 Compatible products for DEDICAM Implant prosthetics

Caution: The following information is a brief description of the products and does not replace the corresponding information for use. These depend on the product and are enclosed with the product or can be viewed or obtained for Camlog products under <https://ifu.camlog.com>.

3.1 Scanbodies for Camlog implant systems

Various scanbodies are available for optical, three-dimensional localization (including determination of the axial inclination and rotational alignment) of Camlog implants or bar abutments in the mouth and the Camlog implant analogs or bar abutment laboratory analogs in the working model.

The geometries of these scanbodies are stored in the DEDICAM CAD libraries. The user must correctly allocate these geometries to the scanbodies used in the mouth or on the model. In other words, it is essential to correctly “match” the two datasets from the scan and the CAD library.

3.1.1 CAMLOG, CONELOG and iSy Scanbodies (at implant level)

The scanbodies specific to the implant system rest on the implant shoulder, thus eliminating the height variance of conical implant-abutment connections inherent to the system.

CAMLOG Scanbodies are available for all CAMLOG Implants with the reduced diameter contact surface to the implant and are therefore also suitable for the Platform Switching option (not for implants with a diameter of 3.3mm).

The scanbodies are available in either PEEK or Titanium alloy.

The ones made of PEEK are for single intraoral use and supplied sterile, each with an sterile abutment screw. For scanning the working model, Camlog recommends the use of the PEEK scanbody, which can be used repeatedly taking its integrity into account.

The ones made of Titanium alloy are for multiple intraoral use and supplied non-sterile, each with a non-sterile abutment screw.

The screw is hand-tightened using the CAMLOG or iSy Screwdriver.

CAMLOG and CONELOG Scanbodies are not suitable for digital acquisition in the design of bridges and bars on CAMLOG and CONELOG Bar abutments. For this purpose, scanning caps are available for bar abutments (see 3.1.3. “Scanning caps for bar abutments (at abutment level)”).

3.1.2 iSy Multifunctional cap (on iSy Implant base)

As an alternative to the iSy Scanbody, iSy Implants can also be localized via the iSy Multifunctional cap. The multifunctional cap is mounted directly on the iSy Implant base, is made of PEEK, is sterile, and included with every iSy Implant and is also available separately.

Because of the high accuracy required, using the iSy Multifunctional cap when scanning for the design of:

- double crowns and for the simultaneous design of abutments and framework or crown (“file-splitting”) is not possible due to the high level of precision required.
- multi-unit structures at implant level are not recommended.

3.1.3 Scanning caps for bar abutments (at abutment level)

Scanning caps for bar abutments are available for the scan required to design bridges and bars on CAMLOG and CONELOG Bar abutments at abutment level.

The scanning caps are available in either PEEK or Titanium alloy and for the two prosthetic platform diameters of 4.3mm (suitable for implants with a diameter of 3.3, 3.8 and 4.3mm) and 6.0mm (suitable for implants with a diameter of 5.0 and 6.0mm).

The ones made of PEEK are for single intraoral use and supplied sterile, each with an sterile prosthetic screw (light blue anodized). For scanning the working model, Camlog recommends the use of the PEEK scanbody, which can be used repeatedly taking its integrity into account.

The ones made of Titanium alloy are for multiple intraoral use and supplied non-sterile, each with an non-sterile prosthetic screw (light blue anodized).

3.1.4 Scanning caps for iSy Multi-unit abutments (at abutment level)

Scanning caps for Multi-unit abutments are available for the scan required to design bridges and bars on iSy Multi-unit abutments at abutment level.

The scanning caps are available in either PEEK or Titanium alloy.

The ones made of PEEK are for single intraoral use and supplied sterile, each with an sterile prosthetic screw (light blue anodized). For scanning the working model, Camlog recommends the use of the PEEK scanbody, which can be used repeatedly taking its integrity into account.

The ones made of Titanium alloy are for multiple intraoral use and supplied non-sterile, each with an non-sterile prosthetic screw (light blue anodized).

3.1.5 Scanbodies for BioHorizons Implants

Various "Internal scanbodies with snap-in function" are available for optical, three-dimensional localization (including determination of the axial inclination and rotational alignment) of BioHorizons Implants with the Internal Hex implant-abutment connection in the mouth and the implant analogs in the working model.

The scanbodies are made of PEEK with a sleeve made of titanium alloy and are supplied non-sterile.

The geometries of these scanbodies are stored in the DEDICAM CAD libraries. The user must correctly allocate these geometries to the scanbodies used in the mouth or on the model. In other words, it is essential to correctly "match" the two datasets from the scan and the CAD library.

3.2 Scanbodies for other implant systems and Multi-unit abutments

DEDICAM Scanbodies are available for the optical, three-dimensional localization (incl. determination of the axial inclination and rotational alignment) of implants from other implant suppliers and of Multi-unit abutments from Camlog or BioHorizons in the mouth as well as the corresponding implant or abutment analogs in the working model.

The scanbodies are made of PEEK with a titanium alloy insert and are supplied non-sterile, each with a pre-assembled titanium alloy fixing screw. The fixing screw is secured against falling out in the scanbody.

Undamaged scanbodies can be used several times, taking into account their integrity.

The fixing screw is activated with the DEDICAM Screwdrivers for further implant systems. The screwdrivers are available in three different lengths.

Only DEDICAM Scanbodies with the additional designation "IO" (article numbers D3013.xxxx or D0068.6417 for the ICX implant system) may be used intraorally.

The geometries of the DEDICAM Scanbodies are stored in the DEDICAM CAD libraries. The user must correctly allocate these geometries to the scanbody used on the model. In other words, it is essential to correctly "match" the two datasets from the scan and the CAD library.

3.3 Abutment screws

For the fixation of DEDICAM Abutments, DEDICAM One-piece healing caps and DEDICAM Bridges, bars and implant base structures at implant level, specific abutment screws for the respective implant system are already supplied with the structure but can be reordered separately.

The following tightening torques of the abutment screws must be compared with the respective information in the implant manufacturer's user information and, if necessary, reduced (not valid for healing caps. In this case, the screws must only be tightened by hand):

Implant provider	Implant system	Platform or implant Ø [mm]	Abutment screw, thread	Abutment screw, tightening torque ¹⁾ [Ncm]
Camlog	CAMLOG	3.3	M1.6	20
		3.8		
		4.3		
		5.0	M2.0	
		6.0		
	CONELOG	3.3	M1.6	
		3.8		
		4.3		
		5.0	M2.0	
	iSy	3.8	M1.6	
		4.4		
		5.0		
BioHorizons	Internal Hex	3.0	M1.8 ²⁾	30
		3.5		
		4.5		
		5.7		
		5.7		

Table 4: Assignment and tightening torque for abutment screws

- 1) Not applicable to healing caps. In this case, the screws must only be tightened by hand.
- 2) Different screws must be used for one-piece and multi-unit structures. A designated screw must be used for angulated screw channels in mesostructures and crowns on bonding bases.

The abutment screws for implant systems from other suppliers are compatible with the screwdrivers of the respective implant manufacturer. The abutment screws are made of titanium alloy and are available non-sterile and individually packaged. Depending on the implant system and platform, some abutment screws have a carbon (DLC) or titanium nitride (TiN) coating.

3.4 Prosthetic screws for CAMLOG and CONELOG Bar abutments

For the fixation of DEDICAM Bridges, bars and implant base structures as well as the scanning cap for bar abutments on CAMLOG and CONELOG Bar abutments, light blue anodized prosthetic screws are available in two thread sizes (see following table).

Implant system	Implant Ø [mm]	Screw thread	Prosthetic screw, torque [Ncm]
CAMLOG	3.3	M1.6	15 (Scanning cap, tightened by hand)
	3.8		
	4.3		
	5.0	M2.0	15 (Scanning cap, tightened by hand)
	6.0		
CONELOG	3.3	M1.6	15 (Scanning cap, tightened by hand)
	3.8		
	4.3		
	5.0	M1.6 (titanium-colored, for bar abutments with article number C2255.50xx)	15 (Scanning cap, tightened by hand)
		M2.0 (for bar abutments with article number C2254.50xx, C2256.50xx, C2257.50xx, C2258.50xx, C2259.50xx)	15 (Scanning cap, tightened by hand)

Table 5: Assignment and tightening torque of prosthetic screws for CAMLOG and CONELOG Bar abutments

The prosthetic screws are made of titanium alloy and are available non-sterile and individually packaged. Activation is via the CAMLOG Screwdrivers, hex.

3.5 Prosthetic screw for Camlog and BioHorizons Multi-unit abutments

For the fixation of DEDICAM bridges, bars and implant base structures on Camlog or BioHorizons Multi-unit abutments, the light blue anodized "Prosthetic screw, multi-unit, regular" is available in three pack sizes of one, five and 25 pieces, independent of the implant diameter. The screw is non-sterile, made of titanium alloy and is tightened with 15Ncm using the "Hex insertion tool" or "Unigrip™ Screwdriver".

Each package of the Bonding base with plastic sleeve, Passive Fit, for Multi-unit abutment as well as the Titanium coping, for Multi-unit abutment includes this screw.

3.6 Prosthetic screws for Nobel Biocare Multi-unit abutments

DEDICAM Prosthetic screws are available in two thread sizes for the fixation of DEDICAM Bridges, bars and implant base structures to Nobel Biocare Multi-unit abutments. The following tightening torques for the prosthetic screws must be compared with the information provided by the abutment manufacturer and reduced, if necessary:

Platform	Thread	Tightening torque [Ncm]
Narrow Platform (NP) 3.5mm	M1.4	15
Regular Platform (RP) 4.3mm		
Wide Platform (WP) 5.0mm	M1.8	15

Table 6: Assignment and tightening torque of prosthetic screws for Multi-unit abutments

3.7 Lab screws for Camlog implant systems

System-related lab screws are available for the fixation of screw-retained DEDICAM Prosthetics and DEDICAM Healing caps to the CAMLOG, CONELOG or iSy Implant analogs or the CAMLOG and CONELOG Bar abutment laboratory analogs in the working model. To avoid confusion with the definitive abutment or prosthetic screws, the lab screws are colored differently (see following table).

Implant system	Implant Ø [mm]	Lab screw, thread	Color code	Tightening torque
CAMLOG	3.3	M1.6	brown	tightened by hand
	3.8			
	4.3			
	5.0	M2.0	brown	
	6.0			
CONELOG	3.3	M1.6	brown	
	3.8			
	4.3			
	5.0	M2.0 (M1.6 for bar abutments with article number C2255.50xx)	brown	

iSy	3.8	M1.6	yellow-green	
	4.4			
	5.0			

Table 7: Assignment and tightening torque lab screws

The lab screws are made of titanium alloy, and are available non-sterile and individually packaged.

3.8 Screwdriver for abutment, prosthetic and lab screws for Camlog implant systems

System-related screwdrivers are available for the fixation of screw-retained DEDICAM Prosthetics and DEDICAM Healing caps on Camlog implants with the respective CAMLOG, CONELOG or iSy Lab or abutment screws or the CAMLOG and CONELOG Prosthetic screws.

The screwdrivers are made of stainless steel, and are available non-sterile and for different clinical situations in different lengths.

For use with angulated screw channels designated screwdrivers, so called "ballpoint screwdrivers" are available. They should only be used in conjunction with CAMLOG, CONELOG or iSy Titanium bases CAD/CAM free, crown, with an angulated screw channel.

3.9 Torque wrench for screws for Camlog implant systems

CAMLOG and iSy Torque wrenches are available for the fixation of screw-retained DEDICAM Prosthetics on Camlog implants with the specified tightening torque for the abutment or prosthetic screws and which can be adjusted continuously up to 30Ncm. The torque wrenches are made of corrosion-resistant stainless steel and are available non-sterile.

3.10 CAMLOG, CONELOG and iSy Titanium bases CAD/CAM, crown

For supporting DEDICAM Mesostructures, crowns and healing caps for bonding bases on a single CAMLOG, CONELOG or iSy Implant, rotationally locked titanium bases CAD/CAM, crown, are available with the implant-abutment connection, which is also used for one-piece abutments.

The CAMLOG Titanium base CAD/CAM, crown is also available in the PS ("Platform Switching") version for implant diameters 3.8, 4.3 and 5.0mm.

The shoulder height of the CAMLOG Titanium base CAD/CAM, crown is 0.3mm in the standard version (0.4mm for implant diameter 3.3mm) and 0.8mm in the PS version. CONELOG and iSy Titanium bases CAD/CAM, crown are available in the gingival heights 0.8 and 2.0mm.

The mesostructure, crown or the healing cap for bonding bases is fixated extra orally with a suitable fixation composite onto the previously blasted with aluminum oxide (grain size 50µm, max. 2bar pressure) and cleaned bonding surface of the titanium base CAD/CAM, crown, and preferably with the help of the included bonding aid.

Titanium bases CAD/CAM, crown, are made of titanium alloy and supplied non-sterile, each with an abutment screw and bonding aid. They are available for all CAMLOG, CONELOG and iSy Implants.

For the manual modeling of DEDICAM Structures on CAMLOG, CONELOG or iSy Titanium bases CAD/CAM, crown, rotationally locked modeling aids are available (see: 3.19 "Modeling aids for CAMLOG, CONELOG or iSy Titanium bases CAD/CAM, crown")

3.11 CAMLOG, CONELOG and iSy Titanium bases CAD/CAM, free, crown

The titanium bases CAD/CAM free enable a vertical shift of the screw channel in the mesostructure or crown of up to 25° in relation to the longitudinal axis. Exception: CONELOG and iSy Titanium base with gingival height 2.0mm is limited up to 15°. In addition, a horizontal deviation of the screw channel of up to 35° relative to the vertical cutout edges for the rotation lock is permissible.

They are available in two different chimney heights and for CAMLOG implants in a version for the Platform Switching option (PS). The bases are supplied with a sandblasted chimney, non-sterile, and each with an abutment screw and a lab screw. When used with an angulated screw access the new unused abutment screw must be placed before bonding and a specially configured Ballpoint screwdriver is required.

3.12 CAMLOG, CONELOG and iSy Titanium bases CAD/CAM, bridge

Non-locking titanium bases CAD/CAM bridge are available for supporting DEDICAM Bridges, bars and implant base structures for bonding bases on CAMLOG, CONELOG or iSy Implants.

The design with a highly shortened implant-abutment connection allows compensation of divergences in the implant axes up to 15° per implant to a joint reference axis.

The CAMLOG Titanium base CAD/CAM, bridge, is not available in the PS ("Platform Switching") version. The shoulder height of all CAMLOG CAD/CAM Titanium bases is 0.4mm. CONELOG and iSy Titanium bases CAD/CAM, bridge are available in gingival heights of 0.8 and 2.0mm.

The bridge, bar implant or base structure is fixated with a suitable fixation composite onto the previously blasted with aluminum oxide (grain size 50 µm, max. 2 bar pressure) and cleaned bonding surface of the titanium base CAD/CAM, bridge on the model and preferably with the help of the included bonding aid.

Titanium bases CAD/CAM, bridge, are made of titanium alloy and supplied non-sterile, each with an abutment screw and bonding aid. They are available for all CAMLOG, CONELOG and iSy Implants.

The restricted indications for the iSy Titanium base CAD/CAM for bridges is to be observed:

Bars with these abutments must be supported on at least 4 implants, bridges and implant base structures on at least 6 implants without mesial or distal extensions. Adequate polygonal support needs to be provided (cross arch).

3.13 Titanium bonding "Passive-Fit" bases and titanium caps without retention for bridges for CAMLOG and CONELOG Bar abutments

"Passive-Fit" titanium bonding bases and titanium caps without retention for bridges are available for the tension-free fixation of DEDICAM Bridges, bars to CAMLOG and CONELOG Bar abutments. They are screwed intraorally onto the bar abutments for fixation in the structure with light blue anodized prosthetic screws. Then the corresponding cavities in the bridge, bar or implant base structure are filled with fixation composite and the structure is placed on the titanium bonding bases or titanium caps until the fixation composite has cured. The seat for the prosthetic screw is located in the titanium bonding base or titanium cap.

"Passive Fit" titanium bonding bases and titanium caps without retention for bridges for CAMLOG and CONELOG Bar abutments are each available in two sizes (see following table).

Implant system	Implant Ø [mm]	"Passive-Fit" titanium bonding bases and titanium caps without retention for bridges; dimension of the prosthetic screw
CAMLOG	3.3	M1.6
	3.8	
	4.3	
	5.0	M2.0
	6.0	
CONELOG	3.3	M1.6
	3.8	
	4.3	
	5.0	M1.6 (for bar abutments with article number C2255.50xx)
		M2.0 (for bar abutments with article number C2254.50xx) C2256.50xx, C2257.50xx, C2258.50xx, C2259.50xx)

Table 8: Assignment of the screw dimension for "Passive-Fit" titanium bonding bases and titanium caps without retention for bridges for CAMLOG and CONELOG Bar abutments to the implant diameter

The titanium bonding bases for "Passive Fit" bar abutments and titanium caps without retention for bridges are made of titanium alloy and are available non-sterile and individually packaged.

When bonding the titanium bonding bases or titanium caps, the plastic screw for bar abutments, hex, can be used as an optional fixation and bonding aid.

3.14 BioHorizons Internal Hybrid Abutment Bases, hexed.

The rotation-locked BioHorizons Internal Hybrid Abutment Base, hexed are available for the support of DEDICAM Meso-/suprastructures and healing caps for bonding bases on a single BioHorizons Tapered internal or Internal Implant.

The bonding bases also have the implant-abutment connection used for one-piece abutments, are made of titanium alloy and non-sterile, have a shoulder height of 0.5mm and are supplied with an abutment screw. Their design enables a vertical shift of the screw channel in the mesostructure or crown of up to 15° in relation to the longitudinal axis. A horizontal deviation of the screw channel relative to the sloped cutout edges is not possible.

The mesostructure, crown or the healing cap for bonding bases is fixated extra orally with a suitable fixation composite onto the previously blasted with aluminum oxide (grain size 50µm, max. 2bar pressure) and cleaned bonding surface of the titanium base.

For the manual modeling of DEDICAM Mesostructures or crowns on the BioHorizons Internal Hybrid Abutment Base, hexed, rotationally locked waxing sleeves are available in packs of 3 pieces.

3.15 BioHorizons Internal Hybrid Abutment Bases, non-hexed

The BioHorizons Internal Hybrid Abutment Bases, non-hexed without antirotational mechanism are available for the support of DEDICAM Bridges, bars and implant base structures for bonding bases on BioHorizons Tapered internal or Internal Implants.

The bonding bases have a strongly shortened, rotationally symmetrical implant-abutment connection, are made of titanium alloy and are non-sterile, have a shoulder height of 0.5mm and are supplied with an abutment screw.

The bridge, bar or implant base structures is fixated with a suitable fixation composite onto the previously blasted with aluminum oxide (grain size 50µm, max. 2bar pressure) and cleaned bonding surface of the Internal Hybrid Abutment Bases, non-hexed on the model.

3.16 BioHorizons Bonding bases and titanium copings for Multi-unit abutments

Bonding bases with plastic sleeve, Passive Fit, for Multi-unit abutment as well as Titanium copings, for Multi-unit abutment are available for the tension-free fixation of DEDICAM Bridges, bars and implant base structures on Camlog or BioHorizons Multi-unit abutments. They are screwed intraorally onto the Multi-unit abutments for fixation in the structure with light blue anodized prosthetic screws. Then the corresponding cavities in the bridge, bar or implant base structure are filled with fixation composite and the structure is placed on the bonding bases/titanium copings until the fixation composite has cured. The seat for the prosthetic screw is located in the bonding base/titanium coping.

The bonding bases/titanium copings are made of titanium alloy, and are available non-sterile and individually packaged with the prosthetic screws.

When bonding the bonding bases/titanium copings, the prosthetic screw, Multi-unit, long can be used as an optional fixation aid (this screw has to be ordered separately).

3.17 Disconnectors for CONELOG and iSy Abutments

Disconnectors are available for loosening a jammed CONELOG or iSy Abutment, healing cap or an iSy Implant base from the lab analog or implant.

After removing the fixing screw from the abutment, healing cap or implant base, the abutment is pressed out of the lab analog or implant by turning the disconnector clockwise. If the conical connection does not come loose this way, the torque wrench (locked setting) can be placed on the disconnector and the fabricated abutment or healing cap can be loosened.

Disconnectors are made of corrosion-resistant stainless steel, and are available non-sterile and for all CONELOG and iSy Abutments and healing caps or iSy Implant bases. The abutment disconnectors are available in different lengths for the different clinical situations.

Implant system	Implant Ø [mm]	Disconnecter, thread
CONELOG	3.3	M1.6
	3.8	
	4.3	
	5.0	M2.0
iSy	3.8	M1.6
	4.4	
	5.0	

Table 9: Assignment of disconnectors to the implant system and diameter

3.18 Modeling aids for CAMLOG, CONELOG or iSy Titanium bases CAD/CAM, crown

For the manual modelling of DEDICAM Mesostructures or crowns on CAMLOG, CONELOG or iSy Titanium bases CAD/CAM, crown, rotationally locked modeling aids are available.

A wax-up is created from an adhesive wax or plastic construction on the modeling aid, that can be digitally scanned and imported into suitable CAD software as a three-dimensional data record for further processing. Please note that scanning a wax-up as a design template does not have the same surface quality for the final product as a purely virtual design.

Modeling aids are made of POM, and are available non-sterile for all CAMLOG, CONELOG or iSy Titanium bases CAD/CAM, crown (not for the “bridge” and “free” versions).

3.19 Locator Abutment with thread M2.0

A Locator Abutment with an M2.0 thread is available for anchoring bar-supported prostheses. A suitable lab kit with a processing replacement male without resilience – as soft tissue does put pressure on the prosthesis –, a spacer/duplication aid and a processing replacement male as spare parts is also available. The existing Locator instruments J2253.0001 and J2253.0002 can be used to install and remove the products (see also CAMLOG or CONELOG or product catalog or the iSy Product instructions for use). The Locator Abutment must be screwed into the bar with a tightening torque of 30Ncm. The use of at least three abutments is recommended for the fixation of a full denture.

4 Clinical benefit

4.1 Implant borne prosthetics and healing caps

Implant borne prosthetic components and healing caps are used for dental prosthetic treatment, which allows the replacement of missing teeth and/or the restoration of masticatory function and the corresponding improvement in quality of life.

4.2 Tooth borne restorations

Tooth borne restorations allow the oral rehabilitation of patients with tooth defects, and/or the restoration of masticatory function and the corresponding improvement in quality of life.

5 Undesirable side-effects

The use of DEDICAM Products is part of dental restoration procedure which may be associated with typical side effects such as:

Temporary side effects:

- pain, swelling, muco-gingival inflammation
- speech difficulties

Prolonged or permanent side effects:

- negative impact on adjacent teeth
- fractures of the restoration
- esthetic problems

6 Intended user and patient population

6.1 Intended user

The products are intended to be used by professional users (dentists, dental/maxillofacial surgeons, dental technicians) only.

6.2 Intended patient population

The products are to be used in skeletally mature patients subject to dental prosthesis treatment. As custom-made devices, they are intended for the sole use of a predefined particular patient exclusively to meet his or her individual conditions and needs. Please make sure that the contraindications and precautions are considered.

7 Preparation

7.1 Handling of packaging

The packaging must be checked before opening for damage and for the expiration date and should only be opened immediately before the products are used.

Should the original packaging be damaged, items may not be returned to the manufacturer.

The prosthetic components are supplied non-sterile unless they are explicitly marked as sterile.

Warning: Products must not be used if the product cannot be identified.

Caution: Products must not be used if the packaging has been damaged or previously opened.

7.2 Handling of non-sterile products and reprocessing

The DEDICAM Products are delivered non-sterile and must be cleaned and disinfected before and after each use on a patient.

In addition, healing caps, impression posts and – if not set after completed healing phase, or if set subcrestal – abutments, implant-supported bridges, bars and implant base structures, must be sterilized.

Warning:

- DEDICAM products intended for single use must not be reused on the patient because safe preparation and/or functional safety cannot be ensured.

Caution:

- All non-sterile packaged products must not be sterilized in the original packaging.
- Plastic components must not be sterilized (only cleaned and disinfected), except for DEDICAM products made out of PEEK or Telio CAD. These products can be sterilized.

7.3 Cleaning, Disinfection and Sterilization

7.3.1 General notes

Thorough cleaning and disinfection is indispensable for effective sterilization.

To ensure effective cleaning, remove the detachable attachments before starting cleaning.

Always observe the applicable legal national regulations and regulations on hygiene relating to dental/medical practices and hospitals. This applies in particular to the guidelines regarding effective prion

inactivation.

Note: every deviation from the processes quoted here must be evaluated carefully by the processor for effectiveness and possible negative consequences. The operator is responsible for ensuring that actual reprocessing achieves the desired results with the equipment, materials and personnel used in the reprocessing facility. This requires validations and routine monitoring.

For sterilizer and disinfectant:

As part of your responsibility for the sterility of the products in application, please observe the following:

- in general, use only adequately validated methods specific to the equipment and product for cleaning/disinfection and sterilization,
- regularly check and service the equipment used (thermal disinfectant, sterilizer) and
- that the validated parameters are complied with in each cycle.

Cleaning and disinfection device

When choosing the cleaning and disinfecting device, care should be taken that:

- the cleaning and disinfecting device has a proven efficacy as a matter of principle (CE marking, conforming to ISO 15883) and is validated specifically for the equipment and product,
- a validated program for thermal disinfection is used (A0 value > 3000 or – for older devices – at least 5 minutes at 90°C/194°F),
- the program used is suitable for the products and contains sufficient rinsing cycles (at least three draining rinsing cycles after the cleaning step or control of the rinsing cycles by conductivity measurement are recommended),
- only deionized water with a low bacterial count (maximum 10 bacteria/ml) and low endotoxin count (maximum 0.25 endotoxin units/ml) is used for rinsing,
- the air used for drying is filtered and
- that the cleaning and disinfecting device is regularly checked and serviced. Use only thermal disinfection that does not require a disinfecting agent.

Cleaning agents for mechanical cleaning and disinfection

Care should be taken that

- these are as a matter of principle suitable for parts made of metals, ceramics and plastics,
- the cleaning agents used are compatible with the devices

Camlog recommends the use of cleaning agents that do not require neutralizing agents.

Cleaning and disinfecting agents for manual cleaning and disinfection

Care should be taken that

- these are suitable for cleaning and disinfecting dental prosthetics,
- these are compatible with the materials of the products to be cleaned and disinfected,
- these are suitable for ultrasonic cleaning (no foam development),
- the disinfectant has a proven efficacy for disinfection (e.g. VAH/DGHM-listed, CE marking),
- cleaning agents and disinfectants are compatible with each other.

The concentrations and application times as well as instructions for post-rinsing specified by the manufacturer must be observed.

Instead of water, suitable cleaning agents or suitable cleaning agents with a disinfectant effect (for performing initial disinfection for personnel protection) can be used for preparation at the point of use (sections 7.3.2.1 and 7.3.3.1). These agents must be formulated and be suitable for the purpose of pre-cleaning soiled prosthetic components (no protein-fixing effect) and must be compatible with the agents used subsequently, as otherwise the effect of subsequent steps may be impaired.

Sterilization

Only steam sterilization methods with fractionated vacuum (forced venting) are permissible.

Care should be taken that:

- the maximum steam sterilization temperature does not exceed 138°C/280°F (products made of Telio CAD: 121°C/250°F),
- the steam sterilizers used have CE marking and comply with the requirements of EN 13060 or EN285,
- only adequately validated methods specific to the equipment and product according to ISO 17665 are used as a matter of principle for sterilization,
- the sterilizers used are regularly checked and serviced and
- that the validated parameters are complied with in each cycle.

Caution:

- Products made of Telio CAD may only be sterilized with the steam sterilization method at 121°C/250°F with a 20 minute hold time.
- One-piece abutments and healing caps made of zirconia may not be sterilized for longer than 6 minutes as well as only twice as a maximum.
- All DEDICAM Products must not be sterilized in the original packaging.
- Plastic parts, e. g. Preci matrices, impression cap and bite registration cap must not be sterilized.

Biotech Dental Digital/ Camlog recommends using disposable sterilization packaging with CE labeling. It must be ensured that the sterilization packaging is suitable for steam sterilization (i.e. constant temperature of at least 142°C/288°F, sufficient vapor permeability) and that the products are adequately protected against mechanical damage.

Two procedures are described in the following for cleaning and disinfection prior to sterilization: a mechanical procedure under section 7.3.2 and a manual procedure under section 7.3.3. Wherever possible, the mechanical procedure described in section 7.3.2 should be employed.

7.3.2 Mechanical cleaning and disinfection

7.3.2.1 Preparation at the point of use

Immediately after use in the patient's mouth, remove coarse contaminations and place the device in a bath with water.

7.3.2.2 Preliminary cleaning (max 2 hours after use)

(after use in the patient's mouth or contamination in the dental laboratory)

1. Remove contaminants from the prosthetic components under running water using a soft brush (no metal bristles or steel wool). The brush must be used exclusively for this purpose only. Brush until no visible contaminants remain. Some prosthetic components have cavities. Remove all residues in these. Contamination should no longer be visible.
2. Rinse all cavities at least 3 times with 20 ml water.

7.3.2.3 Transport of the pre-cleaned parts

Transport the pre-cleaned parts in a contamination-proof container.

7.3.2.4 Cleaning and disinfecting the pre-cleaned products in the cleaning and disinfection device

Worksteps

1. Use a thermal disinfectant.
2. Use a cleaning agent.
3. Place the pretreated prosthetic components in the thermal disinfectant using a small parts basket. The prosthetic components are not to come into contact with one another.
4. Start the program.
5. Remove the prosthetic components from the thermal disinfectant after program end.

7.3.2.5 Drying

Dry the prosthetic components. Use a loose lint free disposable cloth and compressed dry air, free of oil and with a low bacterial count. We also recommend the use of a sterile filter.

7.3.2.6 Inspection

Inspect the prosthetic components for signs of corrosion, surface damage, chipping or contamination. Damaged prosthetic components cannot be used. Re-clean and disinfect prosthetic components which are still contaminated.

7.3.2.7 Packaging

Pack prosthetic components for sterilization promptly. We recommend using disposable sterilization packaging with CE labeling. It must be ensured that the sterilization packaging is suitable for steam sterilization (i.e. constant temperature of at least 141°C/286°F, sufficient vapor permeability) and that the products are adequately protected against mechanical damage.

7.3.3 Manual cleaning and disinfection

7.3.3.1 Preparation at the point of use

Immediately after use in the patient's mouth, remove coarse contaminations and place the devices in a bath with water

7.3.3.2 Preliminary cleaning (maximal 2 hours after use)

1. Remove contaminants from the prosthetic components under running water using a soft brush (no metal bristles or steel wool). The brush must be used exclusively for this purpose only. Brush until no visible contaminants remain. Some prosthetic components have cavities. Remove all residues in these. Contamination should no longer be visible.
2. Rinse all cavities at least 3 times with 20 ml water.

7.3.3.3 Transport of the pre cleaned parts:

Transport the pre-cleaned parts in a contamination-proof container.

7.3.3.4 Cleaning

Caution: Difficult to access areas such as indentations and cavities require special attention during cleaning.

Worksteps:

1. All parts must be placed completely covered in the cleaning agent bath for the prescribed exposure time. The parts are not to come into contact with one another during the exposure time.
2. Brush the parts at the beginning of the exposure time for at least 30 seconds each (do not use metal bristles or steel wool). Brush all accessible cavities additionally for at least 30 seconds each.
3. Rinse all cavities at least 3 times with a syringe with at least 20 ml cleaning agent each at the beginning of the exposure time.
4. At the end of the exposure time, brush all prosthetic components for 30 seconds under running water, additionally brush all accessible cavities for 30 seconds each.

7.3.3.5 Rinsing after cleaning

Worksteps:

1. Rinse all parts completely under running water for at least one minute each. Also rinse areas that have limited access, e.g. cavities of the devices
2. Rinse all cavities 3 times with a syringe with at least 20 ml water each.

7.3.3.6 Disinfection

Worksteps:

1. Place all parts completely covered in the disinfectant bath for the exposure time specified by the disinfectant manufacturer. The parts are not to come into contact with one another

during the exposure time.

2. Rinse all cavities at least once with a syringe with at least 20ml disinfectant each at the beginning of the exposure time.

7.3.3.7 Rinsing after disinfection

Deionized water with a low bacterial count (max. 10 bacteria/ml) and low endotoxin (max. 0.25 endotoxin units/ml) is to be used for these rinses.

Worksteps:

1. Rinse all parts completely under running water for at least one minute each. Also rinse areas that have limited access, e.g. cavities.
2. Rinse all cavities 3 times with a syringe with at least 20ml water each.

7.3.3.8 Drying

Dry the prosthetic components. Use a loose lint free disposable cloth and compressed dry air, free of oil and with a low bacterial count. We also recommend the use of a sterile filter.

7.3.3.9 Inspection

Inspect the prosthetic components for signs of corrosion, surface damage, chipping or contamination. Damaged prosthetic components cannot be used. Re-clean and disinfect prosthetic components which are still contaminated.

7.3.3.10 Packaging

Pack prosthetic components for sterilization promptly. We recommend using disposable sterilization packaging with CE labeling. It must be ensured that the sterilization packaging is suitable for steam sterilization (i.e. constant temperature of at least 141°C/286°F, sufficient vapor permeability) and that the products are adequately protected against mechanical damage.

7.3.3.11 Sterilization

The following parameters apply: (fractionated vacuum/forced venting only):
at least 4 minutes at 134°C/273°F.

For devices made out of Telio CAD:
20 minutes at 121°C/250°F.

Caution:

- Products made of Telio CAD may only be sterilized with the steam sterilization method at 121°C/250°F with a 20-minute hold time.
- One-piece abutments and healing caps made of zirconia may not be sterilized for longer than 6 minutes as well as only twice as a maximum.
- All DEDICAM Products must not be sterilized in the original packaging.
- Plastic parts, e.g. Preci matrices, impression cap and bite registration cap must not be sterilized.

The recommendations of national Guidelines have to be observed.

Biotech Dental Digital/ Camlog recommends a minimum drying time of 30 minutes for each of the cycles described above.

8 Material

Products	Material			
One-piece abutments and healing caps made of titanium	Titanium alloy Ti-6Al-4V ELI according to ASTM F136 & DIN EN ISO 5832-3; Chemical composition (wt. %):			
	Al	5.5 - 6.5	C	< 0.4
	V	3.5 - 4.5	N	< 0.4
	Fe	≤ 0.3	H	< 0.4
	O	< 0.4	T	Balance

Crowns, bridges, bars, implant base structures and suprastructures for bars made of titanium	Titanium alloy Ti-6Al-4V ELI according to ASTM F136 & DIN EN ISO 5832-3; Chemical composition (wt. %): <table><tr><td>Al</td><td>5.5 - 6.5</td><td>C</td><td>≤ 0.08</td></tr><tr><td>V</td><td>3.5 - 4.5</td><td>N</td><td>≤ 0.05</td></tr><tr><td>Fe</td><td>≤ 0.3</td><td>H</td><td>≤ 0.015</td></tr><tr><td>O</td><td>≤ 0.2</td><td>T</td><td>Balance</td></tr></table>	Al	5.5 - 6.5	C	≤ 0.08	V	3.5 - 4.5	N	≤ 0.05	Fe	≤ 0.3	H	≤ 0.015	O	≤ 0.2	T	Balance
Al	5.5 - 6.5	C	≤ 0.08														
V	3.5 - 4.5	N	≤ 0.05														
Fe	≤ 0.3	H	≤ 0.015														
O	≤ 0.2	T	Balance														
Crowns, bridges, bars, implant base structures and suprastructures for bars made of CoCr	CoCr alloy Type 4 according to ISO 22674; Chemical composition (wt. %): <table><tr><td>Co</td><td>61.0%</td></tr><tr><td>Cr</td><td>27.5%</td></tr><tr><td>W</td><td>8.5%</td></tr><tr><td>Si</td><td>1.6%</td></tr><tr><td>Other components: C, Fe, Mn</td><td><1.0%</td></tr></table> Warning: Contains cobalt, which is considered a hazardous substance. During machining or thermal processing of CoCr material (e.g. milling, grinding) dust and fumes can be produced which could entail a health risk if inhaled. Therefore, suitable measures must be applied to reduce dust and vapor exposure, such as the use of extraction and ventilation systems or, for short-term use, suitable personal protective equipment for respiratory air filtration. Caution: Using different alloys in the patient's mouth can lead to galvanic effects.	Co	61.0%	Cr	27.5%	W	8.5%	Si	1.6%	Other components: C, Fe, Mn	<1.0%						
Co	61.0%																
Cr	27.5%																
W	8.5%																
Si	1.6%																
Other components: C, Fe, Mn	<1.0%																
One-piece abutments and healing caps, crowns, meso-structures for titanium bases and bridges made of IPS e.max ZirCAD LT	Pre-sintered, yttrium-stabilized, monochromatic (LT, MT) or polychromatic(MT Multi) zirconia, Type II/Class 5 according to ISO 6872; Chemical composition (wt. %): <table><tr><td>Y₂O₃</td><td>4.5 – 6.0 (MT/MT Multi: 6.5 – 8.0)</td></tr><tr><td>HfO₂</td><td>≤ 5.0</td></tr><tr><td>Al₂O₃</td><td>≤ 1.0</td></tr><tr><td>Other oxids</td><td>≤ 1.0</td></tr><tr><td>ZrO₂</td><td>balance</td></tr></table>	Y ₂ O ₃	4.5 – 6.0 (MT/MT Multi: 6.5 – 8.0)	HfO ₂	≤ 5.0	Al ₂ O ₃	≤ 1.0	Other oxids	≤ 1.0	ZrO ₂	balance						
Y ₂ O ₃	4.5 – 6.0 (MT/MT Multi: 6.5 – 8.0)																
HfO ₂	≤ 5.0																
Al ₂ O ₃	≤ 1.0																
Other oxids	≤ 1.0																
ZrO ₂	balance																
Crowns, meso-structures for titanium bases, bridges, veneers, inlays and onlays made of IPS e.max CAD	Lithium disilicate glass-ceramic; Chemical composition (wt %): <table><tr><td>SiO₂</td><td>57.0 – 80.0</td></tr><tr><td>Li₂O</td><td>11.0 – 19.0</td></tr><tr><td>K₂O</td><td>0.0 – 13.0</td></tr><tr><td>P₂O₅</td><td>0.0 – 11.0</td></tr><tr><td>ZrO₂</td><td>0.0 – 8.0</td></tr><tr><td>ZnO</td><td>0.0 – 8.0</td></tr><tr><td>Al₂O₃</td><td>0.0 – 5.0</td></tr><tr><td>MgO</td><td>0.0 – 5.0</td></tr></table>	SiO ₂	57.0 – 80.0	Li ₂ O	11.0 – 19.0	K ₂ O	0.0 – 13.0	P ₂ O ₅	0.0 – 11.0	ZrO ₂	0.0 – 8.0	ZnO	0.0 – 8.0	Al ₂ O ₃	0.0 – 5.0	MgO	0.0 – 5.0
SiO ₂	57.0 – 80.0																
Li ₂ O	11.0 – 19.0																
K ₂ O	0.0 – 13.0																
P ₂ O ₅	0.0 – 11.0																
ZrO ₂	0.0 – 8.0																
ZnO	0.0 – 8.0																
Al ₂ O ₃	0.0 – 5.0																
MgO	0.0 – 5.0																
Crowns, bridges and healing caps made of Telio CAD	Polymethylmethacrylate (PMMA); Chemical composition (wt %): <table><tr><td>PMMA</td><td>balance</td></tr><tr><td>TiO₂</td><td>≤ 0.50</td></tr><tr><td>Other pigments</td><td>≤ 0.35</td></tr></table> Caution: Wearing period maximum 12 months.	PMMA	balance	TiO ₂	≤ 0.50	Other pigments	≤ 0.35										
PMMA	balance																
TiO ₂	≤ 0.50																
Other pigments	≤ 0.35																

One-piece healing caps and impression posts made of PEEK	Polyetheretherketon (PEEK) Caution: Wearing period maximum 180 days.
--	--

Table 10: Material data

9 Storage conditions and Disposal

The product in its undamaged, original packaging are not subject to any special considerations for storage or handling (during transport and storage). It is recommended to keep the packaged products dry, out of direct sunlight and at room temperature to maintain an optimal visual appearance of the packaging.

Products to be disposed of must be treated and decontaminated as dental surgery waste in compliance with the relevant regulations.

10 General warnings and precautions

- Improper procedures when using products of the implant systems from CAMLOG Biotechnologies GmbH can lead to failures such as implant loss, local bone loss or unsatisfactory esthetic results. The products must only be used by dentists, surgeons and dental technicians. Every patient must be thoroughly examined and evaluated with regard to his/her radiographic, psychological and physical status. This includes the condition of their teeth and deficits in the related hard and soft tissues that might jeopardize the final outcome. Close collaboration between dentists, surgeons and dental technicians is essential for successful treatment. The implant systems from Camlog and the corresponding procedures are developed and clinically tested by experts in the field. Detailed information about the choice of suitable implants, prosthetic components, treatment planning and use of the products from Camlog is contained in the user information available at <https://ifu.camlog.com> or www.camlog.com from CAMLOG Biotechnologies GmbH. In addition, Camlog regularly offers courses or technical consultations regarding the use of its products. Your local Camlog distributor will be glad to advise you.
- Since safe application of Camlog products requires specialized knowledge and skills in oral implant dentistry, these products are sold only to dentists/surgeons and dental laboratories or on their prescription.
- If needed, technical advice/training on the proper use of products from Camlog is provided orally, in writing, by electronic media and/or by demonstration. The delivered information represents the scientific and technological State of the Art at the point in time the products are placed on the market. This does not exempt the user from the responsibility of personally testing the products for suitability for the intended purposes, indications and procedures. Handling and use of the products take place outside the control of Biotech Dental Digital SAS/ Camlog and are the direct responsibility of the user. All liability for damages resulting from such use is disclaimed by Biotech Dental Digital SAS/ CAMLOG Biotechnologies GmbH.
- The use of incompatible components can impair the function and safety of DEDICAM products. Biotech Dental Digital SAS/ CAMLOG Biotechnologies GmbH will provide neither warranty nor replacement if incompatible components are used.
- Implant-supported DEDICAM Products are – in as far as specified – only intended for specific implants/abutments and implant/abutment diameters. Their use with other implants/abutments or diameters can lead to the mechanical failure of system components, tissue damage or unsatisfactory aesthetic results.

- The implant systems from Camlog are each part of a comprehensive treatment concept and must be used only with the pertinent original parts and tools complying with the recommendations and instructions for use provided by the manufacturer. The components of the respective implant systems are matched precisely to one another. The use of third-party components may affect the function and safety of the system. Biotech Dental Digital SAS/ CAMLOG Biotechnologies GmbH neither offer warranty nor provide replacement services when non-system components are used.
- Drills and instruments are dedicated for specific implant systems and implant diameters. Use with different implant systems and/or implant diameters can lead to tissue injury or esthetically unsatisfactory results. Only the system-relevant drills, instruments and prosthetic components or those approved by Biotech Dental Digital SAS/ CAMLOG Biotechnologies GmbH may be used with the respective implant systems.
- Because of the small sizes involved, it may happen that a product is swallowed and/or aspirated. Aspiration can lead to dyspnea and in the worst case to asphyxiation. For this reason, the products should be secured appropriately to prevent them from being swallowed or aspirated during intraoral use.
- Where indications are listed for a particular product, it should be noted that any indications that are not listed are in fact contraindicated.
- Within the framework of our sales and delivery conditions, we guarantee the high quality of our products.
- Not all Camlog products and services are available in all countries.
- Damaged and/or corroded products must not be used.
- Unidentifiable products or products the function of which is doubtful, e.g., by poor readability of the markings and/or inscriptions, must not be used but are to be replaced.
- In case of hypersensitivity to any of the materials as listed in section 8 Material, refrain from using a restoration with the respective products.
- Custom-made devices are made in accordance with a written prescription from a health care professional only who is providing patient specific design characteristic under his or her responsibility.

Reporting in case of a serious incident

Any serious incident that has occurred in relation to the Biotech Dental Digital device must be reported to the manufacturer (Biotech Dental Digital SAS) and the competent authority of your country.

A serious incident means any malfunction or deterioration in the characteristics or performance of the implant or restoration that indirectly led, might have led or might lead to the temporary or permanent serious deterioration of the patient's state of health or the patient's death.

Trademarks and copyright

Protected trade names (trademarks) are not always specially indicated. The absence of such an indication does not mean that this name is NOT a trademark. The document including all its parts is protected by copyright. You may download the content regarding the intended use, but changes to or reproduction of the content are not permitted. Any exploitation beyond the narrow limits of the copyright act is not permitted without prior written approval of CAMLOG Biotechnologies GmbH and/or Biotech Dental Digital SAS is subject to legal sanctions.

11 MRI safety information

Caution:

DEDICAM Products made of titanium alloy and CoCr alloy have not been tested for safety and compatibility in the MR environment, such as heating, migration, or image artifacts. The safety of MR examinations in patients following restorations with metallic DEDICAM Products is unknown. Scanning a patient who has this device may result in patient injury.

Non-clinical testing has demonstrated the “Dental Implants CAMLOG, CONELOG and iSy” are MR conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 Tesla and 3 Tesla, **with**
- Maximum spatial field gradient of 25,700 G/cm (257 T/m)
- Maximum **force product of 438,000,000 G²/cm (438 T²/m)**
- **Calculated** maximum **head** averaged (HA) specific absorption rate (**SAR**) of 2.4 W/kg, **1.5 Tesla** (Normal Operating Mode)

2.3 W/kg, **3 Tesla** (Normal Operating Mode)

Under the scan conditions defined above, the worst-case implants for RF-induced heating are expected to produce a maximum temperature rise of less than

6.0 °C/10.8 °F (**2.4 W/kg, 1.5 Tesla**) **RF-related temperature increase**

6.0 °C/10.8 °F (**2.3 W/kg, 3 Tesla**) **RF-related temperature increase**

after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends approximately 23.3mm from “Dental Implants CAMLOG, CONELOG and iSy” when imaged with a gradient echo pulse sequence and a 3 Tesla MR system.

12 Further information

The latest version of this instruction for use is available at <https://ifu.biotech-dental-digital.com>.

Further information on all DEDICAM Products can be found in the current DEDICAM Product catalog. This document is available at www.camlog.com or from your local distributor.

Basic UDI-DI numbers are not applicable as DEDICAM Products are custom-made devices.

13 Color-Coding for the CAMLOG, CONELOG and BioHorizons Implant

systems according to Implant / Platform Ø

CAMLOG

Implant Ø	3.3mm	3.8mm	4.3mm	5.0mm	6.0mm
Color	gray	yellow	red	blue	green

CONELOG

Implant Ø	3.3mm	3.8mm	4.3mm	5.0mm
Color	gray	yellow	red	blue

iSy

Implant Ø	3.8 mm	4.4 mm	5.0 mm
Color	yellow green		

BioHorizons Tapered Internal and Internal implants

Platform Ø	3.0mm	3.5mm	4.5mm	5.7mm
Color	gray	yellow	green	blue

14 Explanation of symbols

The following symbols per ISO 15223-1 may be present on the device labelling or in information accompanying the device. Please refer to the device labelling or accompanying information for the applicable symbols.

	Consult instructions for use
	Caution, observe the warning notices
	Medical Device
	Article number
	Serial number
	Non-sterile
	Date of manufacture
	Do not reuse
	Do not use if package is damaged
	Manufacturer
	Contains hazardous substances

Table 11: Explanation of symbols on device labels

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BIOTECH DENTAL DIGITAL SAS
CAMLOG Vertriebs GmbH