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# Instructions for Use - Biotech Dental Digital Products (Prosthetic Phases)

**English**

---

This Instructions for use contains information on

## Contents

|        |   |    |
|--------|---|----|
| 1.     | General Notes on Product Design .....   | 42 |
| 2.     | Product description, including Intended Purpose, indications/contraindications, and application ..... | 42 |
| 2.1.   | One-piece titanium and cobalt-chromium abutments for BIOTECH DENTAL implant systems<br>42             |    |
| 2.1.1. | Product description .....   | 42 |
| 2.1.2. | Intended purpose.....   | 43 |
| 2.1.3. | Indications .....   | 43 |
| 2.1.4. | Warning.....  | 43 |
| 2.1.5. | Application.....  | 43 |
| 2.2.   | One-piece abutments for other implant systems.....  | 45 |
| 2.2.1. | Product Description .....   | 45 |
| 2.2.2. | Intended Purpose .....  | 46 |
| 2.2.3. | Indications .....   | 46 |
| 2.2.4. | Warning.....  | 46 |
| 2.2.5. | Application.....  | 46 |
| 2.3.   | Mesostructures and crowns for bonding to Titanium base (on single implants) .....                     | 48 |
| 2.3.1. | Product Description .....   | 48 |
| 2.3.2. | Intended Purpose .....  | 48 |
| 2.3.3. | Indications .....   | 48 |
| 2.3.4. | Warning.....  | 48 |
| 2.3.5. | Application.....  | 49 |
| 2.4.   | Bridges, bars, and implant abutments at the implant level (screw-retained) .....                      | 50 |
| 2.4.1. | Product Description .....   | 50 |
| 2.4.2. | Intended Purpose .....  | 51 |
| 2.4.3. | Indications .....   | 51 |
| 2.4.4. | Warning.....  | 51 |
| 2.4.5. | Application.....  | 52 |
| 2.5.   | Bridges, bars, and implant-supported abutment structures (screw-retained) .....                       | 54 |
| 2.5.1. | Product Description .....   | 54 |
| 2.5.2. | Intended Purpose .....  | 54 |
| 2.5.3. | Indications .....   | 55 |
| 2.5.4. | Warning.....  | 55 |
| 2.5.5. | Application.....  | 55 |
| 2.6.   | Bridges, bars, and implant abutments for adhesive (screw-retained) bases .....                        | 56 |
| 2.6.1. | Product Description .....   | 57 |
| 2.6.2. | Intended Purpose .....  | 57 |
| 2.6.3. | Indications .....   | 58 |
| 2.6.4. | Warning.....  | 58 |
| 2.6.5. | Application.....  | 58 |
| 2.7.   | Crowns and bridges (bonded) .....   | 61 |
| 2.7.1. | Product description.....  | 61 |
| 2.7.2. | Intended Purpose .....  | 61 |
| 2.7.3. | Indications .....   | 61 |
| 2.7.4. | Contraindications.....  | 62 |

|        |  |    |
|--------|--|----|
| 2.7.5. | Application .....  | 62 |
| 2.8.   | Inlays, onlays, and partial crowns .....                                 | 63 |
| 2.8.1. | Product description .....  | 63 |
| 2.8.2. | Intended purpose .....   | 63 |
| 2.8.3. | Indications .....  | 63 |
| 2.8.4. | Warning .....  | 63 |
| 2.8.5. | Application .....  | 63 |
| 3.     | Products compatible with BIOTECH DENTAL DIGITAL implant prostheses ..... | 64 |
| 4.     | Clinical Benefits .....  | 64 |
| 4.1.   | Implant Prosthetics .....  | 64 |
| 4.2.   | Dental Restorations .....  | 64 |
| 5.     | Contraindications .....  | 64 |
| 6.     | Adverse side effects .....   | 64 |
| 7.     | Intended Purposer and patient population .....                           | 65 |
| 7.1.   | Intended User .....  | 65 |
| 7.2.   | Target patient population .....  | 65 |
| 8.     | Preparation .....  | 65 |
| 8.1.   | Handling of Packaging .....  | 65 |
| 8.2.   | Handling of Non-Sterile Products and Reprocessing .....                  | 65 |
| 8.3.   | Cleaning, Disinfection, and Sterilization .....                          | 65 |
| 8.3.1. | General notes .....  | 65 |
| 8.3.2. | Mechanical cleaning and disinfection .....                               | 67 |
| 8.3.3. | Manual cleaning and disinfection .....                                   | 67 |
| 9.     | Material .....   | 69 |
| 10.    | Storage Conditions and Disposal .....                                    | 71 |
| 11.    | Warnings and General Precautions .....                                   | 71 |
| 12.    | MRI Safety Information .....   | 72 |
| 13.    | More information .....   | 72 |
| 14.    | Explanation of symbols .....   | 73 |
| 15.    | Contact .....  | 73 |

# 1. General Notes on Product Design

- All BIOTECH DENTAL DIGITAL products described in this Instructions for use are custom-made devices designed virtually using appropriate dental CAD software. They are manufactured individually for each patient at the BIOTECH DENTAL DIGITAL manufacturing center and are therefore intended for single-use only.
- Custom-made devices are manufactured under the responsibility of a healthcare professional and in accordance with their written prescription. The prescription provides patient-specific design characteristics.
- Unless otherwise specified, for implant-supported structures, the BIOTECH DENTAL DIGITAL CAD library specific to the CAD system must be used. These structures may only be used for implants, prefabricated abutments, or Tibase attachments assigned via the CAD library.
- The design of healing caps, impression abutments, abutments, bars, and implant-supported crowns/bridges must follow the emergence profile to avoid excessive pressure on the gingiva.
- The design of abutments, bars, and implant-supported crowns/bridges must be adapted to the emergence profile to avoid excessive pressure on the gingiva.

The following applies to all BIOTECH DENTAL DIGITAL products:

- Prosthetic and physiological principles must be observed during design. In particular, the lower mechanical load-bearing capacity of implants with reduced diameter must be taken into account.

## 2. Product description, including Intended Purpose, indications/contraindications, and application

### **2.1. One-piece titanium and cobalt-chromium abutments for BIOTECH DENTAL implant systems**

#### **2.1.1. Product description**

BIOTECH DENTAL DIGITAL one-piece abutments for implants serve as customized connecting elements for cemented or adhesive-retained crowns with a single dental implant and for bridges with two or more dental implants.

BIOTECH DENTAL DIGITAL custom prosthetic components are transmucosal components, custom-designed for prosthetic reconstruction. These devices are available in TA6V Titanium and Cobalt-Chrome. Each of these devices is a custom-made medical device intended for a single patient. A device is therefore tailored to a specific patient and may under no circumstances be transferred to a third party. These prostheses have a lifespan, according to the state of the art, of 15 years.

The abutments are manufactured:

- Either by machining from a titanium or cobalt-chromium alloy blank with a prefabricated implant-abutment interface.
- Or by laser sintering from titanium or cobalt-chromium alloy metal powder.

The design requires scanning using the scanbodies stored in the BIOTECH DENTAL DIGITAL CAD library.

Each abutment includes a corresponding abutment screw. These can also be ordered individually.

All possible prosthetic restorations depending on the implant used are listed in tables available on the ordering website (scanbody references, analogs, screws, Titanium base, library reference, etc.). You can view them in the "Documents" section.

Abutments and abutment screws are non-sterile and intended for single use. Abutments are available in rough or polished finishes. For abutments, the polished surface is limited to the subgingival section.

#### **Materials:**

##### **Laser-sintered prostheses**

|                                 |  |
|---------------------------------|--|
| Titanium alloy (Ti-6Al-4V ELI): | CTE (25 to 500 °C): $\sim 10.16 \times 10^{-6} \text{ K}^{-1}$ |
| Chromium-cobalt alloy (Cr-Co)   | CTE (20 to 600 °C): $\sim 14.4 \times 10^{-6} \text{ K}^{-1}$  |

## Machined prostheses

|                                 |  |
|---------------------------------|--|
| Titanium alloy (Ti-6Al-4V ELI): | CTE (25 to 500 °C): $\sim 10.3 \times 10^{-6} \text{ K}^{-1}$                                      |
| Chromium-Cobalt Alloy (Cr-Co)   | CTE (20 to 600 °C): $\sim 14.0 \times 10^{-6} \text{ K}^{-1} / 14.4 \times 10^{-6} \text{ K}^{-1}$ |

### 2.1.2. Intended purpose

The abutments are destined for functional and/or aesthetic rehabilitation of partially or completely edentulous patients, as components of dental prostheses in combination with endosseous implant in the maxilla and/or the mandible.

### 2.1.3. Indications

The abutments are used for prosthetic reconstruction on the lower or upper jaw to restore partially or completely toothless sections. More specifically, for:

- Single-tooth gaps,
- Partially edentulous jaws with multiple missing teeth or
- Toothless jaws.

### 2.1.4. Warning

- Compliance with the clinical situations described in section 2.1.5.1 “Design” is required.
- Restorations with an implant length-to-restoration length ratio greater than 1:1.25.
- Primary abutments with splint.
- Double crowns on implants with a diameter of 3.3 mm.
- Single-unit restorations with open-ended bridges.
- Extension attachments.
- Direct laser welding.
- Casting technique.
- Modification of the implant-abutment connection, including the bearing surface of the implant shoulder.
- Direct full-ceramic veneer
- Single-tooth crown restorations for implants with diameters of 3.0 and 3.3 mm outside the region of the upper lateral incisors and the lower central and lateral incisors.

### 2.1.5. Application

#### 2.1.5.1. Design

The following requirements must be met when designing the abutment:

- Angle of the prosthetic axis relative to the screw axis for abutments up to a maximum of 30°
- The wall thickness in load-bearing areas is at least 0.6 mm
- The radius of the curves is at least 0.5 mm
- The construction height above the screw head for the abutment is at least 0.2 mm.
- Minimum apical diameter according to Table 1


Figure 1: Design limits of the abutment

#### Caution:

- To avoid notch effects and the risk of abutment fracture, abrupt transitions in the geometry must be

avoided. If necessary, the abutment must be discarded, or any features that could cause issues must be removed by polishing under microscopic examination.

- The lower mechanical load-bearing capacity of implants with reduced diameter must be taken into account.

#### **2.1.5.2. Post-treatment**

##### **Caution:**

- The contact surfaces between the abutment and the implant must not be abrasive, sandblasted, or mechanically treated. In this context, please note that discoloration of the titanium alloy can occur at temperatures above 300°C. Heating above 800°C causes structural changes resulting in reduced strength.
- The surface quality and wall thickness must not be altered or reduced in load-bearing areas.

If post-processing of the abutment is necessary, please note that:

- The implant-abutment connection, including the bearing surface of the implant shoulder, must be protected using an implant analog.
- With implant systems, fixation is performed using the corresponding laboratory screw
- That the design specifications are met (see 2.1.5.1 “Design”)

#### **2.1.5.3. Transfer to the clinician**

The post-treated abutment must be thoroughly cleaned of processing residues and any other contamination using an appropriate method.

The specific abutment screw supplied with the abutment must be clearly labeled and handed over to the recipient in its unopened original packaging.

#### **2.1.5.4. Insertion**

The final insertion of the abutment into the implant is performed after cleaning, disinfection, and the recommended sterilization (required for abutments if they are not placed after the healing phase is complete, or if they are placed subcrestally) as follows:

1. Remove the screw from the implant cover or the healing cap.
2. Clean the internal configuration of the implant and the implant shoulder according to the implant type.
3. Ensure that the contact surfaces between the implant and the abutment are absolutely clean and, in particular, that bone debris and soft tissue are completely removed.
4. Insert the abutment into the implant.
5. Screw in the new, unused abutment screw supplied with the abutment. Use the implant system screwdriver corresponding to this purpose.
6. Tighten the abutment screws for securing the abutments to the torque specified by the implant manufacturer in Ncm.

**Warning:**

- Laboratory screws must not be used on the patient.
- The abutment must be permanently secured using new, unused abutment screws.
- Do not mix abutment screws from different implant systems.

**Caution:**

- If the tightening torque exceeds the recommended value, this may result in damage to or destruction of the screw and/or damage to the abutment and/or the implant.
  - If the tightening torque is lower than the recommended value, this may cause the abutment to loosen.
7. Ensure that the implant fits properly. A follow-up X-ray is recommended to ensure that the abutment is correctly positioned on the implant. No soft tissue should be trapped, and there should be no visible gap between the abutment and the implant.
  8. Retighten to the same torque after at least 5 minutes.

**Warning:**

- The screw channel must be sealed with a suitable, removable material to secure the single hand-tightened abutment screw.
  - First, cover the screw head with a fully removable material such as cotton or gutta-percha.
9. In the case of immediate prosthetic restoration of the implant: the edges of the wound are closed tightly with atraumatic suture material. The sutures should not be tied too tightly. The gingiva must fit snugly against the abutment or healing cap in the case of transgingival healing.

**2.2. One-piece abutments for other implant systems****2.2.1. Product Description**

BIOTECH DENTAL DIGITAL one-piece abutments for implant systems from selected other manufacturers serve as custom connection elements for cemented or adhesive-retained crowns with a single dental implant and for bridges with two or more dental implants.

They can only be used for the implant assigned via the CAD library. Selection of compatible implant systems.

The abutments and healing caps are manufactured from a titanium alloy blank with a prefabricated implant-abutment interface and are not color-coded.

Abutments are manufactured:

- Either by machining from a titanium or cobalt-chromium alloy blank with a prefabricated implant-abutment interface.
- Or by laser sintering from titanium or cobalt-chromium metal powder

The design requires scanning using the scanbodies stored in the BIOTECH DENTAL DIGITAL CAD library.

Each abutment includes a corresponding abutment screw packaged separately. These can also be ordered individually.

The abutments and abutment screws are non-sterile and intended for single use. The abutments are available in rough or polished finishes. For the abutments, the polished surface is limited to the subgingival section.

**Material:****Laser-sintered prostheses**

|                                 |  |
|---------------------------------|--|
| Titanium alloy (Ti-6Al-4V ELI): | CTE (25 to 500 °C): $\sim 10.16 \times 10^{-6} \text{ K}^{-1}$ |
| Chromium-cobalt alloy (Cr-Co)   | CTE (20 to 600 °C): $\sim 14.4 \times 10^{-6} \text{ K}^{-1}$  |

**Machined prostheses**

|                                 |  |
|---------------------------------|--|
| Titanium alloy (Ti-6Al-4V ELI): | CTE (25 to 500 °C): $\sim 10.3 \times 10^{-6} \text{ K}^{-1}$                                    |
| Chromium-Cobalt Alloy (Cr-Co)   | CTE (20 to 600 °C): $\sim 14.0 \times 10^{-6} \text{ K}^{-1}/14.4 \times 10^{-6} \text{ K}^{-1}$ |

### 2.2.2. Intended Purpose

The abutments are destined for functional and/or aesthetic rehabilitation of partially or completely edentulous patients, as components of dental prostheses in combination with endosseous implant in the maxilla and/or the mandible.

### 2.2.3. Indications

The abutments are used for prosthetic reconstruction on the lower or upper jaw to restore partially or completely toothless sections. More specifically, for:

- single-tooth gaps,
- partially edentulous jaws with multiple missing teeth, or
- edentulous jaws.

### 2.2.4. Warning

- Clinical situations that do not allow compliance with design specifications (see 2.2.5.1 “Design”)
- Restorations with an implant length-to-restoration length ratio greater than 1:1.25
- Primary abutments with splint
- Double crowns on implants with a diameter of 3.3 mm
- Single-unit restorations with a free-end pontic
- Extension attachments
- Direct laser welding
- Casting technique.
- Modification of the implant-abutment connection, including the bearing surface of the implant shoulder
- Direct full-ceramic veneer

### 2.2.5. Application

#### 2.2.5.1. Design

The following requirements must be met when designing the abutment:

- Angle of the prosthetic axis relative to the screw axis for abutments up to a maximum of 30°
- The wall thickness in load-bearing areas is at least 0.6 mm.
- The radius of the curves is at least 0.5 mm.
- The construction height above the screw head is at least 0.2 mm for the abutment and at least 1.0 mm for the healing cap.
- Minimum apical diameter according to Table 1


Figure 2: Abutment design limits

**Caution:**

To prevent notch effects and the risk of abutment fracture, abrupt transitions in the individual geometry should be avoided. If necessary, the abutment should be discarded, or any features that could cause issues should be removed by polishing under microscopic examination.

The lower mechanical load-bearing capacity of implants with a reduced diameter must be taken into account.

**2.2.5.2. Post-treatment****Caution:**

- The contact surfaces between the abutment and the implant must not be abraded or mechanically treated. In this context, it should be noted that discoloration of the titanium alloy can occur at temperatures above 300°C. Heating to over 800°C results in structural changes and a reduction in strength.
- The surface quality and wall thickness must not be altered or reduced in load-bearing areas.

If further processing of the abutment is necessary, please take the following into account:

- The implant-abutment connection, including the bearing surface of the implant shoulder, must be protected by an implant analog
- Fixation is performed using the corresponding laboratory screw.
- Ensure that the design specifications are met (see 2.2.5.1 "Design")

**2.2.5.3. Return to the clinician**

The abutment must be thoroughly cleaned of processing residues and other contaminants using an appropriate method.

The specific abutment screw included with the abutment must be clearly labeled and delivered to the recipient in its unopened original packaging.

**2.2.5.4. Insertion**

The final insertion of the abutment into the implant is performed after cleaning, disinfection, and recommended sterilization (required for healing caps and abutments if they are not placed after a complete healing phase, or if they are placed in a subcrestal position) as follows:

10. Remove the implant cover screw or healing cap.
11. Clean the internal configuration of the implant and the implant shoulder according to the implant type.
12. Ensure that the contact surfaces between the implant and the abutment are perfectly clean and, in particular, that bone debris and soft tissue are completely removed.
13. Insert the abutment into the implant.
14. Screw in the new, unused BIOTECH DENTAL DIGITAL abutment screw supplied with the abutment. Use the screwdriver provided by the corresponding implant manufacturer for this purpose.
15. Tighten the abutment screw to the specified torque; the value should be compared to the corresponding specification in the implant manufacturer's instructions for use and reduced if necessary.

**Warning:**

- If the torque exceeds the recommended value, this may result in damage to or destruction of the screw and/or damage to the abutment and/or the implant.
  - If the torque is lower than the recommended value, this may cause the abutment to loosen.
16. Ensure that the implant is properly seated. A follow-up X-ray is recommended to verify that the abutment is correctly positioned on the implant. No soft tissue should be trapped, and there should be no visible gap between the abutment and the implant.
  17. Retighten to the same torque after at least 5 minutes.
  18. In the case of immediate prosthetic restoration of the implant: the wound edges are closed tightly using atraumatic suture material. The sutures must not be too tight. The gingiva must fit snugly against the abutment in the case of transgingival healing.

## 2.3. Mesostructures and crowns for bonding to Titanium base (on single implants)

### 2.3.1. Product Description

BIOTECH DENTAL DIGITAL mesostructures and crowns for bonding to Titanium base are used to replace lost tooth structure on a single implant. The mesostructures can also be used as abutments for bridges. They are bonded to a Titanium base compatible with the implant system and a crown. The structures feature a continuous screw hole and internal geometry corresponding to a bonding base assigned in the CAD system.

The mesostructures and crowns are designed virtually using appropriate CAD software and the BIOTECH DENTAL DIGITAL CAD library or freely available, unencrypted CAD libraries. They may only be used for bonding to a Tibase assigned via the CAD library.

Mesostructures and crowns for bonding to Tibase are supplied non-sterile. Structures designed for Tibase are supplied with the Tibase and abutment screw. Structures with an angled screw channel are supplied with the angled screw.

The frameworks are available for individual characterization in various tooth shades. Zirconia is also available in different degrees of translucency and, optionally, with a natural color gradient from dentin to incisal.

Temporary restorations with a maximum lifespan of 12 months can be fabricated using a tooth-colored plastic material (polymethyl methacrylate = PMMA).

Zirconia mesostructures and crowns for Tibase are delivered pre-sintered.

Crowns for adhesive bases are particularly suitable if the screw channel is located occlusally or palatally/lingually.

#### Material:

|   |  |
|---|--|
| IPS e.max ZirCAD (Zirconia) Prime & LT              | CTE (25 to 500°C): $10.5 \times 10^{-6} \text{K}^{-1}$   |
| IPS e.max ZirCAD (Zirconia) Prime<br>Prime Esthetic | CTE (25 to 500°C): $10.4 \times 10^{-6} \text{K}^{-1}$   |
| Polymethyl methacrylate (PMMA)                      | CTE: N/A   |
| IPS e.max CAD (Lithium disilicate)                  | CTE (100 to 400°C): $10.2 \times 10^{-6} \text{K}^{-1}$<br>CTE (100 to 500°C): $10.5 \times 10^{-6} \text{K}^{-1}$ |

CATALOG: All possible prosthetic restorations depending on the implant used are listed in tables available on the ordering website (scanbody references, analogs, screws, Titanium bases, library references, etc.). You can view them in the "Documents" section.

### 2.3.2. Intended Purpose

The mesostructures are destined for functional and/or aesthetic rehabilitation of partially or completely edentulous patients, as components of dental prostheses in combination with endosseous implant in the maxilla and/or the mandible.

Crowns are intended for functional and/or aesthetic rehabilitation of partially edentulous patients, as components of dental prostheses in the maxilla and/or the mandible.

### 2.3.3. Indications

The mesostructures are used for prosthetic reconstruction on the lower or upper jaw to restore partially or completely toothless sections.

Crowns are used for the replacement of partially edentulous jaws, dental defects or compromised teeth with a lack of tooth substance.

Specifically for:

- Single-tooth gaps,
- Partially edentulous jaws with multiple missing teeth, or
- Edentulous jaws.

### 2.3.4. Warning

- Clinical situation that does not allow compliance with design specifications (see 2.5.5.1 "Design")
- Reuse
- Primary abutments
- Restorations with a length-to-implant-length ratio greater than 1:1.25
- Restorations on Titanium base designed for multi-unit restorations and other incompatible third-party Titanium base
- Double crowns on reduced-diameter implants
- Single-tooth restorations with open-ended pontics
- Extension attachments on mesostructures and crowns
- A wear period exceeding 12 months for PMMA restorations

### 2.3.5. Application

#### 2.3.5.1. Design

The following requirements must be met when designing mesostructures and crowns for bonding to Tibase:

- Angle of the prosthetic axis relative to the screw axis up to 20°.
- Tibase CAD/CAM without crown:
- Vertical offset of the screw channel up to 25° relative to the longitudinal axis. Exception: Tibase with a gingival height of 2.0 mm is limited to 15°.
- Horizontal offset of the screw channel up to 35° relative to the vertical edges of the cutout for rotation locking
- Wall thickness:

|   | Mesostructure/crown<br>Zirconia | PMMA crown | Mesostructure/crown<br>IPS e.max CAD |
|---|---------------------------------|------------|--------------------------------------|
| Minimum wall thickness - circular around the screw channel [mm] | 0.7                             | 0.8        | 0.5                                  |
| Minimum wall thickness - circular on the base shoulder [mm]     | 0.5                             | 0.5        | 0.5                                  |
| Minimum wall thickness - incisal/occlusal [mm]                  | 0.7                             | 1.5        | --                                   |

**Table 1: Wall thicknesses of mesostructures and crowns for bonding to Tibase**

The lower mechanical load-bearing capacity of implants with reduced diameters must be taken into account.

#### 2.3.5.2. Post-treatment

PMMA crowns for Tibase cementations can be modified to optimize the emergence profile by applying a suitable composite or by adjusting the wall thicknesses to the minimum permissible levels through grinding.

**Caution:** Surface quality and wall thickness must not be altered or reduced in load-bearing areas.

The material manufacturers' instructions for use must be followed; they are available at:

|                                 |   |
|---------------------------------|---|
| Zirconia:                       | <a href="http://www.ivoclar.com/en_gb/eifu">www.ivoclar.com/en_gb/eifu</a>                                      |
| PMMA (polymethyl methacrylate): | <a href="https://www.white-peaks-dental.com/en/downloads/">https://www.white-peaks-dental.com/en/downloads/</a> |
| Lithium disilicate              | <a href="http://www.ivoclar.com/en_gb/eifu">www.ivoclar.com/en_gb/eifu</a>                                      |

After preparing the surface according to the respective specifications of the manufacturers of the materials and mounting composites, the structures are permanently bonded extraorally to the Titanium bases prepared for bonding, ensuring proper alignment and, if necessary, using bonding aids. Special attention must be paid to following these instructions: when using an with an angled screw, the new, unused abutment screw must be inserted before bonding.

For bonding, rotate the abutment on the titanium base until the anti-rotation features engage, then gently press it onto the titanium base. Excess bonding material must be removed immediately.

**Caution:** Contact between the cement and the gum is known to cause inflammation.

### **Appropriate surface conditioners and bonding composites**

To bond the mesostructures or crown to a Titanium base, we recommend using products that are suitable and validated for this purpose. Follow the manufacturer's instructions. Biotech Dental Digital SAS assumes no liability for the durability of the bond.

#### **2.3.5.3. Delivery to the clinician**

The mesostructure or crown of the Titanium base must be thoroughly cleaned of processing residues and other impurities using an appropriate procedure after bonding to the Titanium base as a hybrid abutment or hybrid crown.

The specific abutment screw supplied with a Titanium base must be clearly identified and delivered to the recipient in its unopened original packaging (does not apply to pre-assembled abutment screws in the case of angled screw channels).

#### **2.3.5.4. Insertion**

The final insertion of the hybrid abutment or hybrid crown into the implant is performed after cleaning, disinfection, and recommended sterilization (required if the abutment is not placed after a complete healing phase, or if it is placed subcrestally) as follows:

1. Remove the implant cover screw or healing cap.
2. Clean the internal configuration of the implant and the implant shoulder according to the implant type.
3. Ensure that the contact surfaces between the implant and the hybrid abutment or hybrid crown are absolutely clean and, in particular, that bone debris and soft tissue are completely removed.
4. Insert the hybrid abutment or hybrid crown into the implant.
5. Screw in the new, unused abutment screw corresponding to the implant system and the implant diameter using the appropriate screwdriver. Note: An angled screw channel requires the use of a specific screwdriver. It is recommended to use the standard screwdriver for straight screw channels.
6. Tighten the abutment screw to the torque specified by the manufacturers of the Titanium base and the implant.

#### **Warning:**

- If the torque exceeds the recommended value, it may result in damage to or destruction of the screw and/or damage to the hybrid abutment or hybrid crown and/or the implant.
  - If the tightening torque is lower than the recommended value, this may cause the hybrid abutment or hybrid crown to loosen.
7. Ensure that the implant fits properly. A follow-up X-ray is recommended to verify that the hybrid abutment or hybrid crown fits the implant properly. No soft tissue should be trapped, and there should be no visible gap between the abutment and the implant.
  8. Hybrid abutment: For hygiene reasons, the screw channel must be sealed with appropriate removable materials. First, cover the screw head with a fully removable material.
  9. Hybrid crown: Cover the screw head with a suitable, fully removable material. Then hermetically seal the screw channel with suitable materials, such as a composite.

## **2.4. Bridges, bars, and implant abutments at the implant level (screw-retained)**

### **2.4.1. Product Description**

BIOTECH DENTAL DIGITAL implant-level bridges are custom-made prosthetic restorations designed to replace lost tooth structure and are screwed onto two or more implants.

BIOTECH DENTAL DIGITAL implant-level bars serve as custom connecting elements for removable prosthetic restorations and are screwed onto two or more implants.

BIOTECH DENTAL DIGITAL implant-supported frameworks serve as custom-made metal support structures for plastic prosthetic dental restorations and prefabricated plastic or ceramic dentures, which are screwed onto two or more implants.

The frameworks are screwed directly onto the implants without additional abutments. They are available for Biotech Dental implant systems as well as for implants from other third-party suppliers.

In the case of conical implant-abutment connections, implant-level bridges, bars, and implant-level base structures are also supported by the implant shoulder.

Anatomical bridges and implant-level bridge frameworks can consist of 2 to 16 units. Bars and implant-level implant frameworks are available for 2 to 10 implants.

The design requires scanning using BIOTECH DENTAL scanbodies or scanbodies for other implant systems.

The structures are non-sterile and intended for single use. Bridges and implant-level abutments have a milled, rough surface; bars are also available polished.

The bridges, bars, and implant-level abutments are supplied with the corresponding abutment screws, packaged separately, in a quantity corresponding to the number of implants to be restored with the structure. Abutment screws can also be ordered individually. The implant-level bridges, bars, and implant abutments are secured to the implants using the supplied abutment screws and in accordance with the specified tightening torque defined by the manufacturer.

The provided CAD files for implant-level bridges, bars, and implant abutments are analyzed for implant axis angulation during fabrication (CAM software), among other things. The implant connections are thus individually optimized based on the overall design and insertion direction, so that guidance and positioning benefit from the necessary freedom on the one hand, and stability is ensured on the other. Thanks to this method of detecting interference contacts and subsequently eliminating them through milling, deviations in implant axes can be compensated for up to 50° depending on the situation and the implant system.

The anatomically reduced bridge framework can be veneered with suitable layering materials or with CAD/CAM-fabricated superstructures.

Thanks to the BIOTECH DENTAL DIGITAL CAD library, bridges and bars can be supplemented with attachments for securing removable restorations, and bars can also be designed with different profiles.

A suitable bar profile is available in the BIOTECH DENTAL DIGITAL CAD library for implant abutments.

**Material:**

**Laser-sintered prostheses**

|                                 |  |
|---------------------------------|--|
| Titanium alloy (Ti-6Al-4V ELI): | CTE (25 to 500 °C): ~ 10.16 x 10 <sup>-6</sup> K <sup>-1</sup> |
| Chromium-cobalt alloy (Cr-Co)   | CTE (20 to 600 °C): ~ 14.4 x 10 <sup>-6</sup> K <sup>-1</sup>  |

**Machined prostheses**

|                            |   |
|----------------------------|---|
| Titanium alloy (Ti-6Al-4V) | CTE (20 to 600°C): ~ 10.3 x 10 <sup>-6</sup> /K <sup>-1</sup> |
| CoCr Alloy                 | CTE (20 to 500°C): ~ 14.0 x 10 <sup>-6</sup> /K <sup>-1</sup> |

**2.4.2. Intended Purpose**

Bridges, bars, and implant abutments are destined for functional and/or aesthetic rehabilitation of partially or completely edentulous patients, as components of dental prostheses in combination with endosseous implant in the maxilla and/or the mandible.

**2.4.3. Indications**

BIOTECH DENTAL DIGITAL Bridges, bars, and implant-level abutments are used for prosthetic reconstruction on the lower or upper jaw to restore partially or completely toothless sections. Specifically for:

- Partially edentulous jaws with multiple missing teeth or
- Edentulous jaws.

**2.4.4. Warning**

- Clinical situation that does not allow compliance with design specifications (see 2.6.5.1 "Design")
- Use on an incompatible implant
- Support of the joint between the tooth (abutment) and the implant

## 2.4.5. Application

### 2.4.5.1. Design

#### **General remarks on direct-to-implant screw-retained prostheses compared to abutment-supported implant prostheses:**

We recommend evaluating each case individually during the design phase. Both restoration options have advantages and disadvantages that should be assessed based on the specific situation.

In the case of direct screw-retained implants (implant-abutment connection modified according to the design), the restoration is supported only by the implant shoulder, and the fixation screw is used for a significant portion of the overall force transmission. Long load-bearing arms and uniaxial distribution of the implants must be viewed critically. To this end, angulations of up to 50° are possible depending on the situation and the implant system, which allows for cost savings on construction components.

The mechanical stability of implant-supported abutment prostheses is considered more advantageous.

The following requirements must be met when designing bridges, bars, and implant abutments at the implant level:

- Minimum wall thickness: 0.5 mm
- Minimum cross-sections of connectors:
  - Bridge element/split crowns, anterior tooth: 6 mm<sup>2</sup>
  - Bridge element/split crowns, posterior tooth: 9 mm<sup>2</sup>
  - Custom-designed bars and implant abutments: 8 mm<sup>2</sup>
- Titanium alloy bridges: maximum of 3 connected pontics; in the case of single support on 2 implants, the maximum distance between implant axes is 30 mm.
- CoCr alloy bridges: maximum of 3 connected pontics in the lateral region and 4 in the anterior region; with single support on 2 implants, the maximum distance between implant axes is 30 mm.
- Bridges: maximum length of abutments; width of one molar, but no more than 10 mm.
- Implant bars and abutments: maximum distance between the two implant axes: 30 mm
- Implant bars and abutments: maximum length of extensions: half the length of the adjacent inter-implant segment, but no more than 10 mm.
- Indications and contraindications for implants, in accordance with the manufacturer's instructions, must be taken into consideration.
- The lower mechanical load-bearing capacity of implants with a reduced diameter must be taken into account.

### 2.4.5.2. Post-treatment

If no information on the surface treatment of bridge frameworks/implant abutments for the selected veneering/construction material is available from the manufacturer, bridge frameworks and implant abutments made of CoCr alloy must be sandblasted with pure aluminum oxide with a grain diameter of 110 µm - 250 µm minimum and a blasting pressure of 2 to 4 bar; titanium alloy bridge frameworks and implant abutments must be sandblasted with pure aluminum oxide with a grain size of 110 µm - 250 µm and a blasting pressure of 2 to 3 bar.

#### **Warning:**

- During the machining or heat treatment of CoCr material (e.g., milling, grinding), dust and fumes may be generated and pose a health risk if inhaled. Therefore, appropriate measures must be taken to reduce exposure to dust and fumes, such as the use of extraction and ventilation systems or, for short-term use, appropriate personal protective equipment for respiratory filtration.

Bridge frameworks and implant abutments must then be rinsed under running water, steam-cleaned, and degreased with ethyl alcohol (for titanium alloys, never use hydrofluoric acid for cleaning). Do not contaminate the structure afterward, for example by touching it with unprotected fingers.

Bridge frameworks are veneered with a coating material compatible with the base material (titanium or CoCr alloy) or with CAD/CAM-fabricated superstructures. The manufacturer's instructions for use must be followed. In this context, it is important to note that discoloration of the titanium alloy can occur at temperatures above 300°C.

#### **Caution:**

- The contact surfaces between the bridge, the bar, or the implant abutment and the implant must not be

abraded or treated.

- Heating titanium to temperatures above 800°C causes structural changes and a reduction in strength.
- The material manufacturers' processing instructions must be followed; they are available at:

- **Laser-sintered prostheses**

|  |   |
|--|---|
| CoCr alloy (Scheftner Starbond CoS Powder 45): | Please contact us by email at the following address: <a href="mailto:service@scheftner.dental">service@scheftner.dental</a> |
| Titanium alloy (TEKMAT™Ti64-53/20-A):          | Please contact us by email at the following address: <a href="mailto:sds@tekna.com">sds@tekna.com</a>                       |

- **Machined prostheses**

|  |   |
|--|---|
| CoCr alloy (Eisenbacher KERA DISC):          | Please contact us by email at the following address: <a href="mailto:info@eisenbacher.de">info@eisenbacher.de</a> |
| Titanium alloy (Eisenbacher Kera TI 5-Disc): | Please contact us by email at the following address: <a href="mailto:info@eisenbacher.de">info@eisenbacher.de</a> |

To avoid aesthetic issues, CoCr restorations covered with a layering material must have the oxide layer removed from the surfaces in contact with soft tissue.

#### **2.4.5.3. Delivery to the clinician**

The structure must be thoroughly cleaned of processing residues and other contaminants using an appropriate method after processing.

The specific abutment screws supplied with the bridges, bars, and implant base structures at the implant level must be delivered to the recipient in the unopened original packaging and clearly identifiable.

#### **2.4.5.4. Insertion**

Bridges, bars, and implant abutments at the implant level are permanently inserted into the implant after cleaning, disinfection, and recommended sterilization (required if the implant was not placed after a complete healing phase, or if it was placed subcrestal) as follows:

1. Remove the implant cover screw or healing cap.
2. Clean the internal configuration of the implant and the implant shoulder.
3. Ensure that the contact surfaces between the implant and the bridge, bar, or implant abutment are absolutely clean and, in particular, that bone debris and soft tissue are completely removed.
4. Insert the bridge, bar, or implant abutment into the implants.
5. Screw in the new, unused abutment screws supplied with the structure. Use the screwdriver for the corresponding implant system for this purpose.
6. Tighten the prosthetic screws in a crosswise pattern using the specified torque; the value should be compared to the corresponding specification in the implant manufacturer's instructions for use and reduced if necessary.

#### **Caution:**

- If the tightening torque exceeds the recommended value, this may result in damage to or destruction of the screw and/or damage to the bridge, the base structure of the bar, or the implant.
  - If the tightening torque is lower than the recommended value, this may result in loosening of the bridge, the bar, or the implant abutment.
1. Ensure that the implant fits properly. It is recommended to take a follow-up X-ray to verify that the bridge, bar, or implant abutment fits correctly into the implant. No soft tissue should be trapped, and there should be no visible gap between the bridge, bar, or implant abutment and the implant.
  2. Retighten to the same torque after at least 5 minutes.
  3. Bridges and implant abutments: Cover the screw head with a fully removable material to protect it. Then, hermetically seal the screw channel with appropriate materials, such as a composite.
  4. Bars: For hygienic reasons, the screw channels must be sealed with appropriate, removable materials. First, cover the screw head with a fully removable material such as cotton or gutta-percha.

## 2.5. Bridges, bars, and implant-supported abutment structures (screw-retained)

### 2.5.1. Product Description

BIOTECH DENTAL DIGITAL abutment-level bridges are custom-made prosthetic restorations designed to replace lost tooth structure and are screwed onto two or more abutments.

BIOTECH DENTAL DIGITAL abutment-level bars serve as custom connecting elements for removable prosthetic restorations and are screwed onto two or more abutments.

BIOTECH DENTAL DIGITAL abutment-level implant frameworks serve as custom-made metal support structures for prosthetic dental restorations made of prosthetic resin and prefabricated prosthetic teeth made of resin or ceramic, and are screwed onto two or more abutments.

The frameworks are screwed occlusally onto prefabricated abutments specifically designed for this application. No intermediate elements requiring bonding are used.

In the case of conical implant-abutment connections, bridges, bars, and implant-level abutments are also supported by the implant shoulder.

Anatomical bridges and abutment-level bridge frameworks can consist of 2 to 16 units. Bars and abutment-level implant frameworks are available for 2 to 10 abutments.

The design of abutment-level bridges, bars, and implant base structures requires a scan to be taken at the abutment level. Scan caps or scan bodies are available for this purpose.

The anatomically reduced bridge framework can be veneered with suitable layering materials or with CAD/CAM-fabricated superstructures.

Thanks to the BIOTECH DENTAL DIGITAL CAD library, bridges and bars can be supplemented with attachments for securing removable restorations, and bars can also be designed with different profiles.

A suitable bar profile is available in the BIOTECH DENTAL DIGITAL CAD library for implant abutments.

The frameworks are non-sterile and intended for single use.

Bridges and implant abutments have a milled surface; bars are also available polished.

Titanium base adapters for bar abutments are available as an option for certain implant systems. With this option, the bridge, bar, or implant abutment is cemented intraorally to the Titanium base adapters screwed onto the abutments or multi-unit abutments. This compensates for any fitting inaccuracies (passive fit) that could lead to stress.

For bridges, bars, and implant abutments, the corresponding prosthetic screws are supplied in a separate package but can also be ordered separately. The structures are secured to the abutments using the torque specified by the manufacturer.

#### Material:

##### Laser-sintered prostheses

|                                 |  |
|---------------------------------|--|
| Titanium alloy (Ti-6Al-4V ELI): | CTE (25 to 500 °C): $\sim 10.16 \times 10^{-6} \text{ K}^{-1}$ |
| Chromium-cobalt alloy (Cr-Co)   | CTE (20 to 600 °C): $\sim 14.4 \times 10^{-6} \text{ K}^{-1}$  |

##### Machined prostheses

|                            |  |
|----------------------------|--|
| Titanium alloy (Ti-6Al-4V) | CTE (20 to 600°C): $\sim 10.3 \times 10^{-6} \text{ K}^{-1}$ |
| CoCr Alloy                 | CTE (20 to 500°C): $\sim 14.0 \times 10^{-6} \text{ K}^{-1}$ |

### 2.5.2. Intended Purpose

Bridges, bars, and implant abutments are destined for functional and/or aesthetic rehabilitation of partially or

completely edentulous patients, as components of dental prostheses in combination with endosseous implant in the maxilla and/or the mandible.

### 2.5.3. Indications

BIOTECH DENTAL DIGITAL bridges, bars, and implant abutments are used for prosthetic reconstruction on the lower or upper jaw to restore partially or completely toothless sections. Specifically for:

- Partially edentulous jaws with multiple missing teeth or
- Edentulous jaws.

### 2.5.4. Warning

- Clinical situation that does not allow compliance with design specifications (see 2.7.5.1 "Design")
- Use on an incompatible abutment
- Support of the joint between the tooth (abutment) and the abutment

### 2.5.5. Application

#### 2.5.5.1. Design

The following requirements must be met when designing bridges, bars, and implant abutments at the abutment level:

- Minimum wall thickness: 0.5 mm
- Minimum cross-sections of connectors:
  - Bridge element/split crowns, anterior tooth: 6 mm<sup>2</sup>
  - Bridge element/split crowns, posterior tooth: 9 mm<sup>2</sup>
  - Custom-designed bars and implant abutments: 8 mm<sup>2</sup>
- Titanium alloy bridges: maximum of 3 connected pontics
- CoCr alloy bridges: maximum of 3 connected pontics in the posterior region and 4 in the anterior region.
- Bridges: maximum length of abutments; width of one molar, but no more than 10 mm.
- Implant bars and abutments: maximum distance between the two implant axes: 30 mm
- Implant bars and abutments: maximum length of extensions: half the length of the adjacent inter-implant segment, but no more than 10 mm.
- Indications and contraindications regarding abutments must be considered in accordance with the manufacturer's instructions.

The lower mechanical load-bearing capacity of implants with reduced diameter must be taken into account. To facilitate the placement of the model's veneered bridges in the sintering furnace, lateral retention pins with a minimum diameter of 1 mm and a maximum length of 5 mm may be designed on the CAD model.

#### 2.5.5.2. Post-processing

If no information on the surface treatment of bridge frameworks or implant abutments for the selected veneering/construction material is available from the manufacturer, CoCr alloy bridge frameworks and implant abutments must be sandblasted with pure aluminum oxide with a grain diameter of 110 µm – 250 µm minimum and a blasting pressure of 2 to 4 bar; bridge frameworks and implant abutments made of titanium alloy with pure aluminum oxide with a grain size of 110 µm – 250 µm and a blasting pressure of 2 to 3 bar.

#### Warning:

During the machining or heat treatment of CoCr material (e.g., milling, grinding), dust and fumes may be generated, which may pose a health risk if inhaled. Therefore, appropriate measures must be taken to reduce exposure to dust and fumes, such as the use of extraction and ventilation systems or, for short-term use, appropriate personal protective equipment for respiratory filtration.

Bridge frameworks and implant abutments must then be rinsed under running water, steam-cleaned, and degreased with ethyl alcohol (for titanium alloys, never use hydrofluoric acid for cleaning). Do not subsequently contaminate the structure by touching it with unprotected fingers.

Bridge frameworks are veneered with the appropriate veneering material for the base material (titanium or CoCr alloy) or with CAD/CAM-fabricated superstructures. The manufacturer's application instructions must be followed. In this context, it should be noted that discoloration of the titanium alloy can occur at temperatures above 300°C.

#### Caution:

- The surface quality and wall thickness must not be altered or reduced in load-bearing areas.
- Heating titanium to temperatures above 800°C causes structural changes and a reduction in strength.
- The material manufacturers' processing instructions must be followed; they are available at
- **Laser-sintered prostheses**

|  |   |
|--|---|
| CoCr alloy (Scheftner Starbond CoS Powder 45): | Please contact us by email at the following address: <a href="mailto:service@scheftner.dental">service@scheftner.dental</a> |
| Titanium alloy (TEKMAT™Ti64-53/20-A):          | Please contact us by email at the following address: <a href="mailto:sds@tekna.com">sds@tekna.com</a>                       |

- **Machined prostheses**

|  |   |
|--|---|
| CoCr alloy (Eisenbacher KERA DISC):          | Please contact us by email at the following address: <a href="mailto:info@eisenbacher.de">info@eisenbacher.de</a> |
| Titanium alloy (Eisenbacher Kera TI 5-Disc): | Please contact us by email at the following address: <a href="mailto:info@eisenbacher.de">info@eisenbacher.de</a> |

To avoid aesthetic issues, CoCr restorations covered with a layering material must have the oxide layer removed from the surfaces in contact with soft tissue.

#### **2.5.5.3. Delivery to the clinician**

The structure must be thoroughly cleaned of processing residues and other contaminants using an appropriate method after processing.

The specific prosthetic screws supplied with bridges, bars, and implant abutment frameworks must be delivered to the recipient in the unopened original packaging and clearly identifiable.

#### **2.5.5.4. Insertion**

Bridges, bars, and implant abutment frameworks are permanently screwed onto the abutment after cleaning, disinfection, and recommended sterilization (required if the abutment was not placed after a complete healing phase, or if it was placed subcrestally) as follows:

11. Remove the protective caps or temporary prostheses.
12. Ensure that the contact surfaces between the implant abutment and the bridge, bar, or implant base structure are perfectly clean and, in particular, that bone debris and soft tissue are completely removed.
13. Place the bridge, bar, or implant abutment on the abutment permanently fixed in the implant.
14. Screw in the new, unused prosthetic screw corresponding to the implant system and implant diameter using the appropriate screwdriver provided by the implant manufacturer.
15. Tighten the prosthetic screws in a crosswise pattern to the specified torque of 15 Ncm.

#### **Warning:**

- If the tightening torque exceeds the recommended value, it may result in damage to or destruction of the screw and/or damage to the abutment and/or the bridge, bar, or implant base structure.
- If the tightening torque is lower than the recommended value, this may cause the bridge, bar, or implant abutment to loosen.

16. Verify that the abutment is securely in place. Soft tissue must not be trapped, and there must be no visible gap between the bridge, bar, or implant base structure and the abutment.

17. Retighten to the same torque after at least 5 minutes.

18. Bridges and implant abutments: cover the screw head with a fully removable material. Then, hermetically seal the screw channel with appropriate materials,

19. such as a composite.

20. Rails: For hygienic reasons, the screw channels must be sealed with appropriate, removable materials. First, cover the screw head with a fully removable material such as cotton or gutta-percha.

## **2.6. Bridges, bars, and implant abutments for adhesive (screw-retained) bases**

### 2.6.1. Product Description

BIOTECH DENTAL DIGITAL bridges for bonding to Titanium base are custom-made artificial prosthetic restorations designed to replace lost tooth structure and are screwed onto two or more implants/abutments.

BIOTECH DENTAL DIGITAL bars for bonding to Titanium base serve as custom-made connecting elements for removable dental prosthetic restorations and are screwed onto two or more implants/abutments.

BIOTECH DENTAL DIGITAL implant base structures for bonding to Titanium base serve as custom-made metal support structures for prosthetic dental restorations made of prosthetic resin and prefabricated dental prostheses made of resin or ceramic, which are screwed onto two or more implants/abutments.

Structures for bonding to Titanium base are screwed onto the implants occlusally.

Anatomical bridges and bridge frameworks on Titanium base can consist of 2 to 16 units. Bars and implant abutments on Titanium base are available for 2 to 10 implants.

Titanium base bridges made of zirconia and plastic (PMMA) are available for individual characterization in various tooth shades. Zirconia is also available in various degrees of translucency and, optionally, with a natural color gradient from dentin to incisal. In cases of hard and/or soft tissue loss, zirconia and plastic bridges can be designed with an optional gingival structure component.

Temporary restorations with a maximum duration of 12 months can be fabricated from tooth-colored plastic (PMMA = polymethyl methacrylate).

The zirconia frameworks are supplied pre-sintered.

The framework of the anatomically reduced bridge can be veneered with suitable layering materials or with CAD/CAM-fabricated superstructures. Implant abutments serve as a metal framework to reinforce dental resin restorations and to support prefabricated prosthetic teeth.

Using the BIOTECH DENTAL DIGITAL CAD library, bridges and bars made of titanium or CoCr alloys can be supplemented with attachments for securing removable restorations, and the bars can also be designed with different profiles. A suitable bar profile is available in the BIOTECH DENTAL DIGITAL CAD library for implant abutments.

The frameworks are non-sterile and intended for single use.

Bridges and implant abutments have a rough, milled surface; bars are also available polished.

Bridges, bars, and implant abutments are delivered with Titanium bases.

#### Material:

##### Laser-sintered prostheses

|                                 |  |
|---------------------------------|--|
| Titanium alloy (Ti-6Al-4V ELI): | CTE (25 to 500 °C): $\sim 10.16 \times 10^{-6} \text{ K}^{-1}$ |
| Chromium-cobalt alloy (Cr-Co)   | CTE (20 to 600 °C): $\sim 14.4 \times 10^{-6} \text{ K}^{-1}$  |

##### Machined prostheses

|  |  |
|--|--|
| Titanium alloy (Ti-6Al-4V)                 | CTE (20 to 600°C): $\sim 10.3 \times 10^{-6} \text{ K}^{-1}$   |
| CoCr Alloy                                 | CTE (20 to 500°C): $\sim 14.0 \times 10^{-6} \text{ K}^{-1}$   |
| IPS e.max ZirCAD (Zirconia) Prime & LT     | CTE (25 to 500°C): $10.5 \times 10^{-6} \text{ K}^{-1}$  |
| IPS e.max ZirCAD (Zirconia) Prime Esthetic | CTE (25 to 500°C): $10.4 \times 10^{-6} \text{ K}^{-1}$  |
| Polymethyl methacrylate (PMMA)             | CTE: N/A   |
| IPS e.max CAD (lithium disilicate)         | CTE (100 to 400°C): $10.2 \times 10^{-6} \text{ K}^{-1}$<br>CTE (100 to 500°C): $10.5 \times 10^{-6} \text{ K}^{-1}$ |

### 2.6.2. Intended Purpose

Bridges, bars, and implant abutments are destined for functional and/or aesthetic rehabilitation of partially or completely edentulous patients, as components of dental prostheses in combination with endosseous implant in the maxilla and/or the mandible.

### 2.6.3. Indications

BIOTECH DENTAL DIGITAL Bridges, bars, and implant abutments for Titanium base are used for prosthetic reconstruction on the lower or upper jaw to restore partially or completely toothless sections. Specifically, for:

- Partially edentulous jaws with multiple missing teeth or
- Edentulous jaws.

### 2.6.4. Warning

- Clinical situation that does not allow compliance with design specifications (see 2.6.5.1 "Design")
- Use with an incompatible bonding base.
- Support of the joint between the tooth (abutment) and the implant or abutment
- Titanium alloy restorations:
  - Bridge restorations with more than 3 connected pontics
- CoCr alloy restorations:
  - Bridge restorations with more than 3 connected pontics in the posterior region and 4 in the anterior region
- Zirconia restorations:
  - Bridge constructions with more than 2 connected pontics
- PMMA restorations:
  - 12-month usage period
  - Bridge constructions with more than 2 connected pontics

### 2.6.5. Application

#### 2.6.5.1. Design

To facilitate the placement of veneered bridges from the model into the sintering furnace, lateral retention pins with a minimum diameter of 1 mm and a maximum length of 5 mm can be designed on the CAD model. The following requirements must be met when designing bridges, bars, and implant abutments for adhesive bases:

- Minimum wall thickness: 0.5 mm
- Minimum cross-sections of connectors:
  - Bridge element/split crowns, anterior tooth: 6 mm<sup>2</sup>
  - Bridge element/split crowns, posterior tooth: 9 mm<sup>2</sup>
  - Custom-designed bars and implant abutments: 8 mm<sup>2</sup>
- Titanium alloy bridges: maximum of 3 connected pontics; in cases of single support on 2 implants, the maximum distance between implant axes is 30 mm.
- CoCr alloy bridges: maximum of 3 connected pontics in the lateral region and 4 in the anterior region; with single support on 2 implants, the maximum distance between implant axes is 30 mm.
- Zirconia or resin bridges: maximum of 2 connected pontics
- Bridges: maximum length of abutments; width of one molar, but no more than 10 mm.
- Implant bars and abutments: maximum distance between the two implant axes: 30 mm
- Implant bars and abutments: maximum length of extensions: half the length of the adjacent inter-implant segment, but no more than 10 mm.
- Indications and contraindications regarding bonding bases must be considered in accordance with the manufacturer's instructions.

The lower mechanical load-bearing capacity of implants with a reduced diameter must be taken into account. To facilitate the placement of model-based veneered bridges in the sintering furnace ( ), lateral retention pins with a minimum diameter of 1 mm and a maximum length of 5 mm may be designed on the CAD model.

#### 2.6.5.2. Post-processing

If no information on the surface treatment of bridge frameworks or implant abutments for the selected veneering/construction material is available from the manufacturer, CoCr alloy bridge frameworks and implant abutments must be sandblasted with pure aluminum oxide with a grain size of at least 110 µm and a blasting pressure of 3 to 4 bar, Bridge frameworks and implant abutments made of titanium alloy must be sandblasted

with pure aluminum oxide with a grain size of 180 µm and a sandblasting pressure of 2 to 3 bar.

**Warning:**

During the machining or heat treatment of CoCr material (e.g., milling, grinding), dust and fumes may be generated and pose a health risk if inhaled. Therefore, appropriate measures must be taken to reduce exposure to dust and fumes, such as the use of extraction and ventilation systems or, for short-term use, appropriate personal protective equipment for respiratory filtration.

Bridge frameworks or implant abutments must then be rinsed under running water, steam-cleaned, and degreased with ethyl alcohol. Do not contaminate the structure afterward, for example by touching it with unprotected fingers.

Bridge frameworks are veneered with a coating material suitable for the base material (titanium alloy, CoCr alloy, or zirconia) or with CAD/CAM-fabricated superstructures. The manufacturer’s processing instructions must be followed. In this context, it should be noted that discoloration of the titanium alloy can occur at temperatures above 300°C.

PMMA bridges can be modified to optimize the emergence profile by applying a suitable composite or

- Taking into account the minimum permitted wall thicknesses
- by grinding them.

**Caution:**

- The surface quality and wall thickness must not be altered or reduced in load-bearing areas.
- Heating titanium to temperatures above 800 °C causes structural changes and a reduction in strength.

The material manufacturers’ processing instructions must be followed; they are available at

○ **Laser-sintered prostheses**

|  |   |
|--|---|
| CoCr alloy (Scheftner Starbond CoS Powder 45): | Please contact us by email at the following address: <a href="mailto:service@scheftner.dental">service@scheftner.dental</a> |
| Titanium alloy (TEKMAT™Ti64-53/20-A):          | Please contact us by email at the following address: <a href="mailto:sds@tekna.com">sds@tekna.com</a>                       |

○ **Machined prostheses**

|  |   |
|--|---|
| CoCr alloy (Eisenbacher KERA DISC):          | Please contact us by email at the following address: <a href="mailto:info@eisenbacher.de">info@eisenbacher.de</a> |
| Titanium alloy (Eisenbacher Kera TI 5-Disc): | Please contact us by email at the following address: <a href="mailto:info@eisenbacher.de">info@eisenbacher.de</a> |
| Zirconia:                                    | <a href="http://www.ivoclar.com/en_gb/eifu">www.ivoclar.com/en_gb/eifu</a>  |
| PMMA (polymethyl methacrylate):              | <a href="https://www.white-peaks-dental.com/en/downloads/">https://www.white-peaks-dental.com/en/downloads/</a>   |
| Lithium disilicate                           | <a href="http://www.ivoclar.com/en_gb/eifu">www.ivoclar.com/en_gb/eifu</a>  |

To avoid aesthetic issues, CoCr restorations covered with a layering material must have the oxide layer removed from the surfaces in contact with soft tissues.

**2.6.5.3. Return to the clinician**

The processed framework must be thoroughly cleaned of processing residues and other contaminants using an appropriate method.

The following items must be sent to the recipient along with the bridge, bar, or implant abutment, in the unopened original packaging and clearly labeled:

**2.6.5.4. Insertion**

Before the final fixation of the structure, the manufacturer’s instructions in the usage information for the bonding bases and the bonding composite must be followed for the intraoral bonding of titanium bonding bases or titanium caps in bridges, bars, or implant abutments.

### **Appropriate surface conditioners and luting composites**

To connect the bridge, bar, or implant abutment structure to the Titanium bases, we recommend using products that are suitable and validated for this purpose. Follow the manufacturer's instructions. Biotech Dental Digital SAS assumes no liability for the durability of the bond.

Bridges, bars, and implant abutments on Titanium base are permanently screwed onto the implants, bar abutments, or multi-unit abutments after cleaning, disinfection, and sterilization (recommended if the bridge is not placed after a complete healing phase, or if it is placed subcrestally), as follows:

1. At the implant level: remove the implant cover screw or the implant abutment and clean the internal configuration and the implant shoulder.
2. At the abutments: remove the protective caps or temporary prostheses.
3. Ensure that the contact surfaces between the implant abutment, bar abutment, or multi-unit abutment and the bridge, bar, or implant base structure are absolutely clean and, in particular, that bone debris and soft tissue are completely removed.
4. Place the bridge, bar, or implant base structure onto the implants, or permanently secure the bar abutments or multi-unit abutments to the implants.
5. Screw in the new abutment or the unused prosthetic screw using the appropriate screwdriver.
6. Tighten the screws in a crosswise pattern to the specified torque.

### **Warning:**

- If the tightening torque exceeds the recommended value, this may result in damage to or destruction of the screw and/or damage to the implant or abutment and/or the bridge, bar, or implant base structure.
  - If the torque is lower than the recommended value, this may cause the bridge, bar, or implant base structure to loosen.
7. Verify that the implants, bar abutments, or multi-unit abutments are securely in place. No soft tissue should be trapped, and there should be no visible gap between the base structure of the bridge, bar, or implant and the implant, bar abutment, or multi-unit abutment.
  8. Retighten to the same torque after at least 5 minutes.
  9. Bridges and implant abutments: Cover the screw head with a fully removable material to protect it. Then, hermetically seal the screw channel with appropriate materials, such as a composite.
  10. Bars: For hygienic reasons, the screw channels must be sealed with appropriate, removable materials. First, cover the screw head with a fully removable material such as cotton or gutta-percha.

## 2.7. Crowns and bridges (bonded)

### 2.7.1. Product description

BIOTECH DENTAL DIGITAL crowns, bridges, and crown/bridge frameworks are cemented onto tooth stumps or custom abutments. They are designed virtually using specialized CAD software and, if compatible with the software used, the BIOTECH DENTAL DIGITAL CAD material library.

Depending on the material selected, bridges or bridge frameworks can consist of 2 to 16 units and can be constructed with up to 3 connected intermediate units (4 in the anterior region when using a CoCr alloy and 2 when using zirconia). Appropriate connection cross-sections must be taken into account.

BIOTECH DENTAL DIGITAL crowns, bridges, and crown-bridge frameworks are available for the selected materials for individual characterization in various tooth shades and degrees of translucency, depending on the indication and material. Zirconia is also available with a natural shade gradient from the dentin to the incisal edge.

Zirconia crowns, bridges, and crown and bridge frameworks are delivered pre-sintered.

The anatomically reduced crown or bridge framework can be veneered with suitable layering materials or with CAD/CAM-fabricated superstructures. In cases of hard and/or soft tissue loss, bridges and bridge frameworks can be designed with an optional gingival structure.

#### Material:

#### Laser-sintered prostheses

|                                 |  |
|---------------------------------|--|
| Titanium alloy (Ti-6Al-4V ELI): | CTE (25 to 500 °C): $\sim 10.16 \times 10^{-6} \text{ K}^{-1}$ |
| Chromium-cobalt alloy (Cr-Co)   | CTE (20 to 600 °C): $\sim 14.4 \times 10^{-6} \text{ K}^{-1}$  |

#### Machined prostheses

|  |  |
|--|--|
| Titanium alloy (Ti-6Al-4V)                 | CTE (20 to 600°C): $\sim 10.3 \times 10^{-6}/\text{K}^{-1}$  |
| CoCr Alloy                                 | CTE (20 to 500°C): $\sim 14.0 \times 10^{-6}/\text{K}^{-1}$  |
| IPS e.max ZirCAD (Zirconia) Prime & LT     | CTE (25 to 500°C): $10.5 \times 10^{-6}/\text{K}$  |
| IPS e.max ZirCAD (Zirconia) Prime Esthetic | CTE (25 to 500°C): $10.4 \times 10^{-6} \text{ K}^{-1}$  |
| Polymethyl methacrylate (PMMA)             | CTE: N/A   |
| IPS e.max CAD (Lithium disilicate)         | CTE (100 to 400°C): $10.2 \times 10^{-6}\text{K}^{-1}$<br>CTE (100 to 500°C): $10.5 \times 10^{-6}\text{K}^{-1}$ |

### 2.7.2. Intended Purpose

BIOTECH DENTAL DIGITAL bridges are destined for functional and/or aesthetic rehabilitation of partially or completely edentulous patients, as components of dental prostheses in combination with endosseous implant in the maxilla and/or the mandible.

BIOTECH DENTAL DIGITAL crowns are intended for functional and/or aesthetic rehabilitation of partially edentulous patients, as components of dental prostheses in the maxilla and/or the mandible.

### 2.7.3. Indications

BIOTECH DENTAL DIGITAL bridges are used for prosthetic reconstruction on the lower or upper jaw to restore partially or completely toothless sections.

BIOTECH DENTAL DIGITAL crowns are used for the replacement of partially edentulous jaws, dental defects or compromised teeth with a lack of tooth substance.

Zirconia restorations have different indications depending on their translucency.

#### **2.7.4. Contraindications**

- Clinical situation that does not allow compliance with design specifications (see 2.9.5.1 "Design")
  - Titanium alloy restorations:
    - Bridge restorations with more than 3 connected pontics
  - CoCr alloy restorations:
    - Bridge restorations with more than 3 connected pontics in the posterior region and 4 in the anterior region
  - Zirconia restorations:
    - Bridge restorations with more than one intermediate element or more than 3 elements
    - Temporary placement

#### **Warning:**

- The use of different alloy groups in the patient's mouth can cause galvanic effects.

#### **2.7.5. Application**

##### **2.7.5.1. Design**

The BIOTECH DENTAL DIGITAL CAD material libraries contain detailed indications and material-specific design data.

In particular, the specified minimum wall thickness must be taken into account during design. Any modification of the cement thickness specified in the material libraries or in the parameter table may result in a poor fit of the restoration.

To achieve optimal results in dental restorations, it is important to follow both general and material-specific preparation guidelines:

- For restorations on stumps: remove sufficient hard tooth structure.
- Blocking undercuts
- Preparation of grooves
- Avoid abrupt transitions
- Observe material-specific transitions

To facilitate the placement of veneered wood bridges and crowns from the model into the sintering furnace, lateral retaining pins with a minimum diameter of 1 mm and a maximum length of 5 mm can be designed using CAD.

##### **2.7.5.2. Post-processing**

If no information regarding the surface condition of bridge and crown frameworks for the selected veneering material is available from the respective manufacturer, CoCr alloy frameworks must be sandblasted with pure aluminum oxide with a grain size of at least 110 µm and a blasting pressure of 3 to 4 bars, and titanium alloy frameworks with pure aluminum oxide with a grain size of 180 µm and a blasting pressure of 2 to 3 bar.

#### **Warning:**

During the machining or heat treatment of CoCr material (e.g., milling, grinding), dust and fumes may be generated and pose a health risk if inhaled. Therefore, appropriate measures must be taken to reduce exposure to dust and fumes, such as the use of extraction and ventilation systems or, for short-term use, appropriate personal protective equipment for respiratory filtration.

The frameworks must then be rinsed under running water, cleaned with a steam jet, and degreased with ethyl alcohol. Do not subsequently contaminate the structure by touching it with unprotected fingers.

The frameworks are coated with a coating material suitable for the material used (titanium, CoCr alloy, or zirconia) or with CAD/CAM-fabricated superstructures. The manufacturer's processing instructions must be followed. In this context, it should be noted that discoloration of the titanium alloy can occur at temperatures above 300°C.

#### **Caution:**

- The surface quality and wall thickness must not be altered or reduced in load-bearing areas.
- Heating titanium to temperatures above 800°C causes structural changes and a reduction in strength. The processing instructions from the manufacturers of the following materials must be followed; they are available at:
  - **Laser-sintered prostheses**

|  |   |
|--|---|
| CoCr alloy (Scheftner Starbond CoS Powder 45): | Please contact us by email at the following address: <a href="mailto:service@scheftner.dental">service@scheftner.dental</a> |
| Titanium alloy (TEKMAT™Ti64-53/20-A):          | Please contact us via email at the following address: <a href="mailto:sds@tekna.com">sds@tekna.com</a>                      |

○ **Machined prostheses**

|  |   |
|--|---|
| CoCr alloy (Eisenbacher KERA DISC):          | Please contact us by email at the following address: <a href="mailto:info@eisenbacher.de">info@eisenbacher.de</a> |
| Titanium alloy (Eisenbacher Kera TI 5-Disc): | Please contact us by email at the following address: <a href="mailto:info@eisenbacher.de">info@eisenbacher.de</a> |
| Zirconia:                                    | <a href="http://www.ivoclar.com/en_gb/eifu">www.ivoclar.com/en_gb/eifu</a>  |
| PMMA (polymethyl methacrylate):              | <a href="https://www.white-peaks-dental.com/en/downloads/">https://www.white-peaks-dental.com/en/downloads/</a>   |

To avoid aesthetic issues, CoCr restorations covered with a layering material must have the oxide layer removed from the surfaces in contact with soft tissues.

**2.7.5.3. Return to the clinician**

After processing, the structures must be thoroughly cleaned of processing residues and other contaminants using an appropriate method.

**2.8. Inlays, onlays, and partial crowns**

**2.8.1. Product description**

BIOTECH DENTAL DIGITAL inlays, onlays, and partial crowns are bonded to the teeth and serve to replace tooth structure. They are available in titanium and cobalt-chrome.

**Material:**

|                                 |  |
|---------------------------------|--|
| Titanium alloy (Ti-6Al-4V ELI): | CTE (25 to 500 °C): $\sim 10.16 \times 10^{-6} \text{ K}^{-1}$ |
| Chromium-cobalt alloy (Cr-Co)   | CTE (20 to 600 °C): $\sim 14.4 \times 10^{-6} \text{ K}^{-1}$  |

**2.8.2. Intended purpose**

BIOTECH DENTAL DIGITAL inlays, onlays, and partial crowns are intended for functional and/or aesthetic rehabilitation of partially edentulous patients, as components of dental prostheses in the maxilla and/or the mandible.

**2.8.3. Indications**

BIOTECH DENTAL DIGITAL inlays, onlays, and partial crowns are used for the replacement of partially edentulous jaws, dental defects or compromised teeth with a lack of tooth substance.

**2.8.4. Warning**

- Very deep subgingival preparations
- Markedly reduced residual dentition
- Temporary placement or trial wear

**2.8.5. Application**

**2.8.5.1. Design**

To achieve optimal results for dental restorations, the general and material-specific preparation guidelines must be followed:

- Removal of a sufficient amount of hard dental tissue
- Elimination of undercuts
- Preparation of grooves
- Avoid abrupt transitions

- Observe material-specific transitions

The BIOTECH DENTAL DIGITAL CAD material libraries contain detailed instructions and material-specific design data.

### 2.8.5.2. *Post-processing*

The material manufacturer’s processing instructions must be followed; they are available at

|  |   |
|--|---|
| CoCr alloy (Scheftner Starbond CoS Powder 45): | Please contact us by email at the following address <a href="mailto:service@scheftner.dental">HYPERLINK</a><br>"mailto:service@scheftner.dental":<br>service@scheftner.dental |
| Titanium alloy (TEKMAT™Ti64-53/20-A):          | Please contact us by email at the following address: <a href="mailto:sds@tekna.com">sds@tekna.com</a>   |

### 2.8.5.3. *Delivery to the clinician*

After processing, the abutments must be thoroughly cleaned of processing residues and other contaminants using an appropriate method.

## 3. Products compatible with BIOTECH DENTAL DIGITAL implant prostheses

### Note:

Compatibility tables for BIOTECH DENTAL DIGITAL implant prostheses, as well as installation guides for the Biotech Dental Digital libraries related to BDD products, are available for download in PDF format on the ordering website. You may contact customer service to request a hard copy, which will be delivered with your order or subsequently within a maximum of seven calendar days following receipt of your request, at no additional cost.

## 4. Clinical Benefits

### 4.1. Implant Prosthetics

Implant prosthetic components are used for dental prosthetic treatments, which allow for the replacement of missing teeth and/or the restoration of masticatory function and the corresponding improvement in quality of life.

### 4.2. Dental Restorations

Dental restorations enable the oral rehabilitation of patients with dental defects, and/or the restoration of masticatory function and the corresponding improvement in quality of life.

## 5. Contraindications

BIOTECH DENTAL DIGITAL custom-made medical devices must not be used in patients with contraindications to dental restoration, including but not limited to:

- Allergy to any of the materials used in the manufacture of the custom dental restoration and identified in the current Instructions for Use (IFU).
- Non-adult patients
- Bone and/or oral soft tissue pathologies or infections
- Pregnancy and breastfeeding
- Poor oral hygiene
- Bruxism, parafunctional habits, malocclusion, and/or temporomandibular joint disorders

In addition to the contraindications listed above, the specific contraindications for the dental implant system used must also be observed.

## 6. Adverse side effects

The use of BIOTECH DENTAL DIGITAL products is part of the dental restoration procedure, which may be associated with typical side effects such as:

- Pain, swelling, mucogingival inflammation
- Difficulty speaking
- Allergic reaction
- Inhalation or ingestion
- Negative impact on adjacent teeth
- Fractures of the restoration
- Aesthetic issues

## **7. Intended Purposer and patient population**

### **7.1. Intended User**

The products are intended for use by professionals (dentists, dental/maxillofacial surgeons, dental technicians) only.

### **7.2. Target patient population**

The products are intended for use in skeletally mature patients undergoing dental prosthetic treatment. As custom-made devices, they are intended for the exclusive use of a specific, predefined patient to meet their individual conditions and needs. Please ensure that contraindications and precautions are taken into account.

## **8. Preparation**

### **8.1. Handling of Packaging**

The packaging must be inspected prior to opening to ensure it is undamaged and does not have an expired date, and should only be opened immediately before using the products.

If the original packaging is damaged, the items cannot be returned to the manufacturer. Prosthetic components are shipped non-sterile, unless explicitly marked as sterile.

**Warning:** The products must not be used if they cannot be identified.

**Caution:** The products must not be used if the packaging has been damaged or previously opened.

### **8.2. Handling of Non-Sterile Products and Reprocessing**

BIOTECH DENTAL DIGITAL products are delivered non-sterile and must be cleaned and disinfected before and after each use on a patient.

**Warning:**

- BIOTECH DENTAL DIGITAL single-use products must not be reused on the patient, as the safety of the preparation and/or functional safety cannot be guaranteed.

**Caution:**

- All non-sterile packaged products must not be sterilized in their original packaging.
- Plastic components should not be sterilized (only cleaned and disinfected), with the exception of BIOTECH DENTAL DIGITAL products made of PMMA. These products can be sterilized.

### **8.3. Cleaning, Disinfection, and Sterilization**

#### **8.3.1. General notes**

Thorough cleaning and disinfection are essential for effective sterilization.

To ensure effective cleaning, remove detachable accessories before beginning the cleaning process.

Always comply with applicable legal and national hygiene regulations for dental/medical practices and hospitals. This applies in particular to guidelines regarding the effective inactivation of prions.

Note: Any deviation from the procedures described here must be carefully evaluated by the person responsible for the treatment to determine its effectiveness and any potential negative consequences. The operator is responsible for ensuring that the actual reprocessing produces the desired results with the equipment, materials, and personnel used in the reprocessing facility. This requires validations and routine

monitoring.

### **For sterilizers and disinfectors:**

As part of your responsibility regarding the sterility of the products used, please observe the following points:

- In general, use only adequately validated methods specific to the equipment and product for cleaning/disinfection and sterilization,
- Regularly check and maintain the equipment used (thermal disinfectant, sterilizer) and
- Ensure that the validated parameters are followed during each cycle.

### **Cleaning and Disinfection Device**

When selecting a cleaning and disinfection device, ensure that:

- The cleaning and disinfection device has proven efficacy (CE marking, compliant with ISO 15883) and is specifically validated for the equipment and product,
- A validated thermal disinfection program is used (A0 value > 3000 or—for older machines—at least 5 minutes at 90°C/194°F),
- The program used is suitable for the products and includes sufficient rinse cycles (it is recommended to perform at least three drain rinse cycles after the cleaning step or to monitor the rinse cycles by measuring conductivity),
- Only deionized water with a low bacterial count (maximum 10 bacteria/ml) and endotoxin content (maximum 0.25 endotoxin units/ml) is used for rinsing,
- The air used for drying is filtered, and the cleaning and disinfection device is regularly inspected and maintained. Use only thermal disinfection that does not require a disinfectant.

### **Cleaning products for mechanical cleaning and disinfection**

Care should be taken to ensure that

- They are generally suitable for metal, ceramic, and plastic parts,
- The cleaning products used are compatible with the devices

### **Cleaning and disinfection agents for manual cleaning and disinfection**

Care should be taken to ensure that

- They are suitable for cleaning and disinfecting dental prostheses,
- They are compatible with the materials of the items to be cleaned and disinfected,
- They are suitable for ultrasonic cleaning (no foaming),
- The disinfectant has proven efficacy for disinfection (e.g., VAH/DGHM list, CE marking),
- The cleaning agents and disinfectants are compatible with each other.

The concentrations, application times, and instructions for subsequent rinsing specified by the manufacturer must be followed.

Instead of water, suitable cleaning agents or suitable cleaning agents with a disinfecting effect (to perform initial disinfection for staff protection) may be used for point-of-use preparation (sections 7.3.2.1 and 7.3.3.1). These agents must be formulated and suitable for the pre-cleaning of soiled prosthetic components (no protein-binding effect) and must be compatible with the agents used subsequently, otherwise the effectiveness of the subsequent steps may be impaired.

### **Sterilization**

Only steam sterilization methods with fractionated vacuum (forced ventilation) are permitted. Care must be taken to ensure that:

- The maximum steam sterilization temperature does not exceed 138°C/280°F (PMMA products: 121°C/250°F),
- The steam sterilizers used bear the CE marking and comply with the requirements of standard EN 13060 or EN 285,
- Only appropriately validated methods specific to the equipment and product, in accordance with ISO 17665, are used for sterilization,
- The sterilizers used are regularly inspected and maintained, and
- Ensure that the validated parameters are adhered to in every cycle.

### **Warning:**

- Products made of PMMA may only be sterilized using the steam sterilization method at 121°C/250°F

- with a holding time of 20 minutes.
- All BIOTECH DENTAL DIGITAL products must not be sterilized in their original packaging.
- Plastic parts must not be sterilized.

Biotech Dental Digital recommends using disposable sterilization pouches that bear the CE mark. Ensure that the sterilization packaging is suitable for steam sterilization (i.e., a constant temperature of at least 142°C/288°F and sufficient steam permeability) and that the products are adequately protected against mechanical damage.

Two procedures are described below for cleaning and disinfection prior to sterilization: a mechanical procedure in Section 7.3.2 and a manual procedure in Section 7.3.3. Whenever possible, the mechanical procedure described in Section 7.3.2 should be used.

### **8.3.2. Mechanical cleaning and disinfection**

#### **8.3.2.1. Preparation at the point of use**

Immediately after use in the patient's mouth, remove coarse contaminants and place the device in a water bath.

#### **8.3.2.2. Preliminary cleaning (within 2 hours of use)**

(After use in the patient's mouth or contamination in the dental laboratory)

1. Remove contaminants from the prosthetic components under running water using a soft brush (no metal bristles or steel wool). The brush must be used exclusively for this purpose. Brush until no visible contaminants remain. Some prosthetic components have cavities. Residues may be present there. Contaminants and residues must no longer be visible.
2. Rinse all cavities at least 3 times with 20 ml of water.

#### **8.3.2.3. Transport of Pre-cleaned Parts**

Transport the pre-cleaned parts in a contamination-proof container.

#### **8.3.2.4. Cleaning and disinfection of pre-cleaned items in the cleaning and disinfection device**

#### **Work steps**

1. Use a thermal disinfectant.
2. Use a cleaning agent.
3. Place the pre-treated prosthetic components in the thermal disinfectant using a small parts basket. The prosthetic components must not come into contact with each other.
4. Start the program.
5. Remove the prosthetic components from the thermal disinfectant after the program has finished.

#### **8.3.2.5. Drying**

Dry the prosthetic components. Use a lint-free disposable cloth and dry compressed air that is oil-free and low in bacteria. We also recommend using a sterile filter.

#### **8.3.2.6. Inspection**

Inspect the prosthetic components for any signs of corrosion, surface damage, chipping, or contamination. Damaged prosthetic components must not be used. Clean and disinfect any prosthetic components that are still contaminated.

#### **8.3.2.7. Packaging**

Pack the prosthetic components promptly for sterilization. We recommend using disposable sterilization packaging bearing the CE mark. Ensure that the sterilization packaging is suitable for steam sterilization (i.e., a constant temperature of at least 141°C/286°F, sufficient steam permeability) and that the products are adequately protected against mechanical damage.

### **8.3.3. Manual cleaning and disinfection**

#### **8.3.3.1. Preparation at the point of use**

Immediately after use in the patient's mouth, remove gross contamination and place the devices in a water bath.

#### **8.3.3.2. Preliminary cleaning (within 2 hours of use)**

1. Remove contaminants from the prosthetic components under running water using a soft brush (no metal bristles or steel wool). The brush must be used exclusively for this purpose. Brush until no visible

contaminants remain. Some prosthetic components have cavities. Remove all residues from these cavities. Contaminants should no longer be visible.

2. Rinse all cavities at least 3 times with 20 ml of water.

#### **8.3.3.3. Transport of Pre-cleaned Parts**

Transport the pre-cleaned parts in a contamination-proof container.

#### **8.3.3.4. Cleaning**

**Caution:** Hard-to-reach areas, such as recesses and cavities, require special attention during cleaning.

##### **Work steps:**

1. All parts must be placed fully submerged in the cleaning solution bath for the prescribed exposure time. Parts must not come into contact with one another during the exposure time.
2. Brush the parts at the start of the exposure time for at least 30 seconds each (do not use metal bristles or steel wool). Also brush all accessible cavities for at least 30 seconds each.
3. Rinse all cavities at least 3 times using a syringe containing at least 20 ml of cleaning agent at the start of the exposure time.
4. At the end of the exposure time, brush all prosthetic components for 30 seconds under running water, and also brush all accessible cavities for 30 seconds each.

#### **8.3.3.5. Rinsing After Cleaning**

1. Thoroughly rinse all parts under running water for at least one minute each. Also rinse areas with limited access, such as the cavities of devices.
2. Rinse all cavities three times using a syringe containing at least 20 ml of water per cavity.

#### **8.3.3.6. Disinfection**

1. Place all parts fully submerged in the disinfectant bath for the exposure time specified by the disinfectant manufacturer. Parts must not come into contact with one another during the exposure time.
2. Rinse all cavities at least once using a syringe containing at least 20 ml of disinfectant at the start of the exposure time.

#### **8.3.3.7. Rinsing after Disinfection**

Deionized water with a low bacterial count (max. 10 bacteria/ml) and endotoxin content (max. 0.25 endotoxin units/ml) must be used for these rinses.

1. Rinse all parts thoroughly under running water for at least one minute each. Also rinse areas with limited access, such as cavities.
2. Rinse all cavities 3 times with a syringe containing at least 20 mL of water per cavity.

#### **8.3.3.8. Drying**

Dry the prosthetic components. Use a lint-free disposable cloth and dry compressed air that is oil-free and low in bacteria. We also recommend using a sterile filter.

#### **8.3.3.9. Inspection**

Inspect the prosthetic components for any signs of corrosion, surface damage, chipping, or contamination. Damaged prosthetic components must not be used. Clean and disinfect any prosthetic components that are still contaminated.

#### **8.3.3.10. Packaging**

Pack the prosthetic components promptly for sterilization. We recommend using disposable sterilization packaging bearing the CE mark. Ensure that the sterilization packaging is suitable for steam sterilization (i.e., a constant temperature of at least 141°C/286°F, sufficient steam permeability) and that the products are adequately protected against mechanical damage.

#### **8.3.3.11. Sterilization**

The following parameters apply: (fractional vacuum/forced evacuation only): at least 4 minutes at 134°C/273°F.

For devices made of PMMA: 20 minutes at 121°C/250°F.

##### **Warning:**

- Products made of PMMA may only be sterilized using the steam sterilization method at 121°C/250°F with a holding time of 20 minutes.
- All BIOTECH DENTAL DIGITAL products must not be sterilized in their original packaging.
- Plastic parts must not be sterilized.

National guidelines must be followed.

Biotech Dental Digital recommends a minimum drying time of 30 minutes for each of the cycles described above

## 9. Material

| Products  | Material  |                     |        |             |
|---|---|---------------------|--------|-------------|
| Machining<br>Monoblock<br>titanium<br>abutments   | Ti-6Al-4V ELI titanium alloy according to ASTM F136 & DIN EN ISO 5832-3;<br>Chemical composition (% by weight): |                     |        |             |
|   | Al  | 5.5–6.5             | C      | ≤ 0.4       |
|   | V   | 3.5–4.5             | N      | ≤ 0.4       |
|   | Fe  | <0.25               | H      | ≤ 0.4       |
|   | O   | ≤ 0.4               | T      | Equilibrium |
| Sintering<br>Abutments, inlays,<br>onlays, partial<br>crowns, crowns,<br>bridges, bars,<br>titanium implant<br>abutments                | Ti-6Al-4V ELI titanium alloy according to ASTM F136 & DIN EN ISO 5832-3;<br>Chemical composition (% by weight): |                     |        |             |
|   | Al  | 6.0–6.5             | C      | ≤ 0.08      |
|   | V   | 3.5–4.5             | N      | ≤ 0.03      |
|   | Fe  | ≤ 0.25              | H      | ≤ 0.012     |
|   | O   | <0.10               | T      | Equilibrium |
| Y   | ≤ 0.005   | Other<br>components | ≤ 0.10 |             |
| Machining - One-<br>piece cobalt-<br>chromium<br>abutments*   | Type 4 CoCr alloy according to ISO 22674;<br>Chemical composition (% by weight):                                |                     |        |             |
|   | Co  | Balance             |        |             |
|   | Cr  | 26.0–30.0           |        |             |
|   | W   | 9.5                 |        |             |
|   | Mo  | 5.0–7.0             |        |             |
|   | Fe  | ≤0.75               |        |             |
|   | Ni  | ≤0.1                |        |             |
| Other components: C, Si, N, Mn  | <2.5  |                     |        |             |
| Sintering<br>Abutments, inlays,<br>onlays, partial<br>crowns, crowns,<br>bridges, bars,<br>implant base<br>structures<br>Chrome-Cobalt* | Type 4 CoCr alloy according to ISO 22674;<br>Chemical composition (% by weight):                                |                     |        |             |
|   | Co  | 59.                 |        |             |
|   | Cr  | 25                  |        |             |
|   | W   | 9.5                 |        |             |
|   | MB  | 3.5                 |        |             |
|   | Si  | 1.0                 |        |             |
|   | Other components: C, N, Fe, Mn  | <1.0                |        |             |
| Bridges, bars,<br>titanium implant<br>abutments   | Ti-6Al-4V ELI Grade 5 titanium alloy per ASTM F136 & DIN EN ISO 5832-3;<br>Chemical composition (% by weight):  |                     |        |             |
|   | Al  | 6.2                 | C      | ≤ 0.4       |
|   | V   | 4                   | N      | ≤ 0.4       |
|   | Fe  | ≤ 0.4               | H      | ≤ 0.4       |
|   | O   | ≤ 0.4               | T      | 89.4        |
| Bridges, bars,<br>and implant<br>abutments made<br>of cobalt-<br>chromium*  | Type 4 CoCr alloy according to ISO 22674;<br>Chemical composition (% by weight):                                |                     |        |             |
|   | Co  | 59.0                |        |             |
|   | Cr  | 25.0                |        |             |
|   | W   | 9.5                 |        |             |
|   | Si  | 1                   |        |             |

|  |   |  |
|--|---|--|
|  | Other components: C, Fe, Mn, N  | <1.0   |
| Mesostructures for titanium-zirconia abutments and bridges | Pre-sintered zirconia, Type II/Class 3 ISO 6872;<br>Chemical composition (% by weight):   |  |
|  | <b>Material / Product</b>   | <b>IPS e.max Zirconium oxide (ZrO<sub>2</sub>)</b> |
|  |   | <b>IPS e.max Zirconium oxide (ZrO<sub>2</sub>)</b> |
|  |   | <b>IPS e.max Zirconium oxide (ZrO<sub>2</sub>)</b> |
|  |   | <b>IPS e.max Zirconium oxide (ZrO<sub>2</sub>)</b> |
|  |   | <b>IPS e.max Zirconium oxide (ZrO<sub>2</sub>)</b> |
|  |   | <b>IPS e.max Zirconium oxide (ZrO<sub>2</sub>)</b> |
| PMMA crowns and bridges                                    | Polymethyl methacrylate (PMMA); Chemical composition (% by weight):   |  |
|  | PMMA  | equilibrium  |
|  | TiO <sub>2</sub>  | ≤ 0.50   |
|  | Other pigments  | ≤ 0.35   |
| <b>Note:</b> Maximum shelf life of 12 months.              |   |  |
| Mesostructures and titanium base crowns in e.max           | SiO <sub>2</sub>  | 57–80%   |
|  | Li <sub>2</sub> O   | 11–19%   |
|  | K <sub>2</sub> O  | 0–13%  |
|  | P <sub>2</sub> O <sub>5</sub>   | 0–11%  |
|  | ZrO <sub>2</sub>  | 0–8%   |
|  | ZnO   | 0–8%   |
|  | Al <sub>2</sub> O <sub>3</sub>  | 0–5%   |
|  | MgO   | 0–5%   |
|  | Coloring oxides (CeO <sub>2</sub> , MnO <sub>2</sub> , V <sub>2</sub> O <sub>5</sub> , Tb <sub>4</sub> O <sub>7</sub> , Er <sub>2</sub> O <sub>3</sub> , La <sub>2</sub> O <sub>3</sub> ) | 0–8%   |

**Table 2: Material data**

**Warning:**

Contains cobalt, which is considered a hazardous substance.

During machining or heat treatment of CoCr material (e.g., milling, grinding), dust and fumes may be generated and pose a health risk if inhaled.

Therefore, appropriate measures must be taken to reduce dust and fumes.

Exposure control measures, such as the use of extraction and ventilation systems or, for short-term use, appropriate personal protective equipment for respiratory filtration.

**Warning:**

The use of different alloys in the patient's mouth may cause galvanic effects.

## 10. Storage Conditions and Disposal

The product in its original, unopened packaging does not require any special storage or handling considerations (during transport and storage). It is recommended to store packaged products in a dry place, away from direct sunlight, and at room temperature to maintain the packaging's optimal visual appearance. Products to be disposed of must be treated and decontaminated as dental surgical waste, in accordance with applicable regulations.

## 11. Warnings and General Precautions

- Before any treatment, a comprehensive clinical and radiographic evaluation of the patient is essential for their safety.
- Improper use of implant system products can lead to complications such as implant failure, local bone loss, or unsatisfactory aesthetic results.
- The contraindications for the respective implant systems must be taken into account. These can be found in the respective instructions for use of the implant systems.
- The products must only be used by dentists, surgeons, and dental technicians with specialized knowledge and skills in oral implantology; these products are sold only to dentists/surgeons and dental laboratories or upon their prescription.
- Custom-made devices are manufactured in accordance with a written prescription issued exclusively by a healthcare professional, who provides the patient-specific design characteristics for the patient under their care. These devices must not be reused under any circumstances to ensure patient safety and treatment efficacy.
- If necessary, technical advice and/or training on the proper use of Biotech Dental Digital products is provided orally, in writing, via electronic media, and/or through demonstrations. The information provided reflects the state of science and technology at the time the products are placed on the market. This does not relieve the user of the responsibility to personally test the products to ensure they are suitable for the intended purposes, indications, and procedures. The handling and use of the products are beyond the control of Biotech Dental Digital SAS and are the direct responsibility of the user. Biotech Dental Digital SAS disclaims all liability for damages resulting from such use.
- The use of incompatible components may impair the function and safety of BIOTECH DENTAL DIGITAL products. BIOTECH DENTAL DIGITAL SAS will provide no warranty or replacement in the event of the use of incompatible components.
- BIOTECH DENTAL DIGITAL implant products are—where specified—intended solely for use with specific implants/abutments and implant/abutment diameters. Their use with other implants/abutments or diameters may result in mechanical failure of system components, tissue damage, or unsatisfactory aesthetic results.
- Drills and instruments are designed for specific implant systems and implant diameters. The use of different implant systems and/or implant diameters may result in tissue damage or unsatisfactory aesthetic results. Only drills, instruments, and prosthetic components that are compatible with the system or approved by Biotech Dental Digital SAS may be used with the respective implant systems.
- Due to their small size, products may be swallowed and/or aspirated. Aspiration can lead to dyspnea and, in the worst case, asphyxiation. Therefore, products must be properly secured to prevent them from being swallowed or aspirated during intraoral use.
- When indications are listed for a specific product, please note that any indications not listed are in fact contraindicated.
- As part of our terms of sale and delivery, we guarantee the high quality of our products.
- Not all BIOTECH DENTAL DIGITAL products and services are available in every country.
- Damaged and/or corroded products must not be used.
- Products that cannot be identified or whose function is questionable, for example due to illegible markings and/or inscriptions, must not be used but must be replaced.

- In case of hypersensitivity to any of the materials listed in Section 8 Materials, refrain from using a restoration with the respective products.

### **Reporting of Serious Incidents**

Any serious incident related to the Biotech Dental Digital device must be reported to the manufacturer (Biotech Dental Digital SAS) and to the competent authority in your country.

A serious incident is defined as any malfunction or alteration in the characteristics or performance of the implant or restoration that has indirectly resulted in, could have resulted in, or could result in serious, temporary, or permanent deterioration of the patient's health or the patient's death.

### **Trademarks and Copyrights**

Protected trade names (trademarks) are not always specifically indicated. The absence of such an indication does not mean that the name is NOT a trademark. This document, including all its parts, is protected by copyright. You may download the content for its Intended Purpose, but you are prohibited from modifying or reproducing it. Any use beyond the narrow limits of copyright law is not permitted without the prior written consent of Biotech Dental Digital SAS and is subject to legal penalties.

## **12. MRI Safety Information**

### **Warning:**

BIOTECH DENTAL DIGITAL products made of titanium alloy and CoCr alloy have not been tested for safety and compatibility in an MRI environment, particularly with regard to heating, migration, or image artifacts. The safety of MRI examinations in patients who have undergone restorations with BIOTECH DENTAL DIGITAL metal products is unknown. Examining a patient wearing this device may result in injury.

## **13. More information**

Product information sheets for Biotech Dental Digital products are available for download in PDF format on the order website. You may contact customer service to have a paper copy included with your order or sent to you later, within a maximum of seven calendar days following receipt of your request, at no additional cost.

Further information on all BIOTECH DENTAL DIGITAL products can be found in the current BIOTECH DENTAL DIGITAL product catalog. This document is available on the ordering website.

Basic UDI-DI numbers do not apply because BIOTECH DENTAL DIGITAL products are custom-made devices.

## 14. Explanation of symbols

The following symbols, in accordance with ISO 15223-1, may appear on the device labeling or in the accompanying information. Please refer to the device labeling or accompanying information for the applicable symbols.


|  |   |
|--|---|
|   | Consult the instructions for use                                    |
|   | Cautin, follow the warnings   |
|   | Medical device  |
|   | Catalogue number  |
|   | Serial number   |
|   | Non-sterile   |
|   | Date of manufacture   |
|   | Do not re-use   |
|   | Do not use if packaging is damaged and consult instructions for use |
|   | Manufacturer  |
|  | Contains hazardous substances                                       |

Table 3: Explanation of symbols on appliance labels

## 15. Contact

### Manufacturer BIOTECH DENTAL DIGITAL SAS

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